Komet Dental
 Gebr. Brasseler GmbH & Co. KG
 Trophagener Weg 25
 32657 Lemgo · Germany
 Telefon +49 (0) 5261 701-0
 Telefax +49 (0) 5261 701-289
 info@kometdental.de
 www.kometdental.de

RepairPost
 © 12/2015, Gebr. Brasseler GmbH & Co. KG
301109 - Rev 12/2015

Stiftnsysteme, Post systems, Systèmes de reconstitution, Sistemas de postes, Sistemi di perni






Abb. 1 / Fig. 1

Abb. 2 / Fig. 2

Abb. 3 / Fig. 3

Abb. 4 / Fig. 4

Abb. 5 / Fig. 5

Trepanbohrer L5 · Trepan bur L5 · Tenon Repair L5 · Trépano L5 · Punta per trapanare L5

| | | | | |
|--|--|--|--|---|
| | | | | Rostfreier Stahl Stainless Steel Acier inoxydable Acero inoxidable Acciaio inossidabile |
| 114.204.2 ○ _{opt.} 1.000 min ⁻¹ ○ _{max.} 20.000 min ⁻¹ | 114.204.2 ○ _{opt.} 1.000 min ⁻¹ ○ _{max.} 20.000 min ⁻¹ | 114.204.2 ○ _{opt.} 1.000 min ⁻¹ ○ _{max.} 20.000 min ⁻¹ | 114.204.2 ○ _{opt.} 1.000 min ⁻¹ ○ _{max.} 20.000 min ⁻¹ | |

Trepanbohrer L7/L9 · Trepan bur L7/L9 · Tenon Repair L7/L9 · Trépano L7/L9 · Punta per trapanare L7/L9

| | | | | |
|--|--|--|--|---|
| | | | | Rostfreier Stahl Stainless Steel Acier inoxydable Acero inoxidable Acciaio inossidabile |
| 113.204.2 ○ _{opt.} 1.000 min ⁻¹ ○ _{max.} 20.000 min ⁻¹ | 113.204.2 ○ _{opt.} 1.000 min ⁻¹ ○ _{max.} 20.000 min ⁻¹ | 113.204.2 ○ _{opt.} 1.000 min ⁻¹ ○ _{max.} 20.000 min ⁻¹ | 113.204.2 ○ _{opt.} 1.000 min ⁻¹ ○ _{max.} 20.000 min ⁻¹ | |

Repair Stift · Repair post · Tenon Repair · Perno Repair · Perno Repair

| | | | |
|-----------|-----------|-----------|--|
| | | | Reintitan, Oberfläche gestrahlt pure titanium, sandblasted surface Titane pur, surface sablée Titanio puro, superficie arenada Titanio puro, superficie sabbiata |
| ② 332L5.2 | ② 332L7.2 | ② 332L9.2 | |

Tab. 1

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Repair root post system designed for teeth with non-removable fragments in the root canal

1. System components
 - See table 1

2. Indications
 - Prosthetic treatment of upper and lower anteriors and premolars where fractured posts (Fig.1) or root filling materials cannot be removed and a new post buildup is necessary (e.g. strategically important teeth, bridge abutments).

3. Contraindications
 - extremely small teeth
 - completely removed fragments
 - only for fragments with a maximum diameter of 1.7 mm

4. Clinical sequence
 Those instruments that are supplied non-sterile have to be prepared prior to first use.
 4.1 Post selection
 - Select the post size adapted to the anatomy of the tooth. Auxiliary means: radiograph of the tooth with millimeter screen and sliding calliper.
 - For stability reasons, the post has to adapt closely to the canal walls and has to be anchored in the root canal as deeply as possible (at least ½ of the fragment length).
 4.2 Post application
 - Enlarge root canal around the post fragment or the root filling which cannot be removed and prepare the auxiliary cavity.
 - For root post 332L5 drill with trepan bur 114.204.2 until stop is reached and for preparation of the auxiliary cavity drill approx. 1 - 2 mm deeper (Fig. 2).
 - If root posts 332L7 or 332L9 are used, the trepanation bore established with bur 114.204.2 has to be extended to the desired length with bur 113.204.2 (Fig. 3).
 - For this purpose the trepan bur 113.204.2 has special depth markings for 7 and 9 mm lengths, respectively.
 Optimal speed: 500 - 1.000 rpm using light pressure and spray coolant
 - Clean and dry root canal
 - Check post seat
 - Cement post with phosphate cement, glass ionomer cement or self-curing fixation composite (Fig. 4).
 For further details, please refer to the instructions for use of the Composite System.
 - Coat the inner and outer surface of the post shank with cement or composite
 - Apply post and hold in position with cotton roll or finger pressure
 - After curing, remove cement in excess
 - Build up core with composite or glass ionomer cement. For further details, please refer to the instructions for use of the Composite System.

5. Further treatment
 For the definite preparation, please observe the following rules:

- The crown stump has to be surrounded in a circumferential manner
 - For all-ceramic crowns we recommend shoulder preparations
 The crown has to surround the post reinforced root in an at least 2 mm wide parallel zone so that forces exerted are distributed over the entire tooth

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Instructions for use
 Please read carefully and retain for future reference.

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- Impression and temporary restoration in the usual manner
 - Manufacture and insert crown
 - Due to the biological, mechanical and static aspects of the definitive prosthesis, the long-term success of the restoration can be improved by replacing the temporary appliance with the definitive prosthesis as soon as possible.

6. Recommendations for disinfection and sterilization
 Those instruments that are supplied non-sterile have to be prepared prior to first use.
 Root posts should be cleaned and disinfected with medical alcohol. Do not use sterile packed instrument inside damaged packaging. In case of sterile packed instruments, observe the date on the “use by” symbol. Root posts are single-use products. The reuse of these products poses a risk of infection and/or the safety of the products can no longer be guaranteed. The instruments are to be treated with anti-corrosive disinfecting and cleaning agents suited for rotary instruments (e. g. Komet® DC1®). Regarding recommendation for use (immersion time, concentration, suitability) of the disinfecting or cleaning agents see manufacturers’ instructions. The instruments can be reprocessed in the thermo disinfecter, provided that the agent used is suitable for rotary instruments (mildly alkaline detergent). Disinfecting and cleaning agents must be rinsed off thoroughly with water and the instruments must be dried carefully (preferably with compressed air according to the recommendations of the Commission for Hospital Hygiene and Infectious Disease Prevention of the Robert Koch Institute). Do not store instruments in wet or humid condition for a prolonged period of time. Make sure that instruments do not come in contact with each other during ultrasonic cleaning. Inspect cleaned instruments visually. Do not reuse damaged or blunt instruments. These have to be discarded. In some cases, trepan burs have to cut the post material as well which can lead to premature blunting of the blades. Damaged and blunt instruments must be separated and should not be reused.

7. Sterilization
 The sterilization of the instruments is carried out in the autoclave at 134°C. Observe the instructions provided by the manufacturer. After sterilization, check instruments for surface defects (corrosion). Do not reuse corroded instruments. The operator of medical products is responsible for ensuring that proper treatment is carried out by qualified personnel, using suitable materials and suited equipment, as recommended by the Commission for Hospital Hygiene and Infectious Disease Prevention of the Robert Koch Institute. Work instructions with regard to the correct reprocessing of instruments according to DIN EN ISO 17664 can be downloaded from our web site www.kometdental.de or requested from the manufacturer Gebr. Brasseler.

8. Storage
 Store sterile packed instruments in their original packaging, away from light and heat in a clean and dry environment. Do not store in the same room as solvents or chemicals.

9. Safety and liability
 See also “Instructions for use and safety recommendations” in the current Komet Dental catalog. The user is responsible for checking the products prior to use as to whether they are suited for the intended purpose. In case of contributory negligence by the user, Gebr. Brasseler partially or totally declines liability for all resulting damages, particularly if these are due to non-observance of our recommendations for use or warnings as well as inadvertent misuse by the user.

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