
Komet Dental
 Gebr. Brasseler GmbH & Co. KG
 Trophagener Weg 25
 32657 Lemgo · Germany
 Telefon +49 (0) 5261 701-0
 Telefax +49 (0) 5261 701-289
 info@kometdental.de
 www.kometdental.de

BKS
 © Gebr. Brasseler GmbH & Co. KG
 301507 - Rev 2016-10

Stiftsysteme, Post systems, Systèmes de reconstitution, Sistemas de postes, Sistemi di perni




Kanalformer
 Root canal reamer
 Instrument à former le canal radicaulaire
 Instrumento para formar el conducto radicular
 Alesatore canalare

Gewindeschneider
 Thread cutter
 Instrument à tarauder
 Instrumento para cortar la rosca
 Maschiatore

Schraubensetzer
 Socket wrench
 Tournevis
 Llave para tuercas
 Chiavetta giravite

Wurzelschraube
 Screw post
 Vis radicaulaire
 Tornillo radicular
 Vite raduculare

Wurzelschraube
 Screw post
 Vis radicaulaire
 Tornillo radicular
 Vite raduculare

Rostfreier Stahl · Stainless steel · Acier inoxydable · Acero inoxidable · Acciaio inossidabile

Reintitan · Pure titanium · Titane pur · Titanio puro · Titanio puro

Tab. 1

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BKS screw post system for the restoration of crowns which are partly or completely destroyed

1. System components
 See table 1

2. Indication
 Prosthetic restoration of teeth with partially or deeply damaged crown.

- Basic requirements:
- remaining substance of 2,0 mm in height around the crown (ferrule design) or surgical crown extension
- successfully concluded endodontic treatment with remaining apical root filling of 4,0 mm
- wall adaptation of the post along its entire length

3. Contraindication
 - funnel-shaped cervical canal lumen, because the post is not adapted to the canal wall along its entire length
 - subgingival caries
 - circular isogingival destruction

4. Clinical sequence
 The instruments are supplied non-sterile and have to be prepared prior to first use
- 4.1 Endodontic treatment completed
 - 4.2 Post selection
 - Adapted to the clinical situation
 For optimal stability, the screw has to be inserted into the canal as deeply as possible - ratio height of buildup to length of anchorage at least 1:1
 - 4.3 Removal of the root filling app. 4 mm short of the apex using pilot drill 183LB (fig. 1)
 - 4.4 Reaming the root canal (Fig. 2)
 - With the root canal reamer 152BKS corresponding to the lumen. Control the desired depth with depth marking rings. Use reamers in ascending order until the desired size is reached.
 - 4.5 Cutting the thread with the respective thread cutter 118BKS (Fig. 3)
 - For lubrication use rinsing liquid. Cut the thread to the depth previously enlarged - observe the depth markings. If this cannot be carried out easily, turn out thread cutter and remove chips; then cut until the pre-enlarged depth is reached.
 - 4.6 Clean and dry root canal
 - 4.7 Control seat of screw by turning in the BKS root screw without cement - Screw-in the BKS screw with socket wrench 119BKS
 - 4.8 Ream the root canal and fill up half of it from the apical end
 - 4.9 Insert root screw (Fig. 4)
 - Turn in until resistance is met (Fig. 5)
 - 4.10 If necessary, shorten screw
 - 4.11 Build up core with plastic material

5. Further treatment
 For the definite preparation, please observe the following rules:
 - The crown stump has to be surrounded in a circumferential manner
 - For all-ceramic crowns we recommend shoulder preparations

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Instructions for use
 Please read carefully and retain for future reference.

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The crown has to surround the post reinforced root in an at least 2 mm wide parallel zone so that forces exerted are distributed over the entire tooth
 - Impression and temporary restoration in the usual manner
 - Manufacture and insert crown
 - Due to the biological, mechanical and static aspects of the definitive prosthesis, the long-term success of the restoration can be improved by replacing the temporary appliance with the definitive prosthesis as soon as possible.

6. Reprocessing
 Those instruments that are supplied non-sterile have to be prepared prior to first use. Root posts should be cleaned and disinfected with medical alcohol. Do not use sterile packed instrument inside damaged packaging. In case of sterile packed instruments, observe the date on the "use by" symbol. Root posts are single-use products. The reuse of these products poses a risk of infection and/or the safety of the products can no longer be guaranteed. The instruments are to be treated with anti-corrosive disinfecting and cleaning agents suited for rotary instruments (e. g. Komet DC1). Regarding recommendation for use (immersion time, concentration, suitability) of the disinfecting or cleaning agents see manufacturers' instructions. The instruments can be reprocessed in the thermo disinfecter, provided that the agent used is suitable for rotary instruments (mildly alkaline detergent). Disinfecting and cleaning agents must be rinsed off thoroughly with water and the instruments must be dried carefully (preferably with compressed air according to the recommendations of the Commission for Hospital Hygiene and Infectious Disease Prevention of the Robert Koch Institute). Do not store instruments in wet or humid condition for a prolonged period of time. Make sure that instruments do not come in contact with each other during ultrasonic cleaning. Inspect cleaned instruments visually. Do not reuse damaged or blunt instruments. These have to be discarded.

Sterilization:
 The sterilization of the instruments is carried out in the autoclave at 134°C. Observe the instructions provided by the manufacturer. After sterilization, check instruments for surface defects (corrosion). Do not reuse corroded instruments. The operator of medical products is responsible for ensuring that proper treatment is carried out by qualified personnel, using suitable materials and suited equipment, as recommended by the Commission for Hospital Hygiene and Infectious Disease Prevention of the Robert Koch Institute. Information on proper reprocessing of instruments according to DIN EN ISO 17664, the present instructions for use, as amended, and more detailed information can be downloaded from our web site www.kometdental.de or requested from the manufacturer Gebr. Brasseler.

7. Storage
 Store sterile packed instruments in their original packaging, away from light and heat in a clean and dry environment. Do not store in the same room as solvents or chemicals.

8. Safety and liability
 See also "Instructions for use and safety recommendations" in the current Komet catalogue. The user is responsible for checking the products prior to use, as to whether they are suited for the intended purpose. In case of contributory negligence by the user, Gebr. Brasseler partially or totally declines liability for all resulting damages, particularly if these are due to non-observance of our recommendations for use or warnings as well as inadvertent misuse by the user.

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