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
OptiPost root post system for restoration of severely destroyed teeth

- System components**
See table 1
- Indication**
Prosthetic restoration of coronally destroyed anteriors, canines and premolars in the upper and lower jaw.
Basic requirements:
- Presence of remaining substance of 2,0 mm in height around the crown (ferrule design) or surgical crown extension required
- Successfully concluded endodontic treatment with remaining apical root filling of 4,0 mm
- Wall adaptation of the post along its entire length
- The lower retention element (up to the first retention groove) can be countersunk in the root cross-section
- Contraindication**
- Funnel-shaped cervical canal lumen, because the post can not be adapted to the canal wall along its entire length
- Subgingival caries
- Circular isogingival destruction
- Clinical sequence**
The instruments are supplied non-sterile and have to be prepared prior to first use.
- Endodontic treatment completed
- Reduction of the clinical crown up to app. 2 mm above the cemento-enamel junction, proximally
- Selection of post corresponding to the tooth to be restored (see chart) and the clinical situation
- In exceptional cases where tooth size and post size do not correspond, another post size of the system can be chosen. For stability reasons, take care that the post offers adequate wall adaptation and that it is anchored in the root canal as deeply as possible (2/3 of the post length).
- Removal of the root filling to the depth required with the pilot bur 183LB, according to the chosen post size.
- Mark the required preparation length on the pilot drill with a stopper.
- Preparing the post canal (Fig. 2)

- For OptiPosts 279.1-3 and 280.3-5 initial step-by-step enlarging of the canal entrance with the corresponding bladed steel pilot drill 29A.204.1-3 and 30A.204.3-5.
- Use the Optipost pilot drill in an intermittent mode and rinse canal in between, for example with 1% NaOCl, followed by distilled water. Dry with paper points.
- Enlarging the post canal (Fig. 3)
- With OptiPost placement drill 27D.204.1-5 and 28D.204.1-5.
- Cement post with phosphate cement, glass ionomer cement or self or dual curing composite (Fig. 4). Important note: Do not use purely light curing composite or purely light curing bonding for adhesive fixation of the post. Prior to applying the composite, it is recommended to cover the post in a suited bonding material to assure adhesion of the composite. Follow the manufacturer's instructions for the Composite System.
- Remove excessive cement
- Build up core with composite

- Further treatment**
For the definite preparation, please observe the following rules:
- The crown stump has to be surrounded in a circumferential manner
- For all-ceramic crowns we recommend shoulder preparations
The crown has to surround the post reinforced root in an at least 2 mm wide parallel zone so that forces exerted are distributed over the entire tooth
- Impression and temporary restoration in the usual manner
- Manufacture and insert crown
- Due to the biological, mechanical and static aspects of the definitive prosthesis, the long-term success of the restoration can be improved by replacing the temporary appliance with the definitive prosthesis as soon as possible.

- Reprocessing**
The instruments which are supplied non-sterile have to be prepared prior to first use. Root posts should be cleaned and disinfected with medical alcohol. Do not use sterile packed instrument inside damaged packaging. In case of sterile packed instruments, observe the date on the "use by" symbol. Root posts are single-use products. The reuse of these products poses a risk of infection and/or the safety of the products can no longer