Komet Dental
Gebr. Brasseler GmbH & Co. KG
Trophagener Weg 25
32657 Lemgo · Germany
Telefon +49 (0) 5261 701-0
Telefax +49 (0) 5261 701-289
info@kometdental.de
www.kometdental.de

ER Heraplat
© Gebr. Brasseler GmbH & Co. KG
302797 - Rev 2016-10

Stiftsysteme, Post systems, Systèmes de reconstitution, Sistemas de postes, Sistemi di perni

Größe Size Taille Tamaño Misura 050	Farbe: gelb Colour: yellow Couleur: jaune Color: amarillo Colore: giallo	Größe Size Taille Tamaño Misura 070	Farbe: orange Colour: orange Couleur: orange Color: naranja Colore: arancione	Größe Size Taille Tamaño Misura 090	Farbe: rot Colour: red Couleur: rouge Color: rojo Colore: rosso	Größe Size Taille Tamaño Misura 110	Farbe: blau Colour: blue Couleur: bleu Color: azul Colore: blu	Werkstoff Material Matériau Material Materiale	

Pilotbohrer · Pilot drill · Foret pilote · Fresa piloto · Fresa pilota

				Rostfreier Stahl Stainless Steel Acier inoxydable Acero inoxidable Acciaio inossidabile
Größe Size Taille Tamaño Misura	Farbe: gelb Colour: yellow Couleur: jaune Color: amarillo Colore: giallo	Größe Size Taille Tamaño Misura	Farbe: orange Colour: orange Couleur: orange Color: naranja Colore: arancione	

Erweiterer · Reamer · Elargisseur · Ensanchador · Allargacanalì

				Rostfreier Stahl Stainless Steel Acier inoxydable Acero inoxidable Acciaio inossidabile
Größe Size Taille Tamaño Misura	Farbe: gelb Colour: yellow Couleur: jaune Color: amarillo Colore: giallo	Größe Size Taille Tamaño Misura	Farbe: orange Colour: orange Couleur: orange Color: naranja Colore: arancione	

Erweiterer, lang · Reamer, long · Elargisseur, long · Ensanchador, largo · Allargacanalì, lungo

				Rostfreier Stahl Stainless Steel Acier inoxydable Acero inoxidable Acciaio inossidabile
Größe Size Taille Tamaño Misura	Farbe: gelb Colour: yellow Couleur: jaune Color: amarillo Colore: giallo	Größe Size Taille Tamaño Misura	Farbe: orange Colour: orange Couleur: orange Color: naranja Colore: arancione	

Universal-Tiefenlehre · Universal depth gauge · Calibre de profondeur universal · Calibre de profundidad universal · Calibre di profondità universale

				Rostfreier Stahl Stainless Steel Acier inoxydable Acero inoxidable Acciaio inossidabile
Größe Size Taille Tamaño Misura	Farbe: gelb Colour: yellow Couleur: jaune Color: amarillo Colore: giallo	Größe Size Taille Tamaño Misura	Farbe: orange Colour: orange Couleur: orange Color: naranja Colore: arancione	

Planschleifer · Root facer · Alésoir à surfacer · Abrasivo para planear · Spinardici diamantata

				Rostfreier Stahl · diamantiert Stainless Steel · diamond coated Acier inoxydable · diamanté Acero inoxidable · diamantado Acciaio inossidabile · diamantato
Größe Size Taille Tamaño Misura	Farbe: gelb Colour: yellow Couleur: jaune Color: amarillo Colore: giallo	Größe Size Taille Tamaño Misura	Farbe: orange Colour: orange Couleur: orange Color: naranja Colore: arancione	

Aufrauinstrument · Roughening instrument · Instrument à rugosifier · Asperizador · Strumento per irruvidire

				Rostfreier Stahl · diamantiert Stainless Steel · diamond coated Acier inoxydable · diamanté Acero inoxidable · diamantado Acciaio inossidabile · diamantato
Größe Size Taille Tamaño Misura	Farbe: gelb Colour: yellow Couleur: jaune Color: amarillo Colore: giallo	Größe Size Taille Tamaño Misura	Farbe: orange Colour: orange Couleur: orange Color: naranja Colore: arancione	

ER Heraplat

				Gold-Platin-Legierung Gold-platinum alloy Alliage en or-platine Aleación de oro-platino Lega in oro-platino
Größe Size Taille Tamaño Misura	Farbe: gelb Colour: yellow Couleur: jaune Color: amarillo Colore: giallo	Größe Size Taille Tamaño Misura	Farbe: orange Colour: orange Couleur: orange Color: naranja Colore: arancione	

ER Heraplat

				Gold-Platin-Legierung Gold-platinum alloy Alliage en or-platine Aleación de oro-platino Lega in oro-platino
Größe Size Taille Tamaño Misura	Farbe: gelb Colour: yellow Couleur: jaune Color: amarillo Colore: giallo	Größe Size Taille Tamaño Misura	Farbe: orange Colour: orange Couleur: orange Color: naranja Colore: arancione	

Tab. 1

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ER Heraplat posts for cast buildups

1. Indication
Prosthetic restoration of teeth with deeply or completely damaged crown.
Basic requirements:
- remaining substance of 2,0 mm in height around the crown (ferrule design) or surgical crown extension
- successfully concluded endodontic treatment with remaining apical root filling of 4,0 mm
- wall adaptation of the post along its entire length

2. Contraindication
- Funnel-shaped cervical canal lumen, because the post is not adapted to the canal wall along its entire length
- Subgingival caries
- Circular isogingival destruction
- Casting-on of non-precious metal alloys for crowns and bridges

3. Proper use
Those instruments that are supplied non-sterile have to be prepared prior to first use.
3.1 Preparatory, general procedure
- Tooth/root has been endodontically treated and prepared according to the indication
- Endodontic treatment
- Case sealing the apical third
- Case partly destroyed teeth prepare the cavity for example with diamond instrument 835KR or 801.
- Severely destroyed teeth prepare the core surface until only sound tooth substance is left and give it a concave shape using ball-shaped diamond instrument 801 (Fig. 1)
- Reduction of the clinical crown to app. 2 mm above the cemento-enamel junction, proximally
- Select the matching Heraplat post size
- Use the radiograph to identify a suitable post diameter and post length according to the anatomic features of the tooth
- Post should extend 4 mm short of the apex
- Removal of the root filling to 4 mm short of the apex with pilot drill 183LB
- Select the post size
- Mark the required preparation length on the pilot bur using stop (Fig. 2)
- Enlarge root canal according to the selected Heraplat post
- With the respective reamer 196; if clinical crown is still available use reamer 196L (Fig. 3).
- Mark the required preparation length on the reamer with depth gauge 74L12
- Reaming should be carried out in an intermittent way with minor pressure until required depth is reached. In the meantime the root canal should be rinsed and the reamer should be cleaned with brush 9791
- Establishing the central inlay cavity (maximum depth 1 - 2 mm)
- For antirotational enhancement and stability prepare a longitudinal oval tapered auxiliary cavity (Fig. 4) with root facer 120D.204.030 or cylindrical diamond instruments, e.g. 835KR
- Check post seat by inserting the selected post

3.2 Direct modellation of the buildup
- Taking impression of the auxiliary cavity
- Insulate auxiliary cavity with silicone or vaseline oil
- Apply polymerizable acrylic material (e.g. Palavit G) on the post in the area of the notches
- Insert post
- Fill auxiliary cavity with modellation acrylics (Fig. 5)
- Remove post and reinsert for control
- If necessary, shorten root post
- Build up core (Fig. 6) with burn-out modelling resin, either freehand or with auxiliary devices (e.g. Frasco-Strip crown)
- Allow acrylics to cure
- Initial preparation of the core
3.2.4 Remove modelled post buildup for further processing in the laboratory and provide root with a temporary post during the laboratory work.

3.3. Temporary treatment
- Temporary treatment during laboratory work with ER TMP posts P75L14, P422L14, P76L14, P77L14. Note: The ER TMP post has been designed in particular for temporary treatment and is not suited for permanent buildups. Note: The post tip of the ER TMP post is 2 mm shorter than the Heraplat post tip. If the post is inserted in the root canal a hollow space of 2mm will remain apically. This hollow space offers the advantage that the tooth can be provided with a medicated inlay by an inflammation needs to be cured.
- Apply ER TMP post with temporary cement
Note: If the patient wears the temporary appliance only for a short time, no fixation agent is necessary. The retention ring on the post shank assures firm seat of the post. For antibacterial protection the post shank can be coated with calcium hydroxide.
- Remove excessive cement
- Fabrication and cementation of temporary crown
- Cast post buildup in the laboratory

Technical data of Heraplat posts:
- Gold-platinum alloy which can be cast on, composition: 61.0% gold, 23.8% platinum, 15.0% palladium, rest rhodium.
- Melting interval: 1.360 - 1.460°C
- Hardness: 250 HV
- Density: 18.0 g/cm³
- Provide the buildup carried out by the dentist with a cast-on and embed

Note: Only silicate or phosphate bonded investments should be used (plaster-bonded investments are not suited).

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Instructions for use

Please read carefully and retain for future reference.

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- Preheat muffle (buildup burns out) according to the instructions of the respective alloy manufacturer
Note: Preheating temperature like alloy to be cast on, max. 750°C
- Preheating temperature: e.g 700°C for Maingold SG, 700°C for Bio Maingold SG, 700°C for Hera SG, (according to the instructions of the respective alloy manufacturer).
- Casting
Note: precious metal alloys with a casting temperature of ≤ 1.200°C can be cast on
Procedure and the proper casting temperature according to the instructions of the manufacturer
- Casting-on temperature: e.g. 1.060°C for Maingold SG, 1.050°C for Bio Maingold SG, 1.000°C for Hera SG, (according to the instructions of the respective alloy manufacturer).
- Non-precious metal alloys cannot be cast on due to the high casting temperature
- Remove cast from the investment (Fig. 7)
3.5. Insert cast (post buildup)
The finished post buildup fabricated in the laboratory is passed on to the dentist
- Remove temporary appliance (temporary crown with ER TMP post)
- Clean root canal and auxiliary cavity
Note: If necessary, remove cement residues from the canal with reamer 196 (or 196L), preferably manually.
- Rinse root canal for example with 1% NaOCl, followed by distilled water. Dry with paper points.
- Seat control
- Insert cast post buildup in the root and check seat
- Eliminate irregularities, if any, with suited instruments
- Mechanical conditioning of the canal wall with roughening instrument 196D
- Turn instrument manually in the root canal without pressure (3 - 5 turns) (Fig. 8)
- Rinse root canal for example with 1% NaOCl, followed by distilled water. Dry with paper points.
- Cement completed post buildup (Fig. 9)
- Mix suited fixation agent to a creamy consistency (e.g. zinc phosphate cement, glass ionomer cement).
- Coat post shank with fixation agent.
- Insert post buildup slowly into the canal or auxiliary cavity with moderate pressure so that excessive cement can escape and tighten by finger pressure.
- After curing remove excessive fixation agent

4. Further treatment
For the definite preparation, please observe the following rules:
- The core has to be circularly surrounded by the crown
- For all-ceramic crowns we recommend the shoulder preparation
The crown has to surround the post reinforced root in an at least 2 mm wide parallel zone so that forces exerted are distributed over the entire tooth
- Impression and temporary restoration in the usual manner
- Fabricate and insert crown
- Due to the biological, mechanical and static aspects of the definitive prosthesis, the long-term success of the restoration can be improved by replacing the temporary appliance with the definitive prosthesis as soon as possible

5. Reprocessing
Those instruments that are supplied non-sterile have to be prepared prior to first use.
Root posts should be cleaned and disinfected with medical alcohol. Do not use sterile packed instrument with damaged packaging. In case of sterile packed instruments, observe the use of the "use by" symbol. Root posts are single-use products. The reuse of these products poses a risk of infection and/or the safety of the products can no longer be guaranteed. The instruments are to be treated with anti-corrosive disinfecting and cleaning agents suited for rotary instruments (e. g. Komet DC1). Regarding recommendation for use (immersion time, concentration, suitability) of the disinfecting or cleaning agents see manufacturers' instructions. The instruments can be reprocessed in the thermo disinfectant, provided that the agent used is suitable for rotary instruments (mildly alkaline detergent). Disinfecting and cleaning agents must be rinsed off thoroughly with water and the instruments must be disinfected (preferably with compressed air according to the recommendations of the Commission for Hospital Hygiene and Infectious Disease Prevention of the Robert Koch Institute). Do not store instruments in wet or humid conditions for a prolonged period of time. Make sure that instruments do not come in contact with each other during ultrasonic cleaning. Inspect cleaned instruments visually. Do not reuse damaged or blunt instruments. These have to be discarded.

Sterilisation
Sterilisation of the instruments is carried out in the autoclave at 134°C. Observe instructions provided by the manufacturer. After sterilization, check instruments for defects on the surface (corrosion). Do not reuse corroded instruments. The operator of medical products is responsible for seeing that proper treatment is carried out by qualified personnel, using the appropriate materials and suited equipment, as recommended by the Commission for Hospital Hygiene and Infectious Disease Prevention of the Robert Koch Institute. Work instructions with regard to proper reprocessing of instruments according to DIN EN ISO 17664, the present instructions for use, as requested, and more detailed information can be downloaded from our web site www.kometdental.de or amended from the manufacturer Gebr. Brasseler.

6. Storage
Store sterile packed instruments in their original packaging, away from light and heat in a clean and dry environment. Do not store in the same room as solvents or chemicals.

7. Safety and liability
See also "Instructions for use and safety recommendations" in the current Komet dental catalogue. The user is responsible to check the products prior to use whether they are suited for the intended purpose. In case of contributory negligence by the user, Gebr. Brasseler partially or totally declines liability for all resulting damages, particularly if these are due to non-observance of our recommendations for use or warnings as well as inadvertent misuse by the user.

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