
 Stiftsysteme, Post systems, Systèmes de reconstitution, Sistemas de postes, Sistemi di perni

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EN

ER ELD posts for cast buildups

1. Indication
 Prosthetic restoration of teeth with deeply damaged or destroyed crown.
 Basic requirements:
 - Remaining substance of 2,0 mm in height around the crown (ferrule design) or surgical crown extension
 - Successfully concluded endodontic treatment with remaining apical root filling of 4,0 mm
 - Wall adaptation of the post along its entire length

2. Contraindication
 - Funnel-shaped cervical canal lumen, because the post is not adapted to the canal wall along its entire length
 - Subgingival caries
 - Circular isogingival destruction
 - Cast-on crowns and bridges made of non-precious metal alloys
 - extremely small teeth or roots
 - casting on precious metal alloys with a melting point of >1.000°C

3. Proper use
 3.1 Preparatory, general procedure
 3.1.1 Tooth/root has been endodontically treated and prepared according to the indication.
 - Endodontic treatment
 - Sealing the apical third
 - With partly destroyed teeth prepare the cavity for example with diamond instrument 835KR or 801.
 - With severely destroyed tooth crowns prepare the core surface until only sound healthy tooth substance is left and give it a concave shape using ball-shaped diamond instrument 801 (Fig. 1)
 - Reduction of the clinical crown to app. 2 mm above the cemento-enamel junction, proximally
 3.1.2 Select the matching ELD post size
 - use the radiograph to identify a suitable post diameter and post length according to the anatomic features of the tooth
 - Post should extend app. 2 - 4 mm short of the apex
 3.1.3 Removal of the root filling to 2 - 4 mm short of the apex with suited drill 183LB corresponding with the selected post size 050, 070, 090 or 110
 - Mark the required preparation length on the pilot bur using stop 9114 (Fig. 2)
 - Engine speed: 1.000 - 5.000 rpm
 3.1.4 Enlarge root canal according to the selected ELD post
 - with the reamer 196; if clinical crown is still available use reamer 196L (Fig. 3).
 - Mark the required preparation length on the reamer with depth gauge 74L12
 - Engine speed: 1.000 - 2.000 rpm
 - Reaming should be carried out in an intermittent way with minor pressure until required depth is reached. In the meantime the root canal should be rinsed and the reamer should be cleaned with brush 9791
 - preferably to be enlarged manually
 3.1.5 Establishing the central inlay cavity (maximum depth 1 - 2 mm)
 - For antirotational enhancement and stabilization prepare a longitudinal oval tapered auxiliary cavity (Fig. 4)
 - If root cross-section is large enough use root facer 120D.204.030
 - Engine speed: 2.000 rpm
 - or with cylinder diamond instruments, e.g. 835KR
 3.1.6 Check post seat by inserting the selected post.
 3.2 Direct modulation of the buildupFor enhancement of retention the ELD post can be provided with vertical and horizontal grooves in the area of the cast-on.
 3.2.1 Taking impression of the auxiliary cavity
 - Insulate auxiliary cavity with silicone or vaseline oil.
 - Apply polymerizable acrylic material (e.g. Palavit G) on the post in the area of the notches
 - Insert post
 - Fill auxiliary cavity with modulation acrylics (Fig. 5)
 - Remove post and reinsert for control
 3.2.2 If necessary, shorten root post
 3.2.3 Build up core (Fig. 6)
 - with burn-out modulation acrylics
 - either freehand or with auxiliary devices (e.g. Frasco-Strip crown)
 - Allow acrylics to cure
 - Initial preparation of the core
 3.2.4 Remove modelled post buildup for further processing in the laboratory and provide root with a temporary post during the laboratory work.
 3.3. Temporary treatment
 3.3.1 Temporary treatment during laboratory work with ER TMP posts P75L14, P422L14, P76L14, P77L14Note: The ER TMP post has been designed in particular for temporary treatment and is not suited for permanent buildups.Note: The post tip of the ER TMP post is 2 mm shorter than the ELD post tip. If the post is inserted in the root canal a hollow space of 2mm will remain apically. This hollow space offers the advantage that the tooth can be provided with a medicated inlay if an inflammation needs to be cured.
 3.3.2 Apply ER TMP post with temporary cement
 Note: If the patient wears the temporary appliance only for a short time, no fixation agent is necessary.

The retention ring on the post shank assures firm seat of the post. For antibacterial protection the post shank can be coated with calcium hydroxide.
 3.3.3 Remove excessive cement
 3.3.4 Fabrication and cementation of temporary crown
 3.4 Cast post buildup in the laboratoryTechnical data of ELD posts:
 - Silver-palladium alloy which alloys can be cast on provided that they have a melting range lower than 1.085°C
 Composition: 52.5% silver, 35.0% palladium, 5.0% platinum, 7.0% copper, 0.5% zinc
 - Melting interval: 1.085 - 1.170°C
 - Hardness: 210 HV
 - Density: 11,1 g/cm³
 3.4.1 Provide the buildup carried out by the dentist with a cast-on and embed.Note: Only silicate or phosphate bonded investments should be used (plaster-bonded investments are not suited).
 3.4.2 Preheat muffle (buildup burns out) according to the instructions of the respective alloy manufacturer.Note: Preheating temperature like alloy to be cast on, max. 700°C
 - Preheating temperature: e.g. 700°C for Degulor M, 700°C for Biocrown IV plus, 700°C for Stabilor IV plus (according to the instructions of the respective alloy manufacturer)
 3.4.3 CastingNote: precious metal alloys with a casting temperature of ≤ 1,000°C can be cast on. Procedure and appropriate melting interval according to the instruction of the respective alloy manufacturer.
 - Melting interval: e.g. 900 - 970°C Degulor M, 900 - 970°C Biocrown IV plus, 870 - 925°C Stabilor IV plus, (according to the instructions of the respective alloy manufacturer).
 - Non-precious metal alloys cannot be cast on due to the high casting temperature.
 3.4.4 Remove cast from the investment (Fig. 7).
 3.5. Insert cast (post buildup).The finished post buildup fabricated in the laboratory is passed on to the dentist.
 3.5.1 Remove temporary appliance (temporary crown with ER TMP post).
 3.5.2 Clean root canal and auxiliary cavity.Note: If necessary, remove cement residues from the canal with reamer 196 (or 196L), preferably manually.
 Rinse canal, for example with 1% NaOCl, followed by distilled water. Dry with paper points.
 3.5.3 Seat control
 - Insert cast post buildup in the root and check seat
 - Eliminate irregularities, if any, with suited instruments.
 3.5.4 Mechanical conditioning of the canal wall with roughening instrument 196D.
 - turn instrument manually in the root canal without pressure (2-3 turns) (Fig. 8).
 3.5.5 Rinse root canal and dry.
 Rinse canal, for example with 1% NaOCl, followed by distilled water. Dry with paper points.
 3.5.6 Cement completed post buildup (Fig. 9).
 - Mix suited fixation agent to a creamy consistency (e.g. zinc phosphate cement, glass ionomer cement).
 - Coat post shank with fixation agent.
 - Insert post buildup slowly into the canal or auxiliary cavity with moderate pressure so that excessive cement can escape and tighten by finger pressure.
 - After curing remove excessive fixation agent.

4. Further treatment
 For the definite preparation, please observe the following rules:
 - The core has to be circularly surrounded by the crown
 - For all-ceramic crowns we recommend the shoulder preparation
 The crown has to surround the post reinforced root in an at least 2 mm wide parallel zone so that forces exerted are distributed over the entire tooth
 - Impression and temporary restoration in the usual manner
 - Fabricate and insert crown
 - Due to the biological, mechanical and static aspects of the definitive prosthesis, the long-term success of the restoration can be improved by replacing the temporary appliance with the definitive prosthesis as soon as possible.

5. Recommendations for maintenance
 All system components are supplied in non-sterile condition. Therefore posts and instruments have to be disinfected, cleaned and sterilized prior to each use. Root posts are single-use products. The reuse of these products poses a risk of infection and/or the safety of the products can no longer be guaranteed. The instruments are to be disinfected with antiseptic and disinfection agents for rotary instruments (e.g. Komet DC1). Regarding recommendations for use (immersion time, concentration, suitability) of the cleaning or disinfection agents see manufacturers' instructions. The instruments can be reprocessed in the thermo disinfectant, provided that the agent used is suitable for rotary instruments (mildly alkaline detergent). Thoroughly remove all residues of cleaning and disinfection agents with water and allow to dry completely (e.g. by means of air stream). Do not store instruments in wet or humid drying conditions for a prolonged period of time. Make sure that instruments do not come in contact with each other during ultrasonic cleaning. Inspect cleaned instruments visually. Any damaged or blunt instruments are to be rejected and must no longer be used.
 Sterilisation:
 Sterilization of the instruments is carried out in the autoclave at 134°C. Observe instructions provided by the manufacturer. After sterilization, check instruments for defects on the surface (corrosion). Do not reuse corroded instruments.
 The operator of medical products is responsible for seeing that proper treatment is carried out by qualified personnel using the appropriate materials and suited equipment. Work instructions with regard to proper