GEBR. BRASSELER GmbH & Co. KG

Trophagener Weg 25 · 32657 Lemgo
Postfach 160 · 32631 Lemgo · Germany

Telefon +49(0)5261 701-0

Telefax +49(0)5261 701-289

0205.303969.DRÖ6

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GB Dentin Scratcher

(instrument with green ring and uncoated tip)

1. Indication

Preparation/roughening of dentine directly prior to the adhesive fixation of a restoration.

2. Contra-indications

- Preparations during which the front end of the instruments is used (for example when opening carious lesions).
- Treatments during the course of which the instrument might jam. Do not use as a lever.

3. Application

- The instrument should rotate at a very low speed ($\varnothing_{rot.}$ 500 rpm) when applied to the dentin surface to be roughened. Roughen uniformly. Thoroughly remove all contaminating substances (for example saliva, residues of imprint material).
- The uncoated tip of the instrument serves to protect prepared surfaces that are not to be roughened (e.g. when roughening a crown core, the crown margin should not be roughened).
- Subsequently apply an adhesive restoration the usual way.
- Always work with sufficient spray cooling (at least 50 ml/min.).
- Avoid excessive contact pressure (>2N).
- Do not exceed the maximum permissible speed indicated on the package.

4. Recommendations for maintenance

The instrument is to be treated with anticorro-

sive disinfecting and cleaning agents for rotary instruments (e.g. Komet DC1). Regarding recommendations for use (immersion time, concentration, suitability) of the disinfecting or cleaning agents please refer to the instructions provided by the manufacturers. Thoroughly remove all residues of disinfecting and cleaning agent with water and allow to dry completely (e.g. by air stream). Do not store instrument in wet or humid condition for a prolonged period of time. Make sure that instruments do not come in contact with each other during ultrasonic cleaning. Inspect cleaned instruments visually. Any damaged or blunt separating strips are to be rejected and discarded.

Sterilisation of the separating strips is carried out in the autoclave (134°C), in the chemiclave (126°C) or in the hot air sterilizer (180°C) according to conventional methods. Observe instructions provided by the manufacturer. After completed sterilization check surfaces for corrosion. Discard any corroded instruments.

5. Safety and liability

Worn and damaged (i.e. bent etc.) instruments are to be discarded and immediately replaced by new ones. The above instructions regarding the handling of the instruments, appropriate cooling and contact pressure are to be strictly observed. Only use the instruments according to their intended purpose. The non-observance of these safety instructions can lead to damage to the power system or injuries. It is the responsibility of the user to check prior to use that the instruments are suitable for the intended purpose. GEBR. BRASSELER excludes all liability or compensation for possible damage caused by non-observance of these safety instructions.

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© 2005, GEBR. BRASSELER GmbH & Co. KG