
ER DentinPost X
 © Gebr. Brasseler GmbH & Co. KG
 304395 - Rev 2018-10

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ER DentinPost X for the restoration of teeth with plastic materials

- System overview
 See table 1
- Indication
 Prosthetic reconstruction of teeth with destroyed crown. General requirements:
 - Wall adaptation of the post along its entire length
 - The lower retention element (up to the first retention groove) can be countersunk in the root cross-section
 - Residual circular substance of 2.0 mm at the crown is available (Ferrule design), otherwise, the crown has to be extended surgically
 - Endodontic treatment completed successfully, with 4 mm of apical root filling left
- Contraindications
 - In case of cervical funnel-shaped canal lumen, wall adaptation of the post along its entire length cannot be achieved.
 - Subgingival caries
 - Circular isogingival destruction
- Clinical sequence/proper use
 - Those instruments that are supplied non-sterile have to be prepared prior to first use
 4.1 General preparatory procedure:
 - The tooth/root has been prepared endodontically according to the indication.
 - Endodontic treatment has been completed
 - Sealing of the apical third
 - Select post according to the tooth to be restored and the clinical situation
 - By means of radiograph, determine the post diameter and length according to the anatomic situation of the tooth
 - The post should end approx. 4 mm short of the apex.
 - For stability reasons, the post has to adapt to the wall and has to be anchored in the root canal as deeply as possible (at least 2/3 of the post length)
 - Removal of the root filling up to approx. 4 mm short of the apex using pilot drill 183LB
 - Select post size 050, 070, 090 or 110
 - By means of a stop, mark the required preparation length on the pilot drill (fig. 1)
 - Prepare the central inlay cavity (maximum depth 1-2 mm)
 - For anti-rotational protection and stabilization of the restoration, prepare an oval, tapered auxiliary cavity (fig. 2)
 - By slightly tilting the root facer
 - With root facer 120D.204.030
 - Prepare inlay cavity with anti-rotational groove (fig. 3)
 - With cylindrical diamond grinding instruments, e.g. 835KR
 - Enlarge root canal according to the selected post
 - Using the corresponding enlarger 196 (fig. 4)
 - Select depth measuring gauge 45L9 or 45L12 to match the root post and place onto the enlarger
 - Enlargement should be carried out in an intermittent way with low pressure until the required depth is reached
 - During this process, the root canal should be rinsed and the reamer should be cleaned with the brush 9791 from time to time
 - The enlargement is preferably done manually
 - Rinse canal, for example with 1% NaOCl, followed by distilled water. Dry with paper points
 - Check the proper fit of the prepared post bed by inserting the chosen post
 - Eliminate irregularities, if any, with enlarger 196
 - If necessary, take a radiograph for control
 - The ER DentinPost X is radiopaque (but the contrast is not as pronounced as with ceramic or metallic posts)
 4.2 Insertion of post and reconstruction
 4.2.1 Mechanical conditioning of the canal wall with roughening instrument 196D:
 - Turn instrument manually in the root canal without applying pressure (2-3 times) (fig. 5)
 4.2.2 Rinse root canal and inlay cavity and dry:
 - Rinse canal, for example with 1% NaOCl, followed by distilled water. Dry with paper points.
 - Clean the post surface: degrease with medical alcohol and dry by means of air jet
 - The post is fixed using an adhesive method
 - Chemical conditioning of the hard dental substance (enamel and dentin). Follow the instructions for the Composite System to be used.
 - Coat post with fixing composite. Follow the instructions for the Composite System.
 - Insert post slowly into the root canal with moderate pressure and fixate by applying pressure with your thumb.

- Allow composite to polymerize and remove excessive material, if necessary.
- Important: The use of purely light curing composite and bonding agents is not authorized for the fixation of the root post.
 4.2.3 Completion of the definite preparation
 - Adaptation of plastic build-up material. Follow the instructions for the Composite System (fig. 6).

- Further treatment
 For the definite preparation, please observe the following rules:
 - The crown stump has to be surrounded in a circumferential manner
 - For all-ceramic crowns we recommend shoulder preparations
 The crown has to surround the post reinforced root in an at least 2 mm wide parallel zone so that forces exerted are distributed over the entire tooth
 - Impression and temporary restoration in the usual manner
 - Manufacture and insert crown
 - Due to the biological, mechanical and static aspects of the definitive prosthesis, the long-term success of the restoration can be improved by replacing the temporary appliance with the definitive prosthesis as soon as possible.

- Preparation
 Those instruments that are supplied non-sterile have to be prepared prior to first use. Root posts should be cleaned and disinfected with medical alcohol. Do not use sterile packed instrument with damaged packaging. In case of sterile packed instruments, observe the date on the "use by" symbol. Root posts are single-use products. The reuse of these products poses a risk of infection and/or the safety of the products can no longer be guaranteed. Posts made of glass fibre reinforced composite are unsuitable for sterilization. They should be disinfected by wiping with a cloth drenched in medical alcohol. The instruments have to be disinfected with anti-corrosive cleaning and disinfecting agents suited for rotary instruments (e. g. Komet DC1). Regarding recommendation for use (immersion time, concentration, suitability) of the cleaning and disinfecting agents see manufacturers' instructions. The instruments can be reprocessed in the thermo disinfector, provided that a mildly agent is used that is suitable for rotary instruments. Cleaning and disinfecting agents must be rinsed off thoroughly with water and the instruments must be dried carefully (preferably with compressed air according to the recommendations of the Commission for Hospital Hygiene and Infectious Disease Prevention of the Robert Koch Institute). Do not store instruments in wet or humid condition for a prolonged period of time. Make sure that instruments do not come in contact with each other during ultrasonic cleaning. Inspect cleaned instruments visually. Do not reuse damaged or blunt instruments. These have to be discarded.

- Sterilization:**
 Sterilization of the instruments is carried out in the autoclave at 134°C. Observe instructions provided by the manufacturer. After sterilization, check instruments for defects on the surface (corrosion). Do not reuse corroded instruments. The operator of medical products is responsible for seeing that proper treatment is carried out by qualified personnel, using the appropriate materials and suited equipment, as recommended by the Commission for Hospital Hygiene and Infectious Disease Prevention of the Robert Koch Institute. Work instructions with regard to proper reprocessing of instruments according to DIN EN ISO 17664, the present instructions for use, as amended, and more detailed information can be downloaded from our website www.kometdental.de or requested from the manufacturer Gebr. Brasseler.

- Storage
 Store sterile packed instruments in their original packaging, away from UV radiation and heat in a clean and dry environment. Do not store in the same room as solvents or chemicals.

- Disposal
 The instruments have to be disposed of in tight, non-breakable and non-pierceable containers (protection from contamination).

- Safety and liability
 See also the instructions for use and safety recommendations in the current Komet dental catalogue. It is the responsibility of the user to check the products prior to use as to whether they are suited for the intended purpose. The user is responsible for the correct use of the instruments. In case of contributory negligence by the user, Gebr. Brasseler GmbH & Co. KG partly or totally declines liability for all resulting damages, particularly due to non-observance of our instructions for use as well as misuse by the user.

- Procedure in case of serious adverse events
 In case of serious adverse events, please proceed in line with the legal regulations in your country.

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