Komet Dental  
Gebr. Brasseler GmbH & Co. KG  
Trophagener Weg 25  
32657 Lemgo · Germany  
Telefon +49 (0) 5261 701-0  
Telefax +49 (0) 5261 701-289  
info@kometdental.de  
www.kometdental.de

AlphaKite  
305934 - Rev 2017-04

Wurzelkanalaufbereitung, Root canal preparation,  
Traitement endodontique, Preparación de canal  
radicular, Preparazione canalare

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AlphaKite System for the mechanical preparation of root canals

#### 1. Indication

The range of use for these instruments applies without

restriction to all morphologies of all types of root canals.

#### 2. Contraindication

Severely calcified canals

#### 3. Clinical sequence:

The general recommendations for a mechanical root canal preparation with NiTi files apply:

- Creation of a straight canal access
- Adherence to the Crown-Down technique or Double-Flare technique
- Gentle work in intermittent mode, without applying pressure
- Only leave instrument in the canal for a short period of time (5 – 8 seconds at most)
- To be used in lubricated canal only
- The instrument must be in rotary mode when inserted into the canal
- Lubricant used must contain EDTA
- Always rinse canal thoroughly (NaOCl, citric acid, chlorhexidine)
- Check patency by inserting a hand instrument, for example a K-file of size 010 or 015
- Monitor and record carefully how many times each instrument has been used (5 – 8 canals at most, when used in extremely curved canals, discard after one use, if necessary)
- Use in a torque-limited endodontic contra-angle, e.g. E-Drive (Komet), KaVo ENDOadvance or Sironiti/Sironiti Air+ or an endodontic motor, e.g. W&H ENTRAN
- 3.1 The use of a suitable torque-limited power system guarantees sufficient safety from torsional fracture. Do not exceed the file's load limit or its maximum fracture moment. Check instruments visually prior to use. Discard damaged or bent instruments. Prior to first use on a patient, the user should acquaint himself with the system by practicing on extracted teeth.
- 3.2 A straight access cavity of sufficient size should be created for a root canal preparation. The canal orifice

is prepared up to the curve with the access reamer AK10(L15 or L19).035 which can be identified by a green spiral. Never insert this file as far as the curve of the canal (risk of fracture or incorrect preparation). The preparation can for example be carried out according to the Crown-Down technique. Subsequent reaming is carried out with a file with a taper of 08 and/or 06, followed by a file with a taper of 04. The preparation in apical direction is carried out with an 02 tapered file. The apical area is then enlarged by one size with a file with a taper of 02. The apex can be further enlarged with an 04 file if this is required by the filling technique used.

#### 3.3 Instrument sequences

Narrow canals, yellow AK files 020: AK10L19.035 green, AK08L25.020 yellow, AK06L25.020 yellow, AK04L25.020 yellow, AK02L25.020 yellow, AK02L25.025.red and optionally once again AK04L25.020 yellow to enlarge to taper 04. Average canals, red AK files 025: AK10L19.035 green, AK08-06-04-02L25.025 red, AK02L25.030.blue, optionally: AK04L25.025 red.

Wide canals, blue AK files 030: AK10L19.035 green, AK08-06-04-02L25.030 blue, AK02L25.035 green, optionally: AK04L25.030 blue  
The files can be used for other techniques as well, for example the Double-Flare technique. A decisive factor for the success of any preparation is the constant probing with a manual file of size 010 or 015, thus not only creating a smooth path for the rotating instruments but also preventing possible blockages during preparation. The canal should also be rinsed with an adequate quantity of rinsing solution after each change of file.

#### 4. Reprocessing

All AlphaKite instruments are supplied in non-sterile condition. They therefore have to be cleaned, disinfected and sterilised prior to each use. Remove silicone stopper from the file prior to reprocessing. The instru-

ments are to be disinfected with cleaning and disinfection agents for rotary instruments (e.g. with Komet DC1). For recommendations for use (immersion time, concentration, suitability) of cleaning and disinfection agents see instructions provided by the manufacturers of these agents. Cleaning and disinfection agents must be rinsed off thoroughly with water and the instruments must be dried carefully (preferably with compressed air according to the recommendations of the Commission for Hospital Hygiene and Infectious Disease Prevention of the Robert Koch Institute). Do not store instruments for a longer period in wet or humid condition. Remove all traces of NaOCL and reduce the contact time to a minimum, i.e. a few minutes at most. Make sure that the instruments do not come in contact with each other during ultrasonic cleaning. The instruments can be reprocessed in a thermo disinfectant, provided that the agent used is suitable for rotary instruments (mildly alkaline detergent). Submit the instruments to a visual inspection. Damaged or blunt instruments must be rejected and their use discontinued.

#### Sterilization

The sterilization of the instruments is carried out in the autoclave at 134°C. Observe the instructions provided by the manufacturer. After sterilization, check instruments for surface defects (corrosion). Do not reuse corroded instruments. The operator of medical products is responsible for ensuring that proper treatment is carried out by qualified personnel, using suitable materials and suited equipment, as recommended by the Commission for Hospital Hygiene and Infectious Disease Prevention of the Robert Koch Institute. Work instructions with regard to proper reprocessing of instruments according to DIN EN ISO 17664, the present instructions for use and more detailed information, as amended, can be downloaded from our web site [www.kometdental.de](http://www.kometdental.de) or requested from the manufacturer Gebr. Brasseler.

5. Safety and liability  
See also the instructions for use and safety recommendations in the current Komet catalogue. It is the responsibility of the user to check the products prior to use as to whether they are suited for the intended purpose. The user is responsible for the correct use of the instruments. In case of contributory negligence by the user, Gebr. Brasseler GmbH & Co. KG partly or totally declines liability for all resulting damages, particularly due to non-observance of our instructions for use as well as misuse by the user.

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© Gebr. Brasseler GmbH & Co. KG