Komet Dental
Gebr. Brasseler GmbH & Co. KG
Trophagener Weg 25
32657 Lemgo · Germany
Telefon +49 (0) 5261 701-0
Telefax +49 (0) 5261 701-289
info@kometdental.de
www.kometdental.de

EasyShape System
306001 - Rev 2016-11

Wurzelkanalaufbereitung, Root canal preparation,
Traitement endodontique, Preparación de canal
radicular, Preparazione canalare

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EasyShape System for the mechanical preparation of
root canals

1. Indication

The range of use for these instruments applies without
restriction to all morphologies of all types of root
canals.

2. Contraindication

Severely calcified canals

3. Appropriate use

The general recommendations for a mechanical root
canal preparation with NiTi files apply:

- Creation of a straight canal access
- Work in “picking motion”, i.e. the instruments per-
form brief up and down movements along the canal
wall

- Only leave instrument in the canal for a short period
of time (5 - 8 seconds at most)

- To be used in lubricated canal only

- The instrument must be in rotary mode when inserted
into the canal

- Lubricant used must contain EDTA

- Always rinse canal thoroughly (NaOCl)

- Check patency by inserting a hand instrument, for
example a K-file of size 010 or 015

- Monitor and record carefully how many times each
instrument has been used (5 - 8 canals at most, when
used in extremely curved canals, discard after one use,
if necessary)

- Use in a torque-limited endomotor or an endodontic
contra-angle.

3.1 The use of a suitable torque-limited power system
guarantees sufficient safety from torsional fracture. Do
not exceed the file's load limit or its maximum fracture

3.1 The use of a suitable torque-limited power system
guarantees sufficient safety from torsional fracture. Do
not exceed the file's load limit or its maximum fracture

moment. Check instruments visually prior to use. Discard damaged or bent instruments. Prior to first use on a patient, the user should acquaint himself with the system by practicing on extracted teeth.

3.2 A straight access cavity of sufficient size should be created for a root canal preparation, for example with the Komet Endo Access Kit 4646. For the preparation of isthmuses use the endodontic round bur H1SML. The precondition for any root canal preparation is always that the patency of the canal has been checked by means of a hand file. The canal to be prepared must be free of any obstructions. The file that corresponds to the natural size of the canal determines the target size of the preparation. The root canal is cleaned, i.e. reamed mechanically to a few sizes beyond that with the ES files.

3.3 EasyShape file sequence

If the canal access zone is to be reamed to a more comfortable size, the access reaming file with taper 10 (AK10L19.035) from the Komet AlphaKite system can be used prior to the ES files. After manual probing, the mechanical preparation with the ES files starts with the white file ES04.015. The instrument is used along its entire working length, provided that this is possible without meeting too much resistance, followed by further sizes in the following order: ES05.020 yellow, ES06.025 red, ES05.030 blue, ES04.035 green and ES04.040 black. The preparation is completed when the file that had previously determined the target size of the preparation has been used, i.e. when the target size has been reached. In case of narrow canals, the target size can already be reached with the ES05.030 blue. The blade configuration of the ES files is particularly sharp. They are used in so-called "picking motion", which means that the files perform brief up and down movements in the root canal. The user will soon notice that the first two ES files 015 and 020 have to be guided with utmost care. All subsequent files can follow this path easily and shape the canal with a very satisfactory

result. A decisive factor for the success of any preparation is the constant probing with a hand file of size 010 or 015, thus not only creating a smooth path for the rotating instruments but also preventing possible blockages during preparation. The canal should also be rinsed with an adequate quantity of rinsing solution after each change of file.

4. Reprocessing

Those instruments that are supplied non-sterile have to be prepared prior to first use. The instruments have to be treated with anti-corrosive disinfecting and cleaning agents suited for rotary instruments (e. g. Komet DC1). Regarding recommendation for use (immersion time, concentration, suitability) of the disinfecting or cleaning agents, see manufacturers' instructions. Disinfecting and cleaning agents must be rinsed off thoroughly with water and the instruments must be dried carefully (for example by air stream). Do not store instruments in a wet or humid condition for a prolonged period of time. Make sure that instruments do not come in contact with each other during ultrasonic cleaning. The instruments can be reprocessed in the thermo disinfectant, provided that a mildly alkaline cleaning and disinfection agent is used that is suitable for rotary instruments. Inspect cleaned instruments visually. Do not reuse damaged or blunt instruments. These have to be discarded.

Sterilization:

The sterilization of the instruments is carried out in the autoclave at 134°C. Observe the instructions provided by the manufacturer. After sterilization, check instruments for surface defects (corrosion). Do not reuse corroded instruments. The operator of medical products is responsible for ensuring that proper treatment is carried out by qualified personnel, using suitable materials and suited equipment, as recommended by the Commission for Hospital Hygiene and Infectious Disease Prevention of the Robert Koch Institute. Work instructions with regard to the correct reprocessing of instruments according to DIN EN ISO 17664 can be

downloaded from our web site www.kometdental.de or requested from the manufacturer Gebr. Brasseler.

5. Safety and liability
See also "Instructions for use and safety recommendations" in the current Komet catalogue. The instruments may only be used for the intended use indicated. Failure to observe these safety recommendations can lead to injury or damage to the power unit. The user is responsible for checking the products prior to use as to whether they are suited for the intended purpose. In case of contributory negligence by the user, Gebr. Brasseler partially or totally declines liability for all resulting damages, particularly if these are due to non-observance of our recommendations for use or warnings as well as inadvertent misuse by the user.

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© Gebr. Brasseler GmbH & Co. KG