Komet Dental
Gebr. Brasseler GmbH & Co. KG
Trophagener Weg 25
32657 Lemgo · Germany
Telefon +49 (0) 5261 701- 0
Telefax +49 (0) 5261 701- 289
info@kometdental.de
www.kometdental.de

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Milling base

1. Product description

This milling base for the dental laboratory consists of two main components and is provided with a chucking device allowing work on cylindrical and slightly tapered laboratory replica, either manually or with the milling device. After removal of the upper part to carry out free- handed work, this can be placed back into its original position which is fixed by the lower part of the milling base.

2. Functional overview - milling base 150.555

- Suitable for cylindrical and slightly tapered laboratory replica/implants with diameters ranging from 1 to 6.5 mm
- Internal height adjuster for individually adjusting the position of the replica
- The milling base can be fixed to a magnetic or model table
- Two- piece construction to allow flexible work (milling device or free- handed)
- Quick removal and repositioning of the chuck, thanks to screw- on locking bolt
- Ergonomically shaped clamping nut
- No risk of injury caused by fast generation of heat
- Metal parts made of stainless steel

3. Start- up and appropriate use

Milling device

- Unscrew clamping nut. Choose the clamping jaw according to the diameter of the laboratory replica or implant and place into the chuck (fig.1). Screw clamping nut back on until slight resistance can be felt.
- Insert laboratory replica/implant into the clamping jaw and clamp gently by fastening the clamping nut.
- Thanks to a through- hole at the underside of the base, the height adjuster can be set to the required height with a slotted screw driver to suit the laboratory replica/implant to be worked on (fig. 2). Then fix the laboratory replica/implant in its final position with the clamping nut.
- Fix the chucking device to the base plate with the locking bolt.
- Use on a magnetic or model table.

Free- handed use

- Loosen locking bolt with 0.5 - 1 turns, pull black handle, twist by 90° and deposit on the rests provided to prevent the spring from snapping back (fig. 3).
- Remove chucking unit (fig. 4) and use for free- handed work.

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Repositioning the chucking unit

- Place chucking unit into the base plate.
- Turn the handle of the locking bolt by 90° towards the grooves (does not yet lock at this stage)
- Turn chucking unit in the base plate until the handle of the locking bolt catches
- Fasten the chucking unit with the locking bolt

4. Maintenance

- The individual elements should be cleaned separately in regular intervals.
- To this end, disassemble the milling base and remove any contamination, such as chips and/or residues of cooling agents and lubricants
- Thoroughly clean the bearing surface of the depression in the base plate and the chuck. Apply a care product containing oil (e.g. Ballistol)
- Clean the clamping jaws and the reception cone in the collet with extra care to avoid damage to the chucking mechanism.
- Whenever cooling spray is used, clean the milling base straight after use and apply a care product containing oil to avoid corrosion.

5. Disassembly/assembly

- Initial situation: complete, with mounted clamping jaw
 - Unscrew clamping nut
 - Remove clamping jaw
 - Unscrew locking bolt
 - Remove chuck
 - Unscrew height adjuster
- To reassemble the milling block, proceed in reverse order

6. Safety and responsibility

The user is responsible for checking the products prior to use to ensure that these are suited for the intended purpose. The user is responsible for the correct use of this milling base. Damaged components, such as the clamping jaws or the locking bolt, have to be discarded immediately and replaced by new ones. Inappropriate use of the milling base can lead to injury. In case of contributory negligence by the user, Gebr. Brasseler party or totally declines liability for all resulting damages, particularly due to non- observance of our recommendations for use or warnings as well as inadvertent misuse by the user.

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© 05/2013, Gebr. Brasseler GmbH & Co. KG