Komet Dental
Gebr. Brasseler GmbH & Co. KG
Trophagener Weg 25
32657 Lemgo · Germany
Telefon+49 (0)5261 701-0
Telefax+49 (0)5261 701-289
info@brasseler.de
www.brasseler.de

Frästechnikpolierer

© Gebr. Brasseler GmbH & Co. KG
306686 - Rev 2018-10

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EN

Multi-step polishing system for polishing precious and non-precious metals in the milling device

1. Indication

Multi-step polishing system for polishing precious metal and non-precious metal.

2. Proper use

- Only ever use properly maintained and cleaned power systems in good working order
 - Chucks must be cleaned in regular intervals
 - Remove all residues of previously used coolants/lubricants (e.g. milling oil) from the work piece. Check polishers prior to use to make sure they are clean
 - Insert polisher to the appropriate depth
 - Prior to use, check that the polisher is concentric and make sure it is safe to use. For your own safety, only use polishers within the scope of their intended use and wear protective glasses at all times
 - Always apply low contact pressure. Excessive pressure and speed lead to increased generation of heat and may destroy the polisher
 - Bring polisher to optimal speed prior to application
 - Best results are achieved if the indicated speed is kept
 - Never exceed the maximum permissible speed
 - Polishing to be carried out dry
- Please make sure that prior to polishing, the surface of the work piece has been ground to a sufficiently smooth level. The shape of the polisher matches that of the work piece (e.g. tapered crown). The polishers are used in several steps.

The dark brown metal polisher (pre-polishing) is merely suitable for use on non-precious metal. The brown metal polisher (polishing) and the green metal polisher (high-shine polishing) are suitable for use on precious metal and non-precious metal.

Depending on the existing surface, the previous polishing step can be omitted.

The shank of polishers dressed at our works is provided with an indication of the cone angle. Non-dressed polishers (without indication on their shanks) with sufficient polishing mass can be dressed to an angle of 0°, 1°, 2°, 4° or 6° by means of the dressing diamond 150.461M (distinguishing mark: blue point). The surface of the polisher should then be smoothed by means of dressing diamond 150.461F (red point). To avoid the formation of grooves, move polisher gently up and down in axial direction during the dressing process. The optimal speed is \odot_{max} 6,000 rpm. The maximum permissible speed is indicated on the label.

Abrasive agent: Silicon carbide (SiC)

3. Safety and liability

Please observe the general instructions for use and safety recommendations in the current Komet dental catalogue. The user is responsible for checking each polisher prior to use to ensure that it is suited for the intended purpose. Bent or out-of-concentric instruments are to be discarded and replaced by new ones immediately. Improper use may lead to damage of the power system or injury. In case of contributory negligence by the user, Gebr. Brasseler partially or totally declines liability for all resulting damages, particularly if these are due to nonobservance of our recommendations for use or warnings as well as inadvertent misuse by the user.

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