	Komet Dental
	Gebr. Brasseler GmbH & Co. KG Trophagener Weg 25 32657 Lemgo · Germany Telefon+49 (0) 526 1 701-0 Telefax+49 (0) 526 1 701-289 info@kometdental.de www.kometdental.de

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EN

Dental bur block with silicone plugs for storing sonic tips and their accessories

1. Indication
Stainless steel bur block for the storage and sterilising sonic tips and their accessories.

2. Contra-indications
Thermal cleaning and disinfection with sonic tips in place.

3. Warnings
The bur block is closed by means of a lid with clamping fit. Risk of injury!Open and close the lid with care to avoid accidental trapping of your fingers. Proceed with due care when manually moving the inner rail into the desired position. Inappropriate use can lead to injury to the user.

4. Recommendations for use
The bur block is designed for the storage and sterilisation of sonic tips.

5. Recommendations for use
- Open lid before filling the bur block.
- The sonic tips are placed into the specially provided silicone plugs in the swivelling instrument bar, shankend first (thread).

6. Maintenance
Clean, disinfect and sterilize bur block prior to first use.
- Information: In compliance with the „Hygiene Requirements for the Reprocessing of Medical Devices“, the preparation of the instruments in the washer/disinfector has been validated by an independent institution.

7. Reprocessing
Warning: The following descriptions of the reprocessing methods (thermal cleaning and disinfection) only apply to empty bur blocks. For recommendations on cleaning and disinfecting sonic tips, please refer to the relevant instructions for use enclosed with the sonic tips.

Manual reprocessing

1. Remove surface contamination thoroughly under running water. Completely remove stubborn contamination with a nylon brush below water level. Finally, rinse under running water. Place bur block in an ultra sonic device or instrument bath filled with a suitable detergent/disinfectant.
2. To avoid stains, rinse off all residues of disinfecting agent thoroughly under running water and dry the bur block (preferably with compressed air according to the recommendations of the Commission for Hospital Hygiene and Infectious Disease Prevention of the Robert Koch Institute).

Mechanical reprocessing

1. Manual pre-cleaning: Remove surface contamination thoroughly under running water. Completely remove stubborn contamination with a nylon brush below water level. Finally, rinse under running water.
2. Place the bur block in the washer/disinfector with its lid open to ensure that the inner surfaces are reached by the spray jet.
3. Put chemical detergent into the washer/disinfector, following the indications on the label and the instructions of the manufacturer of the washer/disinfector.
4. Start the Vario TD programme including thermal disinfection. Thermal disinfection takes place allowing for the A_w value and observing national provisions (prEN/ISO 15883).

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Instructions for use

Please read carefully and retain for future reference.

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EN

5. On completion of the cycle remove bur block and inserts from the washer/disinfector and dry (preferably with compressed air as recommended by the Commission for Hospital Hygiene and Infectious Disease Prevention of the Robert Koch Institute).

6. Visual examination to ensure that the bur block is clean and undamaged. If after mechanical reprocessing there are still visible residues of contamination, repeat the cleaning and disinfecting process until no visible contamination is left.
According to the recommendations of the Commission for Hospital Hygiene and Infectious Disease Prevention of the Robert Koch Institute, reprocessing should preferably be carried out mechanically. Work instructions with regard to proper reprocessing of instruments according to DIN EN ISO 17664, the present instructions for use and more detailed information, as amended, can be downloaded from our web site www.kometdental.de or requested from the manufacturer Gebr. Brasseler.

Sterilisation in the autoclave (with or without inserted sonic tips):
Warning: Prior to sterilisation of bur blocks with sonic tips in place ensure that all sonic tips to be sterilised have undergone a manual or mechanical reprocessing cycle according to the instructions of the manufacturer. For recommendations for care and maintenance as well as correct reprocessing please refer to the instructions of the manufacturer of the sonic tips. For sterilisation, the bur block has to be sealed in suitable packaging. Make sure that the packaging/foil is suitable for the chosen method of sterilisation. The packaging must be large enough to ensure that there is no pressure on the seal.

For this procedure, the following instructions have to be observed:
Steam sterilisation using a vacuum process at 134 °C in a device that complies with the provisions of DIN EN 13060; with validated processes.
- Fractionated pre-vacuum (type B)
- Sterilisation temperature: 134 °C
- Hold time: at least 5 minutes (full cycle)
- Drying time: at least 10 minutes

In order to prevent staining and corrosion, the steam must be free of particles. The recommended limits for particle contents in feed water and condensed steam are defined by standard DIN EN 13060. Make sure not to exceed the maximum capacity of the sterilizer when sterilizing several instruments. Follow the instructions of the device manufacturer.

Universally valid notes:

Observe the legal provisions regarding the reprocessing of medical products valid in your country (e.g. www.rki.de). The manufacturer confirms that the above detailed reprocessing methods are suitable for preparing the above named instrument group to enable their reuse. The person carrying out the reprocessing is responsible for ensuring that the applied method is carried out with appropriate equipment, materials and trained personnel at the reprocessing site and that it actually achieves the desired result. To guarantee this, routine controls of the validated mechanical and/or manual preparation methods are necessary. Any deviation from the above detailed process must be carefully checked by the operator to ensure effectiveness and to avoid possible adverse consequences.

8. Liability
It is the responsibility of the user to check the products prior to use to ensure that they are suited for the intended purpose. In case of contributory negligence by the user, Gebr. Brasseler partially or totally declines liability for all resulting damages, particularly due to non-observance of our recommendations for use or warnings as well as inadvertent misuse by the user.

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© Gebr. Brasseler GmbH & Co. KG