






























Rundbohrer · Round bur · Fraise boule · Fresa redonda · Rosetta  H141.104.023 6.000 min ⁻¹ 100.000 min ⁻¹					Hartmetall Tungsten carbide Carbure de tungstène Carburo di tungsteno
Rundbohrer · Round bur · Fraise boule · Fresa redonda · Rosetta  Ø 1,8 mm KI141A.205.018 1.000 min ⁻¹ 6.000 min ⁻¹		 Ø 2,3 mm KI141A.205.023 1.000 min ⁻¹ 6.000 min ⁻¹			Rostfreier Stahl Stainless Steel Acier inoxydable Acero inoxidable Acciaio inossidabile
Rundbohrer · Round bur · Fraise boule · Fresa redonda · Rosetta  K160A.104.023 6.000 min ⁻¹ 100.000 min ⁻¹					ATZ-Keramik Ceramics Céramique Cerámica Ceramica
Dreikantbohrer · Triangular bur · Fraise triangulaire · Fresa triangular · Fresa triangolare  186A.204.001 800 min ⁻¹ 100.000 min ⁻¹					Rostfreier Stahl Stainless Steel Acier inoxydable Acero inoxidable Acciaio inossidabile
Dreikantbohrer · Triangular bur · Fraise triangulaire · Fresa triangular · Fresa triangolare  Ø 1,5 mm KI186L16.204.015 1.000 min ⁻¹ 6.000 min ⁻¹					Rostfreier Stahl Stainless Steel Acier inoxydable Acero inoxidable Acciaio inossidabile
Pilotbohrer · Pilot drill · Fraise pilote · Fresa piloto · Fresa pilota  210L16.204/205.020 1.000 min ⁻¹ 6.000 min ⁻¹		 210L16.204/205.028 1.000 min ⁻¹ 6.000 min ⁻¹			Rostfreier Stahl Stainless Steel Acier inoxydable Acero inoxidable Acciaio inossidabile
Pilotbohrer · Pilot drill · Fraise pilote · Fresa piloto · Fresa pilota  K210L16.204.020 1.000 min ⁻¹ 6.000 min ⁻¹		 K210L16.204.028 1.000 min ⁻¹ 6.000 min ⁻¹	 K210L20.205.035 1.000 min ⁻¹ 6.000 min ⁻¹	 K210L20.205.042 1.000 min ⁻¹ 6.000 min ⁻¹	ATZ-Keramik Ceramics Céramique Cerámica Ceramica
Erweiterungsbohrer · Enlarger · Elargisseur · Ensanchador · Fresa allargatrice  EWB595.204.035 300 - 600 min ⁻¹ 6.000 min ⁻¹		 EWB595.204.040 300 - 600 min ⁻¹ 6.000 min ⁻¹	 EWB595.204.045 300 - 600 min ⁻¹ 6.000 min ⁻¹	 EWB595.204.055 300 - 600 min ⁻¹ 6.000 min ⁻¹	Rostfreier Stahl Stainless Steel Acier inoxydable Acero inoxidable Acciaio inossidabile
Erweiterungsbohrer · Enlarger · Elargisseur · Ensanchador · Fresa allargatrice  EWB595K.204.035 300 - 600 min ⁻¹ 6.000 min ⁻¹		 EWB595K.204.040 300 - 600 min ⁻¹ 6.000 min ⁻¹	 EWB595K.204.045 300 - 600 min ⁻¹ 6.000 min ⁻¹	 EWB595K.204.055 300 - 600 min ⁻¹ 6.000 min ⁻¹	ATZ-Keramik Ceramics Céramique Cerámica Ceramica
Gewindebohrer · Screw tap · Fraise fileté · Macho de roscar · Maschiatore  GWB596.204.035 30 min ⁻¹ 6.000 min ⁻¹		 GWB596.204.043 30 min ⁻¹ 6.000 min ⁻¹	 GWB596.204.050 30 min ⁻¹ 6.000 min ⁻¹	 GWB596.204.060 30 min ⁻¹ 6.000 min ⁻¹	Rostfreier Stahl Stainless Steel Acier inoxydable Acero inoxidable Acciaio inossidabile
Gewindebohrer · Screw tap · Fraise fileté · Macho de roscar · Maschiatore  GWB596K.204.035 30 min ⁻¹ 6.000 min ⁻¹		 GWB596K.204.043 30 min ⁻¹ 6.000 min ⁻¹	 GWB596K.204.050 30 min ⁻¹ 6.000 min ⁻¹	 GWB596K.204.060 30 min ⁻¹ 6.000 min ⁻¹	ATZ-Keramik Ceramics Céramique Cerámica Ceramica
Senker · Countersink bur · Fraise · Fresa avellanadora · Fresa countersink  SEK597.204.035 300 - 600 min ⁻¹ 6.000 min ⁻¹		 SEK597.204.043 300 - 600 min ⁻¹ 6.000 min ⁻¹	 SEK597.204.050 300 - 600 min ⁻¹ 6.000 min ⁻¹	 SEK597.204.060 300 - 600 min ⁻¹ 6.000 min ⁻¹	Rostfreier Stahl Stainless Steel Acier inoxydable Acero inoxidable Acciaio inossidabile
Senker · Countersink bur · Fraise · Fresa avellanadora · Fresa countersink  SEK597K.204.035 300 - 600 min ⁻¹ 6.000 min ⁻¹		 SEK597K.204.043 300 - 600 min ⁻¹ 6.000 min ⁻¹	 SEK597K.204.050 300 - 600 min ⁻¹ 6.000 min ⁻¹	 SEK597K.204.060 300 - 600 min ⁻¹ 6.000 min ⁻¹	ATZ-Keramik Ceramics Céramique Cerámica Ceramica
Schleimhautstanze · Tissue Punch · Clé pour la mise à nu de la vis de cicatrisation · Punzón de tejido · Mucotomo  222.644.028		 222.644.035			Rostfreier Stahl Stainless Steel Acier inoxydable Acero inoxidable Acciaio inossidabile
Schleimhautstanze · Tissue Punch · Clé pour la mise à nu de la vis de cicatrisation · Punzón de tejido · Mucotomo  SHS594.204.037 30 min ⁻¹ 800 min ⁻¹		 SHS594.204.042 30 min ⁻¹ 800 min ⁻¹	 SHS594.204.050 30 min ⁻¹ 800 min ⁻¹	 SHS594.204.060 30 min ⁻¹ 800 min ⁻¹	Rostfreier Stahl Stainless Steel Acier inoxydable Acero inoxidable Acciaio inossidabile
Trepanbohrer · Trepan · Trépan · Trépano · Fresa carotatrice  227A.204.032 800 - 1.000 min ⁻¹ 6.000 min ⁻¹		 227A.204.040 800 - 1.000 min ⁻¹ 6.000 min ⁻¹	 227A.204.050 800 - 1.000 min ⁻¹ 6.000 min ⁻¹	 227A.204.060 800 - 1.000 min ⁻¹ 6.000 min ⁻¹	Rostfreier Stahl Stainless Steel Acier inoxydable Acero inoxidable Acciaio inossidabile
Richtungs- und Tiefenindikator · Direction and depth indicator · Indicateur de direction et de profondeur · Indicador de dirección y profundidad · Indicatore di direzione e profondità  RIN598.000.035		 RIN598.000.042	 RIN598.000.050		Rostfreier Stahl Stainless Steel Acier inoxydable Acero inoxidable Acciaio inossidabile
Richtungs- und Tiefenindikator · Direction and depth indicator · Indicateur de direction et de profondeur · Indicador de dirección y profundidad · Indicatore di direzione e profondità  RIN598K.000.035		 RIN598K.000.042	 RIN598K.000.050		ATZ-Keramik Ceramics Céramique Cerámica Ceramica
Tiefenstopp · Depth stop · « Stops » de profondeur · Tope de profundidad · Arresto  599.000.L8 599.000.L10 599.000.L12 599.000.L14 599.000.L16					Rostfreier Stahl Stainless Steel Acier inoxydable Acero inoxidable Acciaio inossidabile
Tiefenstopp · Depth stop · « Stops » de profondeur · Tope de profundidad · Arresto  599K.000.L8 599K.000.L10 599K.000.L12 599K.000.L14 599K.000.L16					Peek-Kunststoff Polyetheretherketone Polyéther éther céton Polyetheretherketone Polietereetherketone
Schlüssel · Handle · Clé · Llave · Chiavetta  155A.000					Reintitan Pure titanium Titane pur Titanio puro Titanio puro
Handrad · Hand wheel · Molette · Rueda manual · Avvitatore a mano  454.000					Titan Titanium Titane Titanio Titanio

Tab. 1

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Instrument Set Implantology

1. Indication

The instrument set is intended for the preparation and enlargement of the implant site. When using a drilling template, make sure to prevent the instrument from jamming to avoid possible fracture of the instrument. The user can choose between metal and ceramic instruments.

2. Contra-indication

Danger to anatomic structures in the region of the intended intervention.

3. Treatment sequence/appropriate use

All system instructions on various surgical techniques can be found in the relevant specialist literature. Prior to the surgical intervention, patients have to be informed on the generally valid precautionary measures and rules of conduct. In addition to the bite-wing photograph and the orthopantomogramme, we recommend a computer tomographic examination to determine the exact position and depth of the drill hole. Examine the adjacent structures to exclude any danger to adjacent structures. Administer a local anaesthetic to the site where the instrument is to be used. Cut and prepare the mucosa and the periosteum according to the surgical procedure for implants. The exact position of the alveolar cavity is determined by means of the round bur (H141, KI141A or K160A) or the triangular bur 186A or KI186L16. The pilot drill 210L16 or L19 or K210Lx is applied at the site previously determined by the round bur. The pilot drill serves to drill the initial hole for the intended preparation. When using drilling templates, make sure to prevent the instrument from jamming to avoid possible fracture of the instrument. The enlarger EWB595 or EWB595K is used for enlarging the cavity in a tapered manner up to its final width. Drill in an intermittent mode, providing constant external cooling with a physiological saline solution. The external irrigation prevents thermal damage to the bone tissue and facilitates the evacuation of bone chips. Make sure that the drill does not jam or get stuck during use (increased risk of fracture). Apply low pressure during the preparation up to the required depth. Work at a speed of 300 – 600 rpm. Do not exceed the recommended speed to prevent fracture of the instrument. The screw tap GWB596 or GWB596K is used to enlarge the cavity in case of hard bones, taking care not to exceed the indicated maximum speed of 30 rpm. Alternatively, the cavity can also be enlarged manually. In case of a compact cortical bone, the cavity can be enlarged with the countersink bur SEK597 or SEK597K in that region. The Manual tissue punch (222.644) or mechanical tissue punch (SHS594) is used to perforate the gingiva at the site where the preparation is to take place, whereas the trepan bur 227A serves to remove enossal implants with a max. total length of 18 mm. If the implant is still completely seated inside the bone, remove approx. 1 – 2 mm of bone substance around the implant to enable the trepan bur to grip the implant.

4. Further treatment

The instruments are placed in the alveolar cavity and the depth is determined by means of the ring-shaped depth marks. To limit the drilling depth, the depth stop 599 or 599K is placed over the enlarger. The tools 538, 155A and 454 allow safe adaptation and use of further elements, such as fixing screws.

5. Care and maintenance

All system components mentioned in these instructions are supplied non-sterile. They therefore have to be disinfected and sterilised prior to each use. Please take note of the instructions concerning disposable instruments, where applicable. Prior to first use, prepare the instruments as described below. Only use brushes with metal-free bristles for pre-cleaning ceramic instruments as these do not leave discolorations caused by abrasion.

Cleaning/Disinfection/Sterilisation

Mechanical cleaning and disinfection

- Place the instruments in the washer/disinfector in such a way that the products are directly hit by the spray jet.
- Put chemical detergent into the washer/disinfector, following the indications on the label and the instructions of the manufacturer of the washer/disinfector.
- Start the Vario TD programme including thermal disinfection. Thermal disinfection takes place allowing for the A_s value and observing national provisions (prEN/ISO 15883).

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- On completion of the cycle remove sonic tip from the washer/disinfector and dry (preferably with compressed air as recommended by the Robert Koch Institute).
- Visual examination to ensure that the instrument is clean and undamaged. In case of residual contamination after mechanical reprocessing, repeat the cleaning and disinfecting process until no visible contamination is left. This does not apply to articles that cannot be reprocessed.

Manual cleaning and disinfection (alternative)

- Place the instruments into the ultrasonic device or instrument bath filled with detergent/disinfectant (closed lid).
- During chemical disinfection in the ultrasonic device or instrument bath, observe the instructions of the manufacturer regarding concentration and immersion time. Be sure to observe the full correct immersion time which does not start until the last instrument or system component has been placed into the bath.
- On completion of the immersion time, rinse instruments thoroughly with suitable water (preferably with demineralised water to avoid residues of lime).
- Dry instruments (preferably with compressed air as recommended by the Robert Koch Institute).
- Visual examination to ensure that the instrument is clean and undamaged. In case of residual contamination, repeat the cleaning and chemical disinfecting process until no visible contamination is left. This does not apply to articles that cannot be reprocessed. According to the recommendations of the Robert Koch Institute, mechanical cleaning and disinfection should be given preference. Information on the reprocessing of instruments in compliance with the standard DIN EN ISO 17664 can be downloaded from our website www.kometdental.de or requested directly from the manufacturer Gebr. Brasseler.

Sterilisation in the autoclave

All instruments are suitable for sterilisation. When sealing the products in the foil, make sure that the packaging is large enough to ensure that there is no pressure on the seal. The system components can be steam sterilised using a vacuum process at 134°C in a cycle according to DIN EN 13060. For this procedure, the following instructions have to be observed: Steam sterilisation using a vacuum process at 134°C in a device that complies with the provisions of DIN EN 13060; with validated processes.
- Fractionated pre-vacuum (type B)
- Sterilisation temperature: 134°C
- Hold time: at least 5 minutes (full cycle)
- Drying time: at least 10 minutes

In order to prevent staining and corrosion, the steam must be free of particles. The recommended limits for particle contents in feed water and condensed steam are defined by the standard DIN EN 13060. Make sure not to exceed the maximum capacity of the sterilizer when sterilizing several instruments. Follow the instructions of the device manufacturer. Do not reuse corroded instruments. Ceramic instruments have to be checked for superficial damages (formation of cracks) after sterilisation. Do not reuse damaged instruments due to the increased risk of fracture. Universally valid note (not applicable to disposable instruments): Observe the legal provisions regarding the reprocessing of medical products valid in your country (e.g. www.rki.de). The manufacturer confirms that the above detailed reprocessing methods are suitable for preparing the above named instrument group to enable their reuse. The operator of medical products is responsible for making sure that reprocessing is carried out by qualified personnel, using the appropriate materials and suited equipment. To guarantee this, routine controls of the validated mechanical preparation methods are required. Any deviation from the above detailed process must be carefully checked by the operator to ensure effectiveness and to avoid possible adverse consequences.

6. Safety and liability:

Worn and damaged instruments are to be discarded and immediately replaced by new ones. The above instructions regarding the handling of the instruments are to be strictly observed. Only use the instruments according to their intended purpose. The non-observance of these safety instructions can lead to injuries.

Liability:

It is the responsibility of the user to check the products prior to use to ensure that they are suited for the intended purpose. In case of contributory negligence by the user, Gebr. Brasseler partially or totally declines liability for all resulting damages, particularly due to non-observance of our recommendations for use or warnings as well as inadvertent misuse by the user.

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