Komet Dental  
 Gebr. Brasseler GmbH & Co. KG  
 Trophagener Weg 25  
 32657 Lemgo · Germany  
 Telefon +49 (0) 5261 701-0  
 Telefax +49 (0) 5261 701-289  
 info@kometdental.de  
 www.kometdental.de

**Endo Interimstand**  
 © Gebr. Brasseler GmbH & Co. KG  
**308496 - Rev 2017-07**

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**Instructions for use**  
 Please read carefully and retain for  
 future reference.

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

Bur block for the interim storage of endodontic instruments

1. Indication  
 Interim storage and cleaning of root canal instruments during the course of the treatment. With integrated length  
 scale.

2. Contra-indication  
 Unsuitable for the reprocessing of root canal instruments.

3. Appropriate use  
 - Preparation: Fill the disinfection solution (for example medical alcohol) into the bur block (Komet 595), no  
 further than the filling edge. Then place the foam insert into the bur block, exerting gentle pressure. It should be  
 flush with the upper edge of the bur block.  
 Use:  
 - For interim storage and cleaning during the treatment, the working parts of the root canal instruments are  
 placed into the foam insert drenched with disinfection solution.  
 - For a quick and precise adjustment of the length, the silicone stops of the instruments are placed in the groove  
 at the bottom of the length scale on the front side of the bur block and adjusted to the required working length.  
 - The perforations with length specifications at the upper edge of the bur block serve to determine the length for  
 example of gutta-percha points.  
 - The foam insert has to be discarded after one use.

4. Reprocessing  
 Those medical devices that are supplied non-sterile have to be prepared prior to first use.  
 The bur block can be disinfected with suitable cleaning and disinfecting agents (e.g. with Komet DC1). The  
 recommendations for use (immersion time, concentration, suitability) of cleaning and disinfecting agents can be

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seen in the information provided by the manufacturers of these agents. Cleaning and disinfecting agents must be  
 rinsed off thoroughly with water and the bur block must be dried carefully (preferably with compressed air  
 according to the recommendations of the Commission for Hospital Hygiene and Infectious Disease Prevention of  
 the Robert Koch Institute). The bur block can be reprocessed in the thermo disinfectant, provided that a suitable  
 agent is used. Make sure that the bur block is placed in the washer/disinfection with the opening facing  
 downwards in order to ensure that all areas of the bur block are reached by the spray jet.

Sterilization:  
 Sterilization is carried out in the autoclave at 134°C. Observe the instructions provided by the manufacturer. For  
 sterilization in the autoclave, the bur block has to be sealed in sterilization foil. After sterilization, the bur block  
 remains sealed in the foil until the next use. After sterilization, allow the bur block to cool down properly to avoid  
 excessive condensation. The single-use foam insert can be treated in the autoclave prior to use. To this end, the  
 foam insert must be sealed in sterilization foil and stored inside the foil until use. After use, the foam insert has  
 to be discarded because the safety of this product in case of reuse can no longer be guaranteed. The operator of  
 medical products is responsible for ensuring that proper treatment is carried out by qualified personnel, using  
 suitable materials and suited equipment, as recommended by the Commission for Hospital Hygiene and  
 Infectious Disease Prevention of the Robert Koch Institute. Information with regard to proper reprocessing  
 according to DIN EN ISO 17664, the present instructions for use and more detailed information, as amended,  
 can be downloaded from our web site www.kometdental.de or requested from the manufacturer Gebr. Brasseler.

5. Safety and responsibility  
 See also "Instructions for use and safety recommendations" in the current Komet catalogue. The user is  
 responsible for checking the products prior to use as to whether they are suited for the intended purpose. The  
 user is responsible for the appropriate use of the instruments. In case of contributory negligence by the user,  
 Gebr. Brasseler partially or totally declines liability for all resulting damages, particularly if these are due to  
 non-observance of our recommendations for use or warnings as well as inadvertent misuse by the user.

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