# F360 Fill | Gebrauchsanweisung

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#### F360 Fill oven

##### Classification

The F360 Fill oven is a medical device, class I, in compliance with MDD (Medical Devices Directive) section IX. The oven must only be used in dental clinics by dental personnel, and only for heating F360 Fill obturators. The product is classified according to DS EN 61010-1- 2001, therefore, it must not be used within reach of patients.

The F360 Fill oven has a heat chamber that consists of 5 printed circuit boards with copper lanes. By applying power to the copper lanes heat is generated. The oven is driven by an external power supply.

The F360 Fill oven is covered by a 2 year warranty.

The F360 Fill oven must only be serviced and repaired by the manufacturer otherwise the warranty will be void.

##### Warning

Do not touch the inside of the heat chamber before the oven has completely cooled off.

##### Technical specifications

Power supply

Input power: 100-240V AC, 50/60Hz

Output power: 9V DC, 1,7A

##### Safety

Auto turn off: The oven turns off automatically after 9 hours

##### Heating times

Obturator heating time is 60 seconds.

Preheating of the oven takes 3-4 minutes depending of the surrounding temperature.

#### 1. Start using the oven

1. Connect the power supply to the DC inlet of the oven and the mains and switch on power (A1)

2. Indicator A emits a red light while the oven preheats (A2). When the oven is ready, the indicator emits a green light, and you will hear a beep. NOTICE: The oven can be left on the entire day.

#### 2. Use

When the oven is ready for use, up to four obturators at a time may be placed in the appropriate slots.

Place the obturators in the slots marked 1 to 4 (A3).

TIP: Place the obturator to be used first in slot #1, the second in #2 etc. Generally you should begin with the smallest size obturator.

1. Activate the timer by pressing the ON icon (A4).

2. Indicator B lights red during the heating cycle.

3. When ready, the indicator changes to green and you will hear a beep (A4).

Let the obturator stay in the oven until you are ready to insert it (max. 15 min.).

4. It is possible to turn off the oven by holding down the ON button for 2 seconds (A4). The oven beeps three times, it indicates it is switched off. Hold down the ON button for 2 seconds to reactivate the oven.

Tip: You do not have to turn off the oven after finished treatment, the heater can be left on standby all day.

##### Indication - light and sound

	Time	Heat	Sound
Preheating			-
Ready			Long beep
Activation			Short beep
Obturator ready			Long beep
Turn off			3 x short beep
Power on			Long beep
Auto turn off	-	-	Long beep

#### 3. Cleaning

When the oven has cooled completely, gutta-percha residue may easily be removed.

The top of the oven may be removed for cleaning by rotating the top counter-clockwise 1 until the top is loose and can be removed 2 (A5).

The top may be washed in a dishwasher and/or autoclave.

Gutta-percha residue may be removed from the heat chamber by the

use of a blunt plastic instrument.

The main part of the oven may be wiped with a slightly damp cloth, using standard cleaning solutions and surface disinfection products.

#### Sterilization guide for the removable top part

WARNINGS	Discharge damaged products
Limitations on reprocessing	250 cycles
INSTRUCTIONS	
Point of use	Remove excess soil with disposable cloth/paper tissue
Preparation for decontamination	No specific requirements
Cleaning: Automated	Not recommended
Cleaning: Manual	Equipment: Detergent, brush, distilled or demineralised water 1. Rinse excess soil from cover 2. Using brush, apply detergent solution to all surfaces 3. Rinse under clear running water
Disinfection	Not recommended
Drying	Do not exceed 134 <sup>o</sup> C
Maintenance, inspection and testing	Discharge damaged products
Packaging	Only pack individually
Sterilization	Vacuum autoclave, minimum 18 minutes at 134 <sup>o</sup> C, pressure 2.2 bar Do not exceed 134 <sup>o</sup> C
Storage	No particular requirements

#### 4. Labelling

Consult instructions

Serial number

Medical Electrical Equipment Directive, DS/EN 60601-1: 2006 Electric shock protection, Type B

Medical device, Directive MDD  
MDD Directive 93/42/EEC  
Medical device, class I

Waste Electrical and Electronic Equipment Directive,  
WEEE Directive 2002/96/EC

Please return the product for recycling.

#### F360 Fill obturators

##### Classification

The F360 Fill obturator is a medical device; class IIA, in compliance with MDD (Medical Devices Directive) section IX. The obturator must only be used by dentists. Read and follow the instructions carefully!

##### General precautions

- The product should only be used by dentists.
- The product should only be used for root canal procedures in teeth.
- A deformed or defective product should never be used, but should be discarded.
- It is the responsibility of the user to check the products prior to use to ensure that they are suited for the intended purpose. The operator is responsible for the correct use of the instruments. In case of contributory negligence by the user, Gebr. Brasseler partially or totally declines liability for all resulting damages, particularly due to the non-observance of our recommendations for use or warnings as well as inadvertent misuse by the user. For information on the correct use of the product, please see the step-by-step guide, section 2.
- The product should be stored in the original packaging.

##### Warnings - obturator

- After heating, the gutta-percha is hot. Avoid any direct contact between the obturator and skin or mucous membranes. The obturator should be inserted directly into the root canal.
- It is always recommended to take a control x-ray. If the x-ray is not satisfactory, take another x-ray from a different angle. If the filling still cannot be satisfactorily documented, remove the obturator.

#### 1. Description of the F360 Fill obturator

A F360 Fill obturator is a carrier of biocompatible plastic, which is coated with thermoplastic gutta-percha. The sizes are 20 to 60 according to ISO standard. The obturator sizes correspond to the ISO standards of root canal files. A single obturator is all that is needed to fully obturate a root canal. The result should always be a perfectly placed root canal filling with a tight apical seal.

A F360 Fill obturator consists of three parts:

##### 1. Plastic handle with a metal insertion pin

The plastic handle is color-coded according to the size of the obturator. The stainless steel insertion pin is 9 mm long and permanently attached to the plastic handle. The metal insertion pin is placed 6 mm down in the plastic carrier. When the metal insertion pin is removed it will leave a hollow part of 6 mm in the plastic carrier.

##### 2. Outer layer of thermoplastic Natural GP™ gutta-percha

The gutta-percha is reversible thermoplastic. This means that it becomes soft and highly adhesive when heated to temperatures above approx. 100°C. When cooled, it returns to its original firm, rigid state. Another heating will once again make the gutta-percha soft and adhesive.

##### 3. Plastic carrier with endo stop

The tapered plastic carrier is sufficiently flexible to negotiate curved canals.

The F360 Fill plastic carrier is slightly tapered. This insures adequate back flow of the heated gutta-percha. The length of the carrier is 24 mm. The 6 mm coronal portion of the carrier is hollow to accommodate the metal insertion pin. A rubber endo stop is provided to indicate the working distance measurement. The plastic carrier and gutta-percha are radiopaque.

#### 2. A step-by-step guide

Before using the F360 Fill obturator, assure that you have achieved the following:

- A properly cleaned and shaped canal – the key to any successful root canal filling. Utilize any technique which provides a smooth tapered form and coronal access wide enough to give room for the instruments.
- A dry root canal.

- An established working length. The plastic carrier is flexible enough to negotiate curved canals. It is highly advised that rotary files with a higher conicity are utilized to create a sufficient taper for the obturator. When using files with a conicity of 2 %, it is advised to combine them with a procedure using Gates Glidden Drills or similar files which can be used to create a sufficient opening in the coronal part of the canal. The use of size verifiers is recommended.

##### STEP 1: Determine the working length and choose the proper obturator size

The working length is measured all the way to the apex (B1). We recommend measuring the working length from the edge of the cusp. The working length may be determined by use of x-ray, file and electronic apex locator. Generally, the appropriate size F360 Fill obturator will be the same size as the last file used at the apex of the canal. TIP: For very narrow and/or highly calcified canals, it may be useful to select one size smaller than the last file used at the apex. If you do not use rotary files with a taper of 4% or more, you should always verify the size of the cleaned root canal with a size verifier before inserting the obturator. Mark your working length on your F360 Fill obturator. When the size verifier has a slightly loose fit in the apical third of the root canal, the rubber stop is placed at the working length. Afterwards, the working length is transferred to the obturator.

##### STEP 2: Heating of the obturator

Place the selected obturator in the F360 Fill oven and activate it (B2, B3). See above for more detailed information on the appropriate use of the oven.

##### STEP 3: Application of sealer

While the obturator is heating, mix any heat resistant sealer (usually a non-eugenol sealer) with long working time, e.g. Komet EasySeal, 9978. Make sure that the root canal is completely dry before applying sealer.

Place a very thin coating of sealer on the wall of the canal. Too much sealer is neither desirable nor necessary. Use a paper point or a F360 Fill Size Verifier to apply a thin layer of sealer on the canal walls. When obturating multiple canals, apply sealer to all canals at once.

#### STEP 4: Insertion of the obturator into the canal

When the oven indicates that the obturator is ready for use, the obturator is removed from the oven. See above for more detailed information on the appropriate use of the oven. Without twisting the handle, immediately insert the obturator into the canal to the working distance, using a firm and steady pressure (B4).

#### STEP 5: Confirming radiograph

Confirm your root canal treatment by taking radiographs.

#### STEP 6: Removal of excess central carrier and gutta-percha

When the gutta-percha has hardened (after 2-3 min.), remove the handle in one of two ways (B5).

1. Lean the handle and metal insertion pin sharply to one side to break off the excess.
2. Stabilize the plastic carrier with a finger and cut away the handle and excess plastic carrier with a small bur, e.g. the Komet gutta-percha cutter GP801L.314.014.

Be careful to prevent vertical movement of the carrier to minimize the chance of disturbing the apical seal. Trim away the extra gutta-percha using a gutta-percha cutter e.g. the Komet GP801L.314.014.

### 3. Technical tips

It is very easy to use the F360 Fill system. Feed-back from experienced users, however, will hopefully make it even easier for you to get started.

#### Obturating multi-rooted teeth

There will almost always be more than enough gutta-percha on each obturator. When obturating the canal, the excess gutta-percha and sealer will backfill to the coronal part of the tooth, possibly blocking the adjacent opening. However, you can avoid this:

- Obturate the shortest canal first.
- Place paper points or size verifiers in the other canals prior to obturation to prevent the following root canals in being blocked by gutta-percha. Remove paper points or size verifiers as each canal is obturated.
- While there is sufficient gutta-percha on each obturator to fill even extremely long canals, wide canals or internally resorbed canals, it may be too much in some cases. If it is obvious that there will be an excess of gutta-percha to fill the canal, use a sharp blade to trim some of the material from the coronal end of the central carrier prior to heating.
- In some cases it can be beneficial to remove the handle/metal insertion pin from the actual plastic carrier, especially in multi-rooted teeth where the extra space will make it easier placing the following obturators in the second and third canal.
- Prior to heating of the obturator(s) and prior to marking the working length, always assure that the handle/metal insertion pin can be separated from the plastic carrier. This is done by twisting the handle/metal insertion pin out of the plastic carrier. Re-insert the handle/metal insertion pin into the plastic carrier.

Mark your working length and proceed to Step 2 - Step 4 (section 2). Stabilize the plastic carrier with a cotton plier while the handle/metal insertion pin is twisted out. Any excess material is removed as described in Step 6 (section 2).

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### Use of sealer

As with any other root filling technique, when using the F360 Fill system, always apply sealer in the canal first. However, there is a significant difference when using F360 Fill. Only apply a very small amount of sealer using a size verifier or a paper point. A root spiral (lentulo) should not be utilized when using obturators as there is a tendency to add too much sealer to the canal. Remember; only apply a very thin layer of sealer on the canal walls.

### 4. Creation of space for the post





Post space creation is accomplished by removing the coronal portion of the obturation. Since the coronal 6 mm of the plastic carrier is hollow, this procedure should be easier with the F360 Fill obturator than with other obturation systems.

### 5. Retreatment of an obturated root canal

Removal of the obturator is relatively simple prior to the gutta-percha cooling. Simply grasp the excess carrier with cotton pliers and remove it. Once the gutta-percha has cooled or for retreatment at a later appointment, removal is easiest achieved by using a Peeso bur, Gates or similar:

1. Place the bur between the canal wall and the plastic carrier, so that the rotational direction will be counter clockwise. Use the instrument at very low speed (250 - 400 rpm) and work slowly around the plastic carrier until you feel resistance. The instrument will grasp between the plastic carrier and the canal wall. The plastic carrier will be coiled out of the canal after a short time.
2. If necessary, a small Hedstrom hand file can be worked down apically along the plastic carrier.
3. When the plastic carrier has been pulled out of the canal, the gutta-percha is removed by traditional means.

### 6. Labelling

-  Further information available
-  For single use only
-  Use by
-  Lot number
-  CE Class I F360 Fill Verifier
-  CE Class IIa F360 Fill Obturator

### 7. Storage

Observe the use-by date.

### 8. Safety and liability

It is the responsibility of the user to check the products prior to use to ensure that they are suited for the intended purpose. The operator is responsible for the correct use of the instruments. In case of contributory negligence by the user, Gebr. Brasseler partially or totally declines liability for all resulting damages, particularly due to the non-observance of our recommendations for use or warnings as well as inadvertent misuse by the user.

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Gebr. Brasseler GmbH & Co. KG
CMS Dental A/S