
 <p>② <b>97502L15</b></p> <p>Extrusionsstift Ø 1,2 mm        Extrusion pin Ø 1,2 mm        Tenon d'extrusion Ø 1,2 mm        Pin de extrusión Ø 1,2 mm        Perno di estrusione Ø 1,2 mm</p> <p>Glasfaserverstärktes Composite        Glass fibre reinforced composite        Composite renforcé avec fibres de verre        Composite reforzado por fibras de vidrio        Composito rinforzato in fibra di vetro</p> <p>(Handlingsteil, TPE Thermoplastisches Elastomer)        (Handling sleeve, TPE Thermoplastic elastomer)        (Tuyau de manipulation, TPE Elastomère thermoplastique)        (Manguito, TPE Elastómero termoplástico)        (Sezione di manipolazione, TPE Elastomero termoplastico)</p>	 <p>② <b>97503</b></p> <p>Linse für Extrusionsstift        Lens for extrusion pin        Élément lenticulaire pour le tenon d'extrusion        Elemento lenticular para el pin de extrusión        Lente per perno di estrusione</p> <p>Polyetheretherketon (PEEK)        Polyetheretherketone (PEEK)        Polyéther éther cétone (PEEK)        Poliéter éter cetona(PEEK)        Polietere etere chetone (PEEK)</p> <p>(Handlingsteil, TPE Thermoplastisches Elastomer)        (Handling sleeve, TPE Thermoplastic elastomer)        (Tuyau de manipulation, TPE Elastomère thermoplastique)        (Manguito, TPE Elastómero termoplástico)        (Sezione di manipolazione, TPE Elastomero termoplastico)</p>	 <p>② <b>97505L25</b></p> <p>Spannstange (Steg)        Bar        Barre        Barra tensora        Barra di serraggio (anima)</p> <p>Glasfaserverstärktes Composite        Glass fibre reinforced composite        Composite renforcé avec fibres de verre        Composite reforzado por fibras de vidrio        Composito rinforzato in fibra di vetro</p>	 <p>② <b>97500.032</b> ② <b>97500.064</b> ② <b>97500.095</b> medium        ② <b>97501.032</b> ② <b>97501.064</b> ② <b>97501.095</b> strong</p> <p>Spanngummi        Dental elastic        Elastique de tension        Goma tensora        Gommino di serraggio</p> <p>Silikon, latexhaltig        Silicone, with latex        Silicone, contient du latex        Silicona, contiene látex        Silikone, con latex</p> 
 <p>○<sub>rot.</sub> 160.000 min<sup>-1</sup>        ○<sub>max.</sub> 300.000 min<sup>-1</sup></p> <p><b>S6881.314.012</b></p>	 <p>○<sub>rot.</sub> 20.000 min<sup>-1</sup>        ○<sub>max.</sub> 450.000 min<sup>-1</sup></p> <p><b>8390.314.016</b></p>	 <p>○<sub>rot.</sub> 160.000 min<sup>-1</sup>        ○<sub>max.</sub> 160.000 min<sup>-1</sup></p> <p><b>953M.314.014</b></p>	 <p>○<sub>rot.</sub> 160.000 min<sup>-1</sup>        ○<sub>max.</sub> 160.000 min<sup>-1</sup></p> <p><b>953AM.314.014</b></p>


Abb./Fig. 1                      Abb./Fig. 2                      Abb./Fig. 3                      Abb./Fig. 4                      Abb./Fig. 5                      Abb./Fig. 6

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**EN**  
**TMC Extrusion with a transversal extrusion pin**

The procedure according to the TMC philosophy (Tissue Master Concept, Dr. Stefan Neumeyer) is largely distinguished by a substantial reduction of the treatment time. To this end, the forces applied significantly exceed those applied during traditional orthodontic treatments. Once the extrusion has been completed, the transversal element (extrusion pin) serves to safely fix the extruded tooth or tooth segment.

**Definition of the term "extrusion"**  
 The terms "orthodontic extrusion" or "forced eruption" describe a load-induced movement of a tooth in direction of its longitudinal axis, out of its alveolus. The movement is orthodontically induced and does not lead to the extraction of the tooth and the consequential separation of tissue structures but to a new coronal position of the tooth thanks to the subsequent, stabilizing reconstruction and formation of tissue.

- Indication**  
 Orthodontic extrusion for the conservation and regeneration of periodontal tissue structures  
 - Extrusion of whole teeth  
 - Extrusion of periodontically diseased teeth  
 - Extrusion of reimplanted resected teeth  
 - Extrusion reimplanted root segments
- Contra-indications**  
 The general and special contra-indications applicable in  
 - orthodontics  
 - dental surgery  
 - oral and maxillofacial surgery  
 have to be observed.
- Clinical sequence**  
 The course of the treatment is always structured to suit the respective case and indication. The structure, i.e. the order of the individual treatment steps is optimally adapted to ensure an effective procedure achieving the aspired result. The treatment steps themselves are always carried out in the same way. The transversal extrusion pin is used as follows:  
 3.1 Positioning of the extrusion pin  
 Provided that enough residual tooth substance is available in the coronal region, a transversal perforation can be drilled in labio/lingual direction with the diamond instrument S6881 at  $\omega_{rot.} 160,000 \text{ rpm}$  in the red contra-angle. Work in intermittent mode with plenty of cooling (at least 50 ml/min). Next, the extrusion pin is roughly shortened to the desired length outside the mouth with spray cooling. The shortened extrusion pin is then inserted into the previously drilled perforation from the lingual side. To this end, the extrusion pin can be placed into the handling sleeve which is then held with tweezers. To ensure that the transversal extrusion pin is securely fixed, a small quantity of light-curing or dual-hardening composite can be applied to the centre part of the pin (e.g. RelyX Unicem). Once the pin is in its final position, the construction is light cured and the handling sleeve can be removed with ease. A little composite is then applied onto the lens inside the handling sleeve. Following that, the handling sleeve with the lenticular element inside is gripped with tweezers and slipped over the labial end of the extrusion pin. The joint is then light cured with a dental halogen lamp for 20 seconds, after which the silicon handling sleeve is pulled off the lens. In the event that the end of the extrusion pin protrudes through the lens, remove any excess with a fine rotary diamond instrument, for example 8390 (at a speed of  $\omega_{rot.} 20,000 \text{ rpm}$  in the red contra-angle) and finish the surface. In case of insufficient residual dental substance or in the absence of a clinical dental crown (fig. 1), the coronal area of the tooth to be extruded can be smoothed with a rotary diamond instrument and provided with a groove (fig. 2). This can be done with the diamond instrument S6881 in the red contra-angle at a speed of  $\omega_{rot.} 160,000 \text{ rpm}$ . Following that, the dental substance is bonded and light cured for approx. 20 seconds. The extrusion pin is then fixed on the conditioned tooth stump with composite (fig. 3), after which the area is light cured with UV light for approx. 20 seconds. Self-etching or dual-hardening composite systems are suitable for this purpose.  
 3.2 Creation of a counter bearing  
 To create a counter bearing, the tie rod (bar) can for example be fixed to the neighbouring teeth with a little composite (fig. 4). To this end, shorten the tie rod by means of a rotary diamond instrument (see 3.1) and fix to the bonded surface of the adjacent teeth with a small quantity of DentinBuild EVO.

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3.3 Once the bar is in place, the dental elastic – which comes in different lengths and thicknesses, for example 97500 (medium) or 97501 (strong) – is positioned at both ends of the extrusion pin (fig. 5). Due to the dental elastic placed across the bar, the extrusion pin, and consequently the root of the tooth, is subjected to tension loads (approx. 20-30 N). To ensure sufficient tension, the dental elastic should be exchanged every day. The dental elastics come in 2 thicknesses and 3 different sizes.

As an alternative, a temporary appliance prepared in the laboratory or the clinical crown – bonded to the adjacent teeth – can also serve as a counter bearing (see fig. 6). The gap can be filled, for example with white orthodontic wax, in order to achieve an aesthetically pleasing look. In case of a particularly extensive extrusion, the fixation of the extrusion pin has to be repeated in a lower position. The procedure is identical to that of the initial fixation.

- Further treatment**  
 After 20-30 days, another appointment should be made to check the function of the dental elastic. If the elastic has become too slack to exert sufficient tension loads, it should be replaced by the next thicker or smaller one. Once the root has moved to the desired extent (several mm), the tooth can be extracted. A grade 3 loosening should be achieved at that point. This can be followed by the further reconstruction, for example the immediate insertion of an implant.
- Maintenance/care**  
 All component parts of the system are delivered non-sterile. The extrusion pin, the lens and the tie rod are made of composite and may only be disinfected with medical alcohol – they must not be sterilised! Place rotary instruments in a cleaning/disinfection tank (Fräsator) filled with a suitable detergent/disinfectant (e.g. Komet DC1) immediately after use in the mouth to prevent drying of residues on the instruments (protein fixation). It is recommended to reprocess the instruments within one hour of use at the very latest. The instruments should be in the cleaning/disinfection tank (Fräsator) when transported to the site where the reprocessing is to take place. The instruments designed for extrusion are disposable. Only the rotary instruments of the system can be reused and are suitable for reprocessing. Follow the instructions enclosed with the instruments.

**Universally valid notes:**  
 Observe the legal provisions regarding the reprocessing of medical products valid in your country (e.g. www.rki.de). The manufacturer confirms that the above detailed reprocessing methods are suitable for preparing the above named instrument group to enable their reuse. The operator of medical products is responsible for seeing that proper treatment is carried out by qualified personnel, using the appropriate materials and suited equipment, as recommended by the Commission for Hospital Hygiene and Infectious Disease Prevention of the Robert Koch Institute. Work instructions with regard to proper reprocessing of instruments according to DIN EN ISO 17664, the present instructions for use and more detailed information, as amended, can be downloaded from our web site www.kometdental.de or requested from the manufacturer Gebr. Brasseler.

- Storage**  
 Store packed instruments away from UV radiation and heat in a clean and dry environment. Do not store in the same room as solvents or chemicals.
- Disposal**  
 The instruments have to be disposed of in tight, non-breakable and non-pierceable containers (protection from contamination).
- Safety and liability**  
 See also "Instructions for use and safety recommendations" in the current Komet Dental catalogue. The user is responsible for checking the products prior to use as to whether they are suited for the intended purpose. The user is responsible for the appropriate use of the instruments. In case of contributory negligence by the user, Gebr. Brasseler partially or totally declines liability for all resulting damages, particularly if these are due to nonobservance of our recommendations for use or warnings as well as inadvertent misuse by the user.
- Procedure in case of serious adverse events**  
 In case of serious adverse events, please proceed in line with the legal regulations in your country.

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