


Hubfeilen Skalpell · Files (scalpels) · Limas (bisturis) · Limes (bistouris) · Lime a scalpello			
 <p>RCAS1C = Initial</p>	 <p>RCAS2/RCAS2C = Advanced</p>	 <p>RCAS3C = Extended</p>	Rostfreier Stahl Stainless Steel Acier inoxydable Acero inoxidable Acciaio inossidabile
Applikationshilfe/Halter für Hubfeilen · Applicator/support · Aplicador/Soporte para limas · Applicateur/support pour limes · Ausilio di applicazione/supporto per lime			
			LCP
97508			

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Files (scalpels) for work on soft tissue

1. Indication
 Scalpels for work on soft tissue during surgical interventions in the oral region, intended for use in a reciprocal contra-angle.

2. Contra-indication
 The instruments must not be used in the proximity of nerves.

3. Clinical sequence
 The scalpels are suitable for cutting soft tissue, for example during an advancement flap reconstruction. The scalpel is firmly inserted into the reciprocal contra-angle by means of the applicator 97508. The applicator is then pulled off the working part of the instrument (fig. 1, 2). All reciprocal contra-angles with a stroke between 0.4 and 1.1 mm and a frequency of 5,000 to 10,000 Hz are suitable, e.g. W&H Synea - Profin® WA-67. The instrument is cooled with an externally supplied sterile isotonic saline solution. The cooling liquid can be supplied by means of a filled syringe (e.g. 20 ml) with a tip/needle sized 0.9x42 mm. The instrument must be guided carefully and precisely in order to make an exact incision. The soft issue is cut according to the requirements of the preparation, by means of vertical, horizontal or angled stroke movements performed by the blade and simultaneous lengthwise movements of the reciprocal contra-angle parallel to the cutting edge. Avoid jamming of the instrument and do not apply excessive contact pressure as this can lead to instrument fracture and consequential injury to the patient. After use, the applicator can be placed onto the working part of the instrument which can then be removed from the power system with protected blades. The instrument is reprocessed with mounted applicator. Attention: An increased stroke/amplitude reduces the sensitivity of the instrument. Avoid excessive stroke frequencies and inappropriate movements – such as allowing the instrument to jam – and do not work without cooling. These factors can lead to the fracture of the instrument and consequently to an increased risk of injury to the patient.

4. Maintenance/care
 All system components are delivered non-sterile and therefore have to be processed prior to first use. Place instruments in a cleaning/disinfection tank (Fräsator) filled with a suitable detergent/disinfectant (e.g. Komet DC1) immediately after use in the mouth to prevent drying of residues on the instruments (protein fixation). It is recommended to reprocess the instruments within one hour of use at the very latest. The instruments should be in the cleaning/disinfection tank when transported to the site where the reprocessing is to take place.

Reprocessing
 Manual pre-cleaning (prior to mechanical and manual reprocessing):
 1.Remove instruments from the cleaning /disinfection tank immediately before reprocessing. Any stubborn contamination should be removed with a nylon brush under running water. To this end, remove the reciprocal file from the applicator/support 97508. After manual pre-cleaning, reinsert the reciprocal instrument into the applicator/support.

Mechanical reprocessing (validated):
 1. Rinse reciprocal instrument inside the applicator thoroughly under running water to prevent any residues of the detergent/disinfectant from getting into the washer/disinfectator.
 2. Unscrew the nozzle on the injector bar of the washer/disinfectator. Then press the applicator with the inserted reciprocal instrument into the silicon holder, using a little force, until it locks. The correct position has been reached when the cylindrical part of the applicator is completely inserted in the silicon holder.
 3. Put chemical detergent into the washer/disinfectator, following the indications on the label and the instructions of the manufacturer of the washer/disinfectator.
 4. Start the Vario TD programme including thermal disinfection. Thermal disinfection takes place allowing for the A0 value and observing national provisions (prEN/ISO 15883).
 5. On completion of the cycle remove instruments from the washer/disinfectator and dry (preferably with compressed air as recommended by the Robert Koch Institute). When drying the inside of the hollow applicator, please make sure to hold the compressed air gun close to the opening, in order to guarantee that sufficient air flows through the instrument.
 6. Visual examination to ensure that the instrument is clean and undamaged. If after mechanical reprocessing there are still visible residues of contamination, repeat the cleaning and disinfecting process until no visible contamination is left.

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Instructions for use
 Please read carefully and retain for future reference.

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5. Sicherheit und Haftung
 Make sure that the packaging is suitable for instruments classed as risk group critical B and the chosen method of sterilisation. Single pack: The packaging must be large enough to ensure that there is no pressure on the seal. Seal instruments in a suitable sterilisation container. For this procedure, the following instructions have to be observed: Steam sterilisation using a vacuum process at 134°C in a device that complies with the provisions of DIN EN 13060; with validated processes.
 - Fractionated pre-vacuum (type B)
 - Sterilisation temperature: 134°C
 - Hold time: at least 5 minutes (full cycle)
 - Drying time: at least 10 minutes
 In order to prevent staining and corrosion, the steam must be free of particles. The recommended limits for particle contents in feed water and condensed steam are defined by standard DIN EN 13060. Make sure not to exceed the maximum capacity of the sterilizer when sterilizing several instruments. Follow the instructions of the device manufacturer.

Manual reprocessing (alternative):
 1. Rinse instrument with detergent/disinfectant. In case of residual contamination, rinse again with detergent/disinfectant.
 2. Remove stubborn contamination with a nylon brush under running water. To this end, remove the reciprocal file from the applicator/support 97508.
 3. Rinse instrument thoroughly with running water.
 4. Visual examination to ensure that the instrument is clean. If there are still visible residues of contamination, repeat the cleaning process until no visible contamination is left.
 5. Place instruments into the ultrasonic device or instrument bath filled with detergent/disinfectant.
 6. During chemical disinfection in the ultrasonic device or instrument bath, observe the instructions of the manufacturer regarding concentration and immersion time. Be sure to observe the full correct immersion time which does not start until the last instrument has been immersed. Attention: do not exceed 45°C (risk of protein coagulation)!
 7. On completion of the immersion time, rinse instrument thoroughly with suitable water (preferably with demineralised water to avoid residues of lime).
 8. Dry instruments (preferably with compressed air as recommended by the Robert Koch Institute). When drying the inside of the hollow applicator, please make sure to hold the compressed air gun close to the opening, in order to guarantee that sufficient air flows through the instrument.
 9. Visual examination to ensure that the instrument is clean and undamaged. If there are still visible residues of contamination, repeat the cleaning and chemical disinfecting process until no visible contamination is left. According to the recommendations of the Robert Koch Institute, the further reprocessing should be carried out mechanically. Information on the reprocessing of instruments in compliance with the standard EN ISO 17664 can be downloaded from our website www.kometdental.de or requested directly from the manufacturer Gebr. Brasseler.

Sterilisation in the autoclave
 Make sure that the packaging is suitable for instruments classed as risk group critical B and the chosen method of sterilisation. Single pack: The packaging must be large enough to ensure that there is no pressure on the seal. Seal instruments in a suitable sterilisation container. For this procedure, the following instructions have to be observed: Steam sterilisation using a vacuum process at 134°C in a device that complies with the provisions of DIN EN 13060; with validated processes.
 - Fractionated pre-vacuum (type B)
 - Sterilisation temperature: 134°C
 - Hold time: at least 5 minutes (full cycle)
 - Drying time: at least 10 minutes
 In order to prevent staining and corrosion, the steam must be free of particles. The recommended limits for particle contents in feed water and condensed steam are defined by standard DIN EN 13060. Make sure not to exceed the maximum capacity of the sterilizer when sterilizing several instruments. Follow the instructions of the device manufacturer.
 After sterilisation, check instruments for corrosion. Do not reuse corroded instruments.

Universally valid notes
 Observe the legal provisions regarding the reprocessing of medical products valid in your country (e.g. www.rki.de). The manufacturer confirms that the above detailed reprocessing methods are suitable for preparing the above named instrument group to enable their reuse. The person carrying out the reprocessing is responsible for ensuring that the applied method is carried out with appropriate equipment, materials and trained personnel at the reprocessing site and that it actually achieves the desired result. To guarantee this, routine controls of the validated mechanical preparation methods are normally necessary. Any deviation from the above detailed process must be carefully checked by the person carrying out the reprocessing to ensure effectiveness and to avoid possible adverse consequences.

5. Safety and liability
 Worn and damaged instruments (blunt or bent instruments etc.) are to be discarded and immediately replaced by new ones. The above instructions regarding the handling of the instruments, cooling and contact pressure are to be strictly observed. Only use the instruments according to their intended purpose. The non-observance of these safety instructions can lead to injuries or damage to the power system. It is the responsibility of the user to check the products prior to use to ensure that they are suited for the intended purpose. In case of contributory negligence by the user, Gebr. Brasseler partially or totally declines liability for all resulting damages, particularly due to non-observance of our recommendations for use or warnings as well as inadvertent misuse by the user.

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