


**PA-Feilen, diamantiert · Periodontal files (diamond coated) · Limes parodontales, diamantées · Limas periodontales, diamantadas · Lime PA, diamantate**

 <p><b>RCAP1 = Surfer</b></p>	 <p><b>RCAP2 = Smoother</b></p>	 <p><b>RCAP3 = Beaver</b></p>	<p>           Rostfreier Stahl            Stainless Steel            Acier inoxydable            Acero inoxidable            Acciaio inossidabile         </p>
---	--	--	--

**Applikationshilfe/Halter für Hubfeilen · Applicator/support for files · Appicateur/support pour limes · Aplicador/Soporte para limas · Ausilio di applicazione/supporto per lime**

 <p><b>97508</b></p>			<p>LCP</p>
--	--	--	------------

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**EN**

Diamond coated files of optimal shape for periodontal treatments

1. Indication  
Mechanical surface treatments and smoothing the surfaces of teeth and roots as part of a periodontal treatment.

2. Contra-indication  
Undermining periodontal defects, especially on curved root surfaces.

3. Clinical sequence  
The periodontal instruments are suitable for removing deposits such as plaque, tartar and concretions, from the surfaces of dental roots, both below and above the gingiva. The instrument is firmly inserted into the reciprocal contra-angle by means of the applicator 97508. The applicator is then pulled off the working part of the instrument (fig. 1, 2). Reciprocal contra-angles with a stroke between 0.4 and 1.1 mm and a frequency of 5,000 to 10,000 Hz are suitable, e.g. W&H Synea - Profin® WA-67. The instrument is guided along the neck of the tooth with careful, controlled movements, with low pressure and supplying constant cooling. The cooling liquid can be supplied by means of a filled syringe (e.g. 20 ml) with a tip/needle sized 0.9x42 mm. After use, the applicator can be placed back onto the working part so that the instrument can be removed from the power system with protected blades. The instrument is reprocessed with mounted applicator. Attention: Avoid excessive stroke frequencies and inappropriate movements – such as allowing the instrument to jam – and do not work without cooling. These factors can lead to the fracture of the instrument and consequently to an increased risk of injury to the patient.

4. Maintenance/care  
All system components are delivered non-sterile and therefore have to be processed prior to first use. Place instruments in a cleaning/disinfection tank (Fräsator) filled with a suitable detergent/disinfectant (e.g. Komet DC1) immediately after use in the mouth to prevent drying of residues on the instruments (protein fixation). It is recommended to reprocess the instruments within one hour of use at the very latest. The instruments should be in the cleaning/disinfection tank when transported to the site where the reprocessing is to take place.

Reprocessing  
Manual pre-cleaning (prior to mechanical and manual reprocessing):  
Remove instruments from the cleaning /disinfection tank immediately before reprocessing. Any stubborn contamination should be removed with a nylon brush under running water. To this end, remove the reciprocal file from the applicator/support 97508. After manual pre-cleaning, reinsert the instrument into the applicator/ support.

Mechanical reprocessing (validated):  
1. Rinse the reciprocal instrument inside the applicator thoroughly under running water to prevent any residues of the detergent/disinfectant from getting into the washer/disinfectant.  
2. Unscrew the nozzle on the injector bar of the washer/disinfectant. Then press the applicator with the inserted reciprocal instrument into the silicon holder, using a little force, until it locks. The correct position has been reached when the cylindrical part of the applicator is completely inserted in the silicon holder.  
3. Put chemical detergent into the washer/disinfectant, following the indications on the label and the instructions of the manufacturer of the washer/disinfectant.  
4. Start the Vario TD programme including thermal disinfection. Thermal disinfection takes place allowing for the A0 value and observing national provisions (prEN/ISO 15883).  
5. On completion of the cycle remove instruments from the washer/disinfectant and dry (preferably with compressed air as recommended by the Robert Koch Institute). When drying the inside of the hollow applicator, please make sure to hold the compressed air gun close to the opening, in order to guarantee that sufficient air flows through the instrument.  
6. Visual examination to ensure that the instrument is clean and undamaged. If after mechanical reprocessing there are still visible residues of contamination, repeat the cleaning and disinfecting process until no visible contamination is left.

Manual reprocessing (alternative):  
1. Rinse instrument with detergent/disinfectant. In case of residual contamination, rinse again with detergent/disinfectant.  
2. Remove stubborn contamination with a nylon brush under running water. To this end, remove the reciprocal file

from the applicator/support 97508.

3. Rinse instrument thoroughly with running water.  
4. Visual examination to ensure that the instrument is clean. If there are still visible residues of contamination, repeat the cleaning process until no visible contamination is left.  
5. Place instruments into the ultrasonic device or instrument bath filled with detergent/disinfectant.  
6. During chemical disinfection in the ultrasonic device or instrument bath, observe the instructions of the manufacturer regarding concentration and immersion time. Be sure to observe the full correct immersion time which does not start until the last instrument has been immersed. Attention: do not exceed 45°C (risk of protein coagulation)!  
7. On completion of the immersion time, rinse instrument thoroughly with suitable water (preferably with demineralised water to avoid residues of lime).  
8. Dry instruments (preferably with compressed air as recommended by the Robert Koch Institute). When drying the inside of the hollow applicator, please make sure to hold the compressed air gun close to the opening, in order to guarantee that sufficient air flows through the instrument.  
9. Visual examination to ensure that the instrument is clean and undamaged. If there are still visible residues of contamination, repeat the cleaning and chemical disinfecting process until no visible contamination is left. According to the recommendations of the Robert Koch Institute, the further reprocessing should be carried out mechanically. Information on the reprocessing of instruments in compliance with the standard EN ISO 17664 can be downloaded from our website [www.kometdental.de](http://www.kometdental.de) or requested directly from the manufacturer Gebr. Brasseler.

Sterilisation in the autoclave  
Make sure that the packaging is suitable for instruments classed as risk group critical B and the chosen method of sterilisation. Single pack: The packaging must be large enough to ensure that there is no pressure on the seal. Seal instruments in a suitable sterilisation container. For this procedure, the following instructions have to be observed: Steam sterilisation using a vacuum process at 134°C in a device that complies with the provisions of DIN EN 13060; with validated processes.  
- Fractionated pre-vacuum (type B)  
- Sterilisation temperature: 134°C  
- Hold time: at least 5 minutes (full cycle)  
- Drying time: at least 10 minutes  
In order to prevent staining and corrosion, the steam must be free of particles. The recommended limits for particle contents in feed water and condensed steam are defined by standard DIN EN 13060. Make sure not to exceed the maximum capacity of the sterilizer when sterilizing several instruments. Follow the instructions of the device manufacturer.  
After sterilisation, check instruments for corrosion. Do not reuse corroded instruments.

Universally valid notes  
Observe the legal provisions regarding the reprocessing of medical products valid in your country (e.g. [www.rki.de](http://www.rki.de)). The manufacturer confirms that the above detailed reprocessing methods are suitable for preparing the above named instrument group to enable their reuse. The person carrying out the reprocessing is responsible for ensuring that the applied method is carried out with appropriate equipment, materials and trained personnel at the reprocessing site and that it actually achieves the desired result. To guarantee this, routine controls of the validated mechanical preparation methods are normally necessary. Any deviation from the above detailed process must be carefully checked by the person carrying out the reprocessing to ensure effectiveness and to avoid possible adverse consequences.

5. Safety and liability  
Worn and damaged instruments (blunt or bent instruments etc.) are to be discarded and immediately replaced by new ones. The above instructions regarding the handling of the instruments, cooling and contact pressure are to be strictly observed. Only use the instruments according to their intended purpose. The non-observance of these safety instructions can lead to injuries or damage to the power system. It is the responsibility of the user to check the products prior to use to ensure that they are suited for the intended purpose. In case of contributory negligence by the user, Gebr. Brasseler partially or totally declines liability for all resulting damages, particularly due to non-observance of our recommendations for use or warnings as well as inadvertent misuse by the user.

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