
bei beschädigter Verpackung nicht verwenden
 (nur zutreffend bei sterilen Instrumenten). Do not
 use in case of damaged packaging (only applies to
 sterile instruments). Ne pas utiliser si l'emballage
 est endommagé (ne s'applique qu'aux instruments
 stériles). No utilizar si el embalaje está dañado (se
 aplica solamente a instrumentos estériles). Non
 utilizzare se l'imballo è danneggiato (riguarda soltanto
 strumenti sterili)

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Diamond coated sonic tips for use in implantology, developed by Dr. Ivo Agabiti

1. Indication
 The diamond coated sonic tips SFS99 provided by Komet are used for oscillating work on the jaw bone during the creation of an implant cavity before inserting the implant. Prior to using the sonic tips, the alveolar crest is opened from crestal direction (for example with a fine rotary bone cutter such as the H254E or the serrated sonic tip SFS101). The alveolar ridge must be sufficiently flexible. The sonic tips SFS99 are then used for creating the implant cavity, starting with a small size which is followed by other sonic tips in ascending order of size, until the required target size is reached.
 The sonic tips may only be used by dentists and medical personnel that are sufficiently trained in the relevant fields of dental surgery.

2. Contraindication
 When a patient is considered for treatment, all general contraindications for dental surgery must be taken into account. All commonly known contraindications in oral maxillo facial surgery (e.g., implantation, jaw surgery) also have to be considered.

3. Warnings
 The useful life of the sonic tips is limited. They therefore have to be used with the utmost care.
 Attention:
 - Risk of fracture! Spontaneous and premature fracture must be expected. The sonic tips may fracture unexpectedly and without warning. In order to reduce the risk of fracture, make sure to stick to the correct power setting and observe all recommendations for use. Never allow the instruments to jam. Do not use instruments for leverage. Avoid excessive contact pressure. Use sonic tips with the utmost care as even minor mechanical damage can provoke the fracture of the instruments.
 - Risk of ingestion! Fracture of the sonic tip leads to the risk of the tip being accidentally swallowed. Take all appropriate steps to prevent this.
 - Premature wear! Work on hard dental substance can cause wear of the diamond coating. Likewise, excessive contact pressure can lead to damage to the working part or the diamond coating. Inspect instruments after every use for signs of wear or damage. Exchange if necessary.
 - Excessive heat generation! Excessive contact pressure and/or the use of worn/blunt instruments can cause hazardous heat generation which increases the risk of thermal necrosis. Excessive contact pressure therefore has to be avoided. Damaged sonic tips (changed shape, corroded surfaces, void spots or defect thread etc.) have to be replaced. The sonic tips are medical devices and as such subject to the applicable legal provisions of the relevant country. In compliance with such regulations, the sonic tips may only be used by properly qualified operators within the scope of the indicated application and according to the valid industrial safety regulations and observing the intended rules of accident prevention. The provisions further stipulate that all equipment used must be free of defects and suitable for the intended application. Furthermore, the patient and third persons must be protected from hazards. Possible cross-contamination caused by the product must be ruled out.
 Safety control:
 Fracture may occur due to the tips being subjected to continuous or sudden stress (for example if the tips are dropped on the floor or through a mechanical change of the original shape). Check the tips prior to each use to make sure they are safe to use. To do this, apply slight pressure with your thumb and index finger. For additional safety, load the tips mechanically with approx. 10 N (1 kg), without function. Always wear protective gloves when checking, inserting or removing the sonic tips to avoid risks of injury and infection. We recommend that the tips be exchanged regularly.

4. Start-up
 The sonic tips are supplied in a non-sterile condition and have to be cleaned, disinfected and sterilised prior to first use.
 The Komet sonic tips SFS99 can be used in the following power systems:
 - in the Komet sonic hand piece SF1LM/SF1LS
 - in the SONICflex™ hand piece made by the co. KaVo Series 2003N/L/X/LX
 Insert and remove sonic tips and cooling adapters with suitable tip changers (i.e. Komet tip changer SF1975). When using the tip changer, make sure that the tip faces towards the semi-radius during insertion (fig. 1). Screw the chosen sonic tip into the hand piece with the sagid ring, turning the tool in a clockwise direction. To remove the tip, place tip changer on the tip and turn anticlockwise. Please observe the directions in the instructions supplied with your sonic hand piece.

5. Proper use
 The sonic tips may only be used with the authorised power systems and permitted power settings. The power setting 1-3 may be used. This applies to both the Komet sonic hand piece SF1LM and the KaVo SONICflex™ 2003N/L/X/LX. In both power systems, the power setting is regulated by means of a regulating ring. Please observe the recommendations in the instructions for use provided with the sonic hand pieces. The sonic tips are at their most effective when used with low contact pressure. Excessive contact pressure may lead to undesirable effects (see section 3 "warnings"). Given that these sonic tips may fracture without warning even when used according to instruction and their designated use, always plan two sonic tips for the operation, i.e. always have a reserve at hand so that the planned operation can still be finished if one tip has broken.
 Provision of cooling
 The surgical sonic tips are cooled with at least 50 ml/min of an externally fed sterile cooling liquid. Make sure to turn off the water supply of the dental treatment unit. It is absolutely essential to thoroughly and continuously irrigate the site at all times during the operation. This can be ensured by means of the cooling adapter SF1979 (fig. 3). One end of the cooling adapter is screwed into the sonic hand piece and the sonic tip is screwed into the other end of the cooling adapter. The flexible tube supplying the cooling agent is placed onto the small tube of the cooling adapter (fig. 2). The sonic tip can be removed from the cooling adapter by means of the mounting wrench 566. Controlled supply of the cooling agent by means of a suitable surgical power system is preferable, for example with the KaVo INTRASurg 1000. There are a number of compression systems for pressure supported transport of the cooling agent available on the market, for example:
 - System OMNIA, manual supply of sterile cooling agents by means of a compression pouch, available in Germany from the company Implants Medical Devices, www.implantis.de
 From Germany: (0180) 5 90 75 90, free fax: (00800) 05 05 51 23
 From Austria: (0049 180) 5 90 75 90, free fax: (00800) 05 05 51 23
 Article no. H2200 Manual cooling system with compression pouch (Omnia no. 30E0008.00), article no. H2204 Manual cooling system with flux control, to be attached to the hand piece (Omnia no. 32F0157). For further information, visit www.omniar.com.

6. Maintenance
 Place sonic tip and/or cooling adapter in a cleaning/disinfection tank (Fräsator) filled with a suitable detergent/disinfectant (e.g. Komet DC1) immediately after use in the mouth to prevent drying of residues on the instruments (protein fixation). Hold sonic tip and/or cooling adapter at an angle when placing them into the liquid to ensure that all hollow spaces are reached. It is recommended to reprocess the instruments within one hour of use at the very latest. The instruments should be in the cleaning/disinfection tank (Fräsator) when transported to the site where the reprocessing is to take place.

7. Reprocessing
 Manual pre-cleaning (prior to mechanical and manual reprocessing):
 1. Remove sonic tip and cooling adapter from cleaning/disinfection tank (Fräsator) immediately before reprocessing.
 2. Always separate the sonic tip from the cooling adapter before the start of the cleaning process with the tip changer to ensure proper reprocessing of both components.
 3. Remove stubborn contamination with a nylon brush under running water, turning the sonic tip and the rinse adapter constantly.
 4. Check patency of the cooling passage by guiding the nozzle cleaner (accessory of the sonic hand piece) through the passage. Replace instrument if the passage is obstructed.
 Mechanical reprocessing (validated):

1. Rinse sonic tip/cooling adapter thoroughly under running water to prevent any residues of the detergent/ disinfectant from getting into the washer/disinfecter.
 2. Unscrew the nozzle on the injector bar of the washer/disinfecter. Then screw the external thread of the rinse adapter SF1978 in the silicon holders on the injector bar. Now screw the sonic tip in the internal thread at the top of the rinse adapter (Fig. 4). When using the rinse adapter SF1979, proceed as follows: Unscrew the nozzle on the injector bar of the washer/disinfecter. Then screw the large external thread of the rinse adapter SF1977 in the silicon holders on the injector bar. Now screw the cooling adapter in the external thread at the top of the rinse adapter. Observe the directions in the instructions for use provided with the rinse adapter.
 3. Put chemical detergent into the washer/disinfecter, following the indications on the label and the instructions of the manufacturer of the washer/disinfecter.
 4. Start the Vario TD programme including thermal disinfection. Thermal disinfection takes place allowing for the A₀ value and observing national provisions (prEN/ISO 15883).
 5. On completion of the cycle remove sonic tip/cooling adapter from the washer/disinfecter and dry (preferably with compressed air as recommended by the Commission for Hospital Hygiene and Infectious Disease Prevention of the Robert Koch Institute). When drying the inside of the instruments, please make sure to hold the compressed air gun close to the cooling passage, in order to guarantee that sufficient air flows through the sonic tip.
 6. Visual examination to ensure that the instrument is clean and undamaged. If after mechanical reprocessing there are still visible residues of contamination, repeat the cleaning and disinfecting process until no visible contamination is left.
 Manual reprocessing (alternative):
 1. Place sonic tip/cooling adapter onto a 10 ml syringe and rinse with detergent/disinfectant. In case of visible contamination emerging from the perforations, rinse again with detergent/disinfectant.
 2. Remove stubborn contamination with a nylon brush under running water, turning the sonic tip/cooling adapter constantly.
 3. Rinse sonic tip/cooling adapter thoroughly with running water.
 4. Visual examination to ensure that the instrument is clean. If there are still visible residues of contamination, repeat the cleaning process until no visible contamination is left.
 5. Place sonic tip/cooling adapter in a suitable container (e.g. Komet ref. 9952) into the ultra sonic device or instrument bath filled with detergent/disinfectant.
 6. During chemical disinfection in the ultrasonic device or instrument bath, observe the instructions of the manufacturer regarding concentration and immersion time. Be sure to observe the full correct immersion time which does not start until the last tip has been immersed. Attention: do not exceed 45°C (risk of protein coagulation)!
 7. On completion of the immersion time, rinse sonic tip/cooling adapter thoroughly with suitable water (preferably with demineralised water to avoid residues of lime). Rinse the perforation of the sonic tip with at least 10 ml of demineralised water by means of a cannula, to ensure that no residues of the detergent are left in the cooling passage.
 8. Dry sonic tip (preferably with compressed air as recommended by the Commission for Hospital Hygiene and Infectious Disease Prevention of the Robert Koch Institute). When drying the inside of the instruments, please make sure to hold the compressed air gun close to the cooling passage, in order to guarantee that sufficient air flows through the sonic tip.
 9. Visual examination to ensure that the instrument is clean and undamaged. If there are still visible residues of contamination, repeat the cleaning and chemical disinfecting process until no visible contamination is left.
 According to the recommendations of the Commission for Hospital Hygiene and Infectious Disease Prevention of the Robert Koch Institute, reprocessing should preferably be carried out mechanically. Work instructions with regard to proper reprocessing of instruments according to DIN EN ISO 17664, the present instructions for use, as amended, and more detailed information can be downloaded from our web site www.kometdental.de or requested from the manufacturer Gebr. Brasseler.
 Sterilisation in the autoclave:
 Make sure that the packaging is suitable for sonic tips ensured that there is no pressure on the seal. In a suitable sterilisation container: Place sonic tips into a suitable sterilisation container (e.g. ref. 9952). For sterilisation, the container has to be sealed in suitable packaging. For this procedure, the following instructions have to be observed:
 Steam sterilisation using a vacuum process at 134°C in a device that complies with the provisions of DIN EN 13060; with validated processes.
 - fractionated pre-vacuum (type B)
 - sterilisation temperature: 134°C
 - hold time: at least 5 minutes (full cycle)
 - drying time: at least 10 minutes
 In order to prevent staining and corrosion, the steam must be free of particles. The recommended limits for particle contents in feed water and condensed steam are defined by standard DIN EN 13060. Make sure not to exceed the maximum capacity of the sterilizer when sterilizing several instruments. Follow the instructions of the device manufacturer.
 Universal valid notes:
 Observe the legal provisions regarding the reprocessing of medical products valid in your country (e.g. www.rki.de). The manufacturer confirms that the above detailed reprocessing methods are suitable for preparing the above named instrument group to enable their reuse. The person carrying out the reprocessing is responsible for ensuring that the applied method is carried out with appropriate equipment, materials and trained personnel at the reprocessing site and that it actually achieves the desired result. To guarantee this, routine controls of the validated mechanical preparation methods are necessary. Any deviation from the above detailed process must be carefully checked by the person carrying out the reprocessing to ensure effectiveness and to avoid possible adverse consequences.
 Accessories:
 Komet tip changer SF1975
 Komet cooling adapter SF1979 + mounting wrench 566 (in the set 4602)
 For reprocessing in the Miele washer/disinfecter: Komet rinse adapter SF1978 for Komet sonic tips and Komet rinse adapter SF1977 for the Komet cooling adapter SF1979
 Komet bur block/sterilisation container 9952

8. Storage
 Store packed instruments away from UV radiation and heat in a clean and dry environment. Do not store in the same room as solvents or chemicals.

9. Disposal
 The instruments have to be disposed of in tight, non-breakable and non-pierceable containers (protection from contamination).

10. Safety and liability
 Only use the instruments according to their intended purpose. The non-observance of these safety instructions can lead to injury or damage to the power system. It is the responsibility of the user to check the products prior to use to ensure that they are suited for the intended purpose. In case of contributory negligence by the user, Gebr. Brasseler partially or totally declines liability for all resulting damages, particularly due to non-observance of our recommendations for use or warnings as well as inadvertent misuse by the user.

11. Procedure in case of serious adverse events
 In case of serious adverse events, please proceed in line with the legal regulations in your country.

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