Komet Dental
Gebr. Brasseler GmbH & Co. KG
Trophagener Weg 25
32657 Lemgo · Germany
Telefon +49 (0) 5261 701-0
Telefax +49 (0) 5261 701-289
info@kometdental.de
www.kometdental.de

Instrumentenständer 97507 und 97607
© 02/2014, Gebr. Brasseler GmbH & Co. KG
309569 - Rev 02/2014

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Dental bur block for receiving ultrasonic tips in a torque wrench

1. Indication

Stainless steel bur block for the storage and sterilisation of ultrasonic tips in a torque wrench.

2. Contra-indications

Thermal cleaning and disinfection of ultrasonic tips in a torque wrench.

3. Warnings

The bur block is closed by means of a lid with clamping fit. Risk of injury!

Open and close the lid with care to avoid accidental trapping of your fingers. Inappropriate use can lead to injury to the user.

4. Recommendations for use

The bur block is designed for the storage and sterilisation of ultrasonic tips.

5. Recommendations for use

- Open lid before filling the bur block.
- The ultrasonic tips inside the torque wrench are placed into the round perforations in the insert inside the bur block. The working parts of the tips have to face upwards.

6. Maintenance

Clean, disinfect and sterilize bur block prior to first use.

7. Reprocessing

Warning: The following descriptions of the reprocessing methods (thermal cleaning and disinfection) only apply to empty bur blocks. For recommendations on cleaning and disinfecting ultrasonic tips, please refer to the relevant instructions for use enclosed with the ultrasonic tips.

Manual reprocessing

- Remove any contamination with a nylon brush under running water. Place bur block in an ultrasonic device or instrument bath filled with a suitable detergent/disinfectant.
- To avoid stains, rinse off all residues of disinfecting agent thoroughly under running water and dry the bur block by air jet.

Mechanical reprocessing

- Manual pre-cleaning: Clean bur block under running water with a nylon brush and rinse thoroughly in order to remove any contamination.
- Remove insert and place separately onto one of the storage trays of the washer/disinfector. Place the bur block in the washer/disinfector with its lid open to ensure that the inner surfaces are reached by the spray jet.
- Put chemical detergent into the washer/disinfector, following the indications on the label and the instructions of the manufacturer of the washer/disinfector.
- Start the Vario TD programme including thermal disinfection. Thermal disinfection takes place allowing for the A₀ value and observing national provisions (prEN/ISO 15883).
- On completion of the cycle remove bur block and insert from the washer/disinfector and dry preferably with compressed air (as recommended by the Commission for Hospital Hygiene and Infectious Disease Prevention).

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Instructions for use

Please read carefully and retain for future reference.

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6. Visual examination to ensure that the bur block is clean and undamaged. If after mechanical reprocessing there are still visible residues of contamination, repeat the cleaning and disinfecting process until no visible contamination is left. Place insert back into the bur block and close the bur block.

According to the recommendations of the Commission for Hospital Hygiene and Infectious Disease Prevention, the further reprocessing should be carried out mechanically. General information on the reprocessing of instruments in compliance with the standard ISO 17664 can be downloaded from our website www.kometdental.de or requested directly from the manufacturer Gebr. Brasseler.

Sterilisation in the autoclave (with or without inserted ultrasonic tips):

Warning: Prior to sterilisation of bur blocks with inserted ultrasonic tips, ensure that all ultrasonic tips to be sterilised have undergone a manual or mechanical reprocessing cycle according to the instructions of the manufacturer. For recommendations for care and maintenance as well as correct reprocessing please refer to the instructions of the manufacturer of the ultrasonic tips. For sterilisation, the bur block has to be sealed in suitable packaging. Make sure that the packaging/foil is suitable for the chosen method of sterilisation. The packaging must be large enough to ensure that there is no pressure on the seal.

For this procedure, the following instructions have to be observed:

- Steam sterilisation using a vacuum process at 134 °C in a device that complies with the provisions of DIN EN 13060; with validated processes.
- Fractionated pre-vacuum (type B)
- Sterilisation temperature: 134 °C
- Hold time: at least 5 minutes (full cycle)
- Drying time: at least 10 minutes

In order to prevent staining and corrosion, the steam must be free of particles. The recommended limits for particle contents in feed water and condensed steam are defined by standard DIN EN 13060. Make sure not to exceed the maximum capacity of the sterilizer when sterilizing several instruments. Follow the instructions of the device manufacturer.

Universally valid notes:

Observe the legal provisions regarding the reprocessing of medical products valid in your country (e.g. www.rki.de). The manufacturer confirms that the above detailed reprocessing methods are suitable for preparing the above named instrument group to enable their reuse. The person carrying out the reprocessing is responsible for ensuring that the applied method is carried out with appropriate equipment, materials and trained personnel at the reprocessing site and that it actually achieves the desired result. To guarantee this, routine controls of the validated mechanical and/or manual preparation methods are necessary. Any deviation from the above detailed process must be carefully checked by the operator to ensure effectiveness and to avoid possible adverse consequences.

8. Liability

It is the responsibility of the user to check the products prior to use to ensure that they are suited for the intended purpose. In case of contributory negligence by the user, Gebr. Brasseler partially or totally declines liability for all resulting damages, particularly due to non-observance of our recommendations for use or warnings as well as inadvertent misuse by the user.

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