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EN

Sonic tips for interproximal cavity preparation

- Indication**
For use as part of restorative dental treatments
Preparation of teeth for treatments with direct fillings and indirect inlays as well as partial crowns and finishing of preparation margins of premolars and molars. The following restorations can be used subsequently:
- direct fillings made of tooth-coloured composites
- ceramic inlays made of pressed ceramics
- ceramic inlays made of CAD/CAM milled ceramics
- partial crowns made of pressed ceramics
- partial crowns made of CAD/CAM milled ceramics
- Contraindication**
There are no known contra-indications related to the instruments.
- Warnings**
The useful life of the sonic tips is limited. They therefore have to be used with the utmost care.
Attention:
- Risk of fracture! Spontaneous and premature fracture must be expected. The sonic tips may fracture unexpectedly and without warning. In order to reduce the risk of fracture, make sure to stick to the correct power setting (see instructions for use) and observe all recommendations for use. Never allow the instruments to jam. Do not use instruments for leverage. Avoid excessive contact pressure. Use sonic tips with the utmost care as even minor mechanical damage can provoke the fracture of the instruments.
- Risk of ingestion! Fracture of the sonic tip leads to the risk of the tip being accidentally swallowed. Take all appropriate steps to prevent this.
- Premature wear! Work on hard dental substance can cause premature wear of the diamond coating. Likewise, excessive contact pressure can lead to damage to the working part or the diamond coating. Inspect instruments after every use for signs of wear or damage. Exchange if necessary.
- Excessive heat generation! Excessive contact pressure and/or worn or blunt instruments can cause hazardous heat generation which increases the risk of thermal necrosis (damage to tissue). Excessive contact pressure therefore has to be avoided. Damaged sonic tips (changed shape, corroded surfaces, void spots or defect thread etc.) have to be replaced.
Incorrect use can lead to the above detailed undesirable effects as well as to injury to the user, the patient or third persons.
Safety control:
Fracture may occur due to the tips being subjected to continuous or sudden stress (for example if the tips are dropped on the floor or through a mechanical change of the original shape). Check the tips prior to each use to make sure they are safe to use. To do this, apply slight pressure with your thumb and index finger (fig. 1). For additional safety, load the tips mechanically with approx. 10 N (1 kg), without function (fig. 2). Always wear protective gloves when checking, inserting or removing the sonic tips to avoid risks of injury and infection. We recommend that the tips be exchanged regularly.
- Start-up**
Those sonic tips that are supplied non-sterile have to be prepared prior to first use.
The sonic tips can be used:
- in the Komet sonic hand piece SF 1LM/S
- in the SONICflex™ hand piece made by the co. KaVo (series 2000N/L/LX or series 2003N/L/X/LX)
- in scalers made by co. W&H (i.e. series Proxeo® ZA-55/L/LM/M/LS and ST ZE-55RM/BC, series Synea® ZA-55/L/LM/M or Alegria® ST ZE-55RM/BC)
- in the SIROAIR L provided by co. Sirona®
Insert and remove sonic tips with the tip changer SF 1975 with torque. Screw the chosen sonic tip into the hand piece with the tip changer, turning the tool in clockwise direction. To remove the tip, place tip changer on the tip and turn anticlockwise.
Note: Please observe the directions in the instructions supplied with your sonic hand piece.
- Proper use**
The sonic tip is used for final shaping of the cavities in the interproximal area and for smoothing the interproximal cavity margins. Remove old restorations with rotary instruments according to the usual methods, if applicable.
- Always carry out the basic preparation with rotary instruments (i.e., prepare the cavity with the instruments from the Expert Set 4562/4562ST provided by Komet).
- At the approximal side, prepare almost up to the adjacent tooth, merely leaving a thin enamel wall.
- Remove the enamel wall with a manual tool or a flame shaped rotary diamond instrument (use a steel matrix to protect the adjacent tooth).
- To prepare the mesial or distal cavity to the correct shape, use one of the above mentioned sonic tips at power level 3.
- Apply low contact pressure. Excessive contact pressure will lessen the oscillation of the sonic tip and reduce its efficiency.
- Start sonic tip with plenty of water spray prior to applying the tip to the tooth.
- The interproximal cavity margin is shaped and smoothed performing vestibular/oral movements. The sonic tip is guided along the cavity margin in mesio/distal direction in order to remove any instable enamel structures.
- For fine finishing of the preparation margin, work at power level 1.
- Do not use power level 2 when working with these sonic tips as this power level is not powerful enough for adequate substance removal during the preparation but too powerful for the fine finishing of the preparation margin.
Finally, continue the treatment of the prepared tooth as usual.
- Maintenance**
Place sonic tips in a cleaning/disinfection tank (Fräsator) filled with a suitable detergent/disinfectant (e.g. Komet DC1) immediately after use in the mouth to prevent drying of residues on the instruments (protein fixation). Hold sonic tips at an angle when placing them into the liquid to ensure that all hollow spaces are reached. It is recommended to reprocess the instruments within one hour of use at the very latest. The instruments should be in the cleaning/disinfection tank (Fräsator) when transported to the site where the reprocessing is to take place.
- Reprocessing**
The tips are supplied non-sterile and have to be prepared prior to first use.
Manual pre-cleaning (prior to mechanical and manual reprocessing):
1. Remove sonic tip from cleaning/disinfection tank (Fräsator) immediately before mechanical reprocessing.
2. Remove surface contamination thoroughly under running water. Completely remove stubborn contamination with a nylon brush below water level, turning the tip constantly. Rinse again under running water.
3. Check patency of the cooling passage by guiding the nozzle cleaner through the passage. Replace instrument if the passage is obstructed.

- Mechanical reprocessing (validated):
1. Rinse sonic tip thoroughly under running water to prevent any residues of the detergent/disinfectant from getting into the washer/disinfectator.
2. Unscrew the nozzle on the injector bar of the washer/disinfectator or pull it out. Then screw or push the external thread of the rinse adapter SF 1978/L in the silicon holders on the injector bar. Now screw the sonic tip in the internal thread at the top of the rinse adapter. (Fig. 3 shows the rinse adapter SF 1978/L). Observe the directions in the instructions for use provided with the SF 1978/L.
3. Put chemical detergent into the washer/disinfectator, following the indications on the label and the instructions of the manufacturer of the washer/disinfectator.
4. Start the Vario TD programme including thermal disinfection. Thermal disinfection takes place allowing for the A₀ value and observing national provisions (prEN/ISO 15883).
5. On completion of the cycle remove sonic tip from the washer/disinfectator and dry (preferably with compressed air as recommended by the Commission for Hospital Hygiene and Infectious Disease Prevention of the Robert Koch Institute). When drying the inside of the instruments, please make sure to hold the compressed air gun close to the cooling passage, in order to guarantee that sufficient air flows through the sonic tip.
6. Visual examination to ensure that the instrument is clean and undamaged. If after mechanical reprocessing there are still visible residues of contamination, repeat the cleaning and disinfecting process until no visible contamination is left.
Manual reprocessing (alternative):
1. Place sonic tip onto a 10 ml syringe and rinse with detergent/disinfectant. In case of visible contamination emerging from the perforation, rinse again with detergent/disinfectant.
2. Remove stubborn contamination with a nylon brush under running water, turning the sonic tip constantly.
3. Rinse sonic tip thoroughly with running water.
4. Visual examination to ensure that the instrument is clean. If there are still visible residues of contamination, repeat the cleaning process until no visible contamination is left.
5. Place sonic tip in a suitable container into the ultra sonic device or instrument bath filled with detergent/disinfectant.
6. During chemical disinfection in the ultrasonic device or instrument bath, observe the instructions of the manufacturer regarding concentration and immersion time. Please note that due to their internal bore, the chemical disinfection of sonic and ultrasonic tips (with internal irrigation) in the ultrasonic bath is always carried out in a 2% solution and takes 10 minutes. Be sure to observe the full correct immersion time which does not start until the last tip has been immersed. Attention: do not exceed 45 °C (risk of protein coagulation!)
7. On completion of the immersion time, rinse sonic tip thoroughly with suitable water (preferably with demineralised water to avoid residues of lime). Rinse the perforation of the sonic tip with at least 10 ml of demineralised water by means of a cannula, to ensure that no residues of the detergent are left in the cooling passage.
8. Dry sonic tip (preferably with compressed air as recommended by the Commission for Hospital Hygiene and Infectious Disease Prevention of the Robert Koch Institute). When drying the inside of the instruments, please make sure to hold the compressed air gun close to the cooling passage, in order to guarantee that sufficient air flows through the sonic tip.
9. Visual examination to ensure that the instrument is clean and undamaged. If there are still visible residues of contamination, repeat the cleaning and chemical disinfecting process until no visible contamination is left.
According to the recommendations of the Commission for Hospital Hygiene and Infectious Disease Prevention of the Robert Koch Institute, reprocessing should preferably be carried out mechanically. Work instructions with regard to proper reprocessing of instruments according to DIN EN ISO 17664, the present instructions for use and more detailed information, as amended, can be downloaded from our web site www.kometdental.de or requested from the manufacturer Gebr. Brasseler.
Thermal disinfection in the autoclave:
Classed as risk group semi-critical B, the sonic tips can be subjected to thermal sterilisation in a steam steriliser, in unwrapped condition inside a suitable container (e.g. 9952).
For this procedure, the following instructions have to be observed:
Steam sterilisation using a vacuum process at 134 °C in a device that complies with the provisions of DIN EN 13060; with validated processes.
- fractionated pre-vacuum (type B)
- sterilisation temperature: 134 °C
- hold time: at least 5 minutes (full cycle)
- drying time: at least 10 minutes
In order to prevent staining and corrosion, the steam must be free of particles. The recommended limits for particle contents in feed water and condensed steam are defined by standard DIN EN 13060. Make sure not to exceed the maximum capacity of the sterilizer when sterilizing several instruments. Follow the instructions of the device manufacturer.

Universally valid notes:
Observe the legal provisions regarding the reprocessing of medical products valid in your country (e.g. www.rki.de). The manufacturer confirms that the above detailed reprocessing methods are suitable for preparing the above named instrument group to enable their reuse. The person carrying out the reprocessing is responsible for ensuring that the applied method is carried out with appropriate equipment, materials and trained personnel at the reprocessing site and that it actually achieves the desired result. To guarantee this, routine controls of the validated mechanical and/or manual preparation methods are necessary. Any deviation from the above detailed process must be carefully checked by the operator to ensure effectiveness and to avoid possible adverse consequences.

8. Storage
Store packed instruments away from UV radiation and heat in a clean and dry environment. Do not store in the same room as solvents or chemicals.

9. Disposal
The instruments have to be disposed of in tight, non-breakable and non-pierceable containers (protection from contamination).

10. Safety and liability
Only use the instruments according to their intended purpose. The non-observance of these safety instructions can lead to injuries or damage to the power system. It is the responsibility of the user to check the products prior to use to ensure that they are suited for the intended purpose. In case of contributory negligence by the user, Gebr. Brasseler partially or totally declines liability for all resulting damages, particularly due to non-observance of our recommendations for use or warnings as well as inadvertent misuse by the user.

11. Procedure in case of serious adverse events
In case of serious adverse events, please proceed in line with the legal regulations in your country.

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© Gebr. Brasseler GmbH & Co. KG