
Komet Dental
Gebr. Brasseler GmbH & Co. KG
Trophagener Weg 25
32657 Lemgo · Germany
Telefon +49 (0) 5261 701-0
Telefax +49 (0) 5261 701-289
info@kometdental.de
www.kometdental.de

Endodontie Schallspitzen

© Gebr. Brasseler GmbH & Co. KG
310450 - Rev 2018-07

| | Wurzelkanalaufbereitung, Root canal preparation, Traitement endodontique, Preparación de canal radicular, Preparazione canalare | | Strahlensterilisiert (Sterile Instruments sind auf dem Verpackungsetikett entsprechend gekennzeichnet), Sterilized using irradiation (Sterile instruments are marked as such on the label of the packaging), Stérilisation par irradiation (Les instruments stériles sont marqués comme tels sur l'étiquette de l'emballage), Esterilización por radiación (Los instrumentos estériles son marcados como tales sobre la etiqueta del empaque), Sterilizzazione per mezzo di radiazione (Gli strumenti sterili sono marcati come tali sull'etichetta dell'imballaggio) |
|--|---|--|---|
|--|---|--|---|

DE

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EN

Sonic tips for orthograde endodontic treatments

- Indication**
Orthograde preparation of the pulp chamber and preparation of the cervical third of the root canal as part of an endodontic treatment. The following tips are currently available for this indication: SF66, SF67, SF68, SF69 and SF70.
Indications of each individual tip:
SF66: Initial preparation of the access cavity, removal of protruding material, straightening and rough preparation of the floor of the chamber. Shape: Large bud;
SF67: Retrieval of fine and calcified canals, opening of the upper canal portions during revision, initial pre-drilling in case of intended post reconstructions and creation of retentions in the canal access area, during preparations instead of Gates-Glidden burs. Shape: tapered 125°;
SF68: Alternative to the SF67 with a more pronounced angle. Shape: tapered 112°;
SF69: Finishing of the access cavity and minimal freshening of the dentin layer, removal of residues of root canal fillings, creation of retentions in the canal access area. Work with brush-like movements. Shape: Small bud;
SF70: Enlargement of long and wide canals, unblocking of fractured instruments, removal of root fillings made of gutta-percha and soft cements. Shape: tapered 122°.
- Contra-indications**
There are no known contra-indications.
- Warnings**
The useful life of the sonic tips is limited. They therefore have to be used with utmost care.
Attention:
- **Risk of fracture!** Spontaneous and premature fracture must be expected. The sonic tips may fracture unexpectedly and without warning. In order to reduce the risk of fracture, make sure to stick to the correct power setting (see instructions for use) and observe all recommendations for use. Never allow the instruments to jam. Do not use instruments for leverage. Avoid excessive contact pressure. Use sonic tips with utmost care as even minor mechanical damage can provoke the fracture of the instruments.
- **Risk of ingestion!** Fracture of the sonic tip leads to the risk of the tip being accidentally swallowed. Take all appropriate steps to prevent this (for example by using a rubber dam).
- **Premature wear!** Work on hard dental substance can cause premature blunting of the blades and/or wear of the diamond coating. Likewise, excessive contact pressure can lead to damage to the working part, the blades or the diamond coating. Inspect instruments after every use for signs of wear or damage. Exchange if necessary.
- **Excessive heat generation!** Excessive contact pressure and/or worn or blunt instruments can cause hazardous heat generation which increases the risk of thermal necrosis (damage to tissue). Excessive contact pressure therefore has to be avoided. Damaged sonic tips (changed shape, corroded surfaces, void spots or defect thread etc.) have to be replaced.
Incorrect use can lead to the above detailed undesirable effects as well as to injury to the user, the patient or third persons.
Safety control:
Fracture may occur due to the tips being subjected to continuous or sudden stress (for example if the tips are dropped on the floor or through a mechanical change of the original shape). Check the tips prior to each use to make sure they are safe to use. To do this, apply slight pressure with your thumb and index finger (fig. 1). For additional safety, load the tips mechanically with approx. 10 N (1 kg), without function (fig. 2). Always wear protective gloves when checking, inserting or removing the sonic tips to avoid risks of injury and infection. We recommend that the tips be exchanged regularly.

- Appropriate use**
Those sonic tips that are supplied non-sterile have to be prepared prior to first use.
The sonic tips can be used:
- In the Komet sonic hand piece SF1LM/S
- In the SONICflex™ hand piece made by the co. KaVo (Series 2000N/L/LX or series 2003N/L/X/LX)
- In scalers made by co. W&H (i.e. series Proxoe® ZA-55/L/LM/M/LS and ST ZE-55RM/BC, series Synea® ZA-55/L/LM/M or Alegra® ST ZE-55RM/BC)
- In the SIROAIR L provided by co. Sirona®
The tip changer with torque SF1975 is used for inserting and removing the sonic tips. Screw the chosen sonic tip into the hand piece with the mentioned tip changer, turning the tool in clockwise direction. To remove the tip, place tip changer on the tip and turn anticlockwise.
Note: Please observe the directions in the instructions supplied with your sonic hand piece.

- Recommendations for use**
- First prepare tooth with rotary instruments the usual way.
- Start up sonic tip prior to application to the tooth, providing plenty of spray cooling.
- When removing large quantities of substance, any of the above named power systems can be used at power setting 2 (this corresponds to 100% capacity of the SONICflex™ 2000).
- For finishing of surfaces and fine trimming, use power setting 1 (this corresponds to 75% capacity of the SONICflex™ 2000).
- The use of power setting 3 is not permitted when working with these sonic tips (see section 3 "warnings").
- Apply low contact pressure only, excessive contact pressure reduces the vibration of the sonic tip and consequently its effectiveness.
- After use of the sonic tips, complete the treatment of the endodontically treated tooth according to the usual methods.

- Maintenance/care**
Place sonic tips in a cleaning/disinfection tank filled with a suitable detergent/disinfectant (e.g. Komet DC1) immediately after use in the mouth to prevent drying of residues on the instruments (protein fixation). Hold sonic tips at an angle when placing them into the liquid to ensure that all hollow spaces are reached. It is recommended to reprocess the instruments within one hour of use at the very latest. The instruments should be in the cleaning/disinfection tank when transported to the site where the reprocessing is to take place.

- Reprocessing**
The instruments are supplied non-sterile and have to be prepared prior to first use.
Manual pre-cleaning (prior to mechanical and manual reprocessing):
1. Remove sonic tip from cleaning/disinfection tank (Fräsator) immediately before reprocessing.
2. Remove surface contamination thoroughly under running water. Completely remove stubborn contamination with a nylon brush below water level, turning the tip constantly.

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Instructions for use

Please read carefully and retain for future reference.

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- Check patency of the cooling passage by guiding the nozzle cleaner through the passage. Replace instrument if the passage is obstructed.

Mechanical reprocessing (validated):

- Rinse sonic tip thoroughly under running water to prevent any residues of the detergent/disinfectant from getting into the washer/disinfecter.
- Unscrew the nozzle on the injector bar of the washer/disinfecter or pull it out. Then screw or push the external thread of the rinse adapter SF1978/L in the silicon holders on the injector bar. Now screw the sonic tip in the internal thread at the top of the rinse adapter (Fig. 3 shows the rinse adapter SF1978L). Observe the directions in the instructions for use provided with the SF1978/L.
- Put chemical detergent into the washer/disinfecter, following the indications on the label and the instructions of the manufacturer of the washer/disinfecter.
- Start the Vario TD programme including thermal disinfection. Thermal disinfection takes place allowing for the A₀ value and observing national provisions (prEN/ISO 15883).
- On completion of the cycle remove sonic tip from the washer/disinfecter and dry (preferably with compressed air according to the recommendations of the Commission for Hospital Hygiene and Infectious Disease Prevention of the Robert Koch Institute). When drying the inside of the instruments, please make sure to hold the compressed air gun close to the cooling passage, in order to guarantee that sufficient air flows through the sonic tip.
- Visual examination to ensure that the instrument is clean and undamaged. If after mechanical reprocessing there are still visible residues of contamination, repeat the cleaning and disinfecting process until no visible contamination is left.

Manual reprocessing (alternative):

- Place sonic tip onto a 10 ml syringe and rinse with detergent/disinfectant. In case of visible contamination emerging from the perforations, rinse again with detergent/disinfectant.
- Remove stubborn contamination with a nylon brush under running water, turning the sonic tip constantly.
- Rinse sonic tip thoroughly with running water.
- Visual examination to ensure that the instrument is clean. If there are still visible residues of contamination, repeat the cleaning process until no visible contamination is left.
- Place sonic tip in a suitable container into the ultra sonic device or instrument bath filled with detergent/disinfectant.
- Please note that, due to their internal orifice, the chemical disinfection of sonic and ultrasonic tips in an ultrasonic bath using a 2% solution always takes 10 minutes. Be sure to observe the full correct immersion time which does not start until the last tip has been immersed. Attention: do not exceed 45°C (risk of protein coagulation)!
- On completion of the immersion time, rinse sonic tip thoroughly with suitable water (preferably with demineralised water to avoid residues of lime). Rinse the perforation of the sonic tip with at least 10 ml of demineralised water by means of a cannula, to ensure that no residues of the detergent are left in the cooling passage.
- Dry sonic tip (preferably with compressed air according to the recommendations of the Commission for Hospital Hygiene and Infectious Disease Prevention of the Robert Koch Institute). When drying the inside of the instruments, please make sure to hold the compressed air gun close to the cooling passage, in order to guarantee that sufficient air flows through the sonic tip.
- Visual examination to ensure that the instrument is clean and undamaged. If there are still visible residues of contamination, repeat the cleaning and chemical disinfecting process until no visible contamination is left.

According to the recommendations of the Commission for Hospital Hygiene and Infectious Disease Prevention of the Robert Koch Institute, reprocessing should preferably be carried out mechanically. Information on proper reprocessing of instruments according to DIN EN ISO 17664, the present instructions for use, as amended, and more detailed information can be downloaded from our web site www.kometdental.de or requested from the manufacturer Gebr. Brasseler.

Sterilisation in the autoclave:

Make sure that the packaging is suitable for sonic tips classed as risk group critical B and the chosen method of sterilisation. Single pack: The packaging must be large enough to ensure that there is no pressure on the seal. In the sterilisation container: Place sonic tips into a suitable sterilisation container (e.g. 9952). For sterilisation, the container has to be sealed in suitable packaging.
For this procedure, the following instructions have to be observed:
Steam sterilisation using a vacuum process at 134°C in a device that complies with the provisions of DIN EN 13060; with validated processes.
- Fractionated pre-vacuum (type B)
- Sterilisation temperature: 134°C
- Hold time: at least 5 minutes (full cycle)
- Drying time: at least 10 minutes
In order to prevent staining and corrosion, the steam must be free of particles. The recommended limits for particle contents in feed water and condensed steam are defined by standard DIN EN 13060. Make sure not to exceed the maximum capacity of the sterilizer when sterilizing several instruments. Follow the instructions of the device manufacturer.

Universally valid notes:

Observe the legal provisions regarding the reprocessing of medical products valid in your country (e.g. www.rki.de). The manufacturer confirms that the above detailed reprocessing methods are suitable for preparing the above named instrument group to enable their reuse. The person carrying out the reprocessing is responsible for ensuring that the applied method is carried out with appropriate equipment, materials and trained personnel at the reprocessing site and that it actually achieves the desired result. To guarantee this, routine controls of the validated mechanical preparation methods are necessary. Any deviation from the above detailed process must be carefully checked by the operator to ensure effectiveness and to avoid possible adverse consequences.

8. Storage

Store sterile packed instruments in their original packaging, away from light and heat in a clean and dry environment. Do not store in the same room as solvents or chemicals.

9. Safety and liability

See also "Instructions for use and safety recommendations" in the current Komet catalogue. The user is responsible for checking the products prior to use as to whether they are suited for the intended purpose. In case of contributory negligence by the user, Gebr. Brasseler partially or totally declines liability for all resulting damages, particularly if these are due to non-observance of our recommendations for use or warnings as well as inadvertent misuse by the user.

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