
	Komet Dental Gebr. Brasseler GmbH & Co. KG Trophagener Weg 25 32657 Lemgo · Germany Telefon+49 (0) 5261 701-0 Telefax+49 (0) 5261 701-289 info@kometdental.de www.kometdental.de	Gebrauchsanweisung Bitte unbedingt beachten und für späteren Gebrauch gut aufbewahren.	Instruções de empleo Leer detenidamente y guardar este ejemplar para uso más tarde.
	Schallspitzen SFS100/101/102 © Gebr. Brasseler GmbH & Co. KG 310452 - Rev 2016-07	Instructions for use Please read carefully and retain for future reference.	Istruzioni d'uso Si prega di leggere attentamente e conservare per un'eventuale consultazione.
		Conseils d'utilisation A respecter scrupuleusement et garder cette copie pour une utilisation ultérieure.	

Kieferchirurgie, Oral surgery, Chirurgie maxillo-faciale, Cirugía maxilofacial, Chirurgia maxillo-facciale


STERILE R Sterilisation durch Bestrahlung (Sterile Instrumente sind auf dem Verpackungsetikett entsprechend gekennzeichnet). Sterilized using irradiation (Sterile instruments are marked as such on the label of the packaging), Stérilisation par irradiation (Les instruments stériles sont marqués comme tels sur l'étiquette de l'emballage), Esterilización por radiación (Los instrumentos estériles son marcados como tales sobre la etiqueta del embalaje), Sterilizzazione per mezzo di radiazione (Gli strumenti sterili sono marcati come tali sull'etichetta dell'imballaggio)

Verwendbar bis (nur zutreffend bei sterilen Instrumenten), Use by (only applies to sterile instruments), A utiliser avant le (ne s'applique qu'aux instruments stériles), Fecha de caducidad (se aplica solamente a instrumentos estériles), Data di scadenza (riguarda soltanto strumenti sterili)

Bei beschädigter Verpackung nicht verwenden (nur zutreffend bei sterilen Instrumenten), Do not use in case of damaged packaging (only applies to sterile instruments), Ne pas utiliser si l'emballage est endommagé (ne s'applique qu'aux instruments stériles), No utilizar si el empaque está dañado (se aplica solamente a instrumentos estériles), Non utilizzare se l'imballo è danneggiato (riguarda soltanto strumenti sterili)

SFS 100 | sagittal
SFS 101 | axial
SFS 102 | gerade/straight

rostfreier Stahl - stainless steel - acier inoxydable - acero inoxidable - acero inoxidabile

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EN
Sonic tips for gentle surgery, developed by Dr. Agabiti

1. Indications
The Komet sonic tips SFS 100, 101 and 102 are designed for oscillating work on bones in addition to rotary work with hand pieces and contra angles. These sonic tips are particularly suited for the fragmentation of bones, alveolotomy, osteotomy – in particular for exposure of impacted teeth –, orthodontic corticotomy and for bone splitting. The saw blade of an SFS100/101/102 tip is provided with sharp blades. These sonic tips may only be used by dentists or medical personnel who are sufficiently qualified in dental implantology and the relevant fields of surgery.

2. Contraindications
When a patient is considered for treatment, all general contraindications for dental surgery must be taken into account. All commonly known contraindications in oral maxillo facial surgery (e.g. implantation, jaw surgery) also have to be considered.

3. Warnings
The useful life of the sonic tips is limited. They therefore have to be used with the utmost care.
Attention:
- Risk of fracture! The sonic tips may fracture unexpectedly and without warning. In order to reduce the risk of fracture, make sure to stick to the correct power setting (see instructions for use) and observe all recommendations for use. Never allow the instruments to jam. Do not use instruments for leverage. Avoid excessive contact pressure. Use sonic tips with the utmost care as even minor mechanical damage can provoke the fracture of the instruments.
- Risk of ingestion! Fracture of the sonic tip leads to the risk of the tip being accidentally swallowed. Take all appropriate steps to prevent this.
- Premature wear! Inspect instruments after every use for signs of wear or damage. Exchange if necessary. Work on hard dental substance can cause moderate wear. Likewise, excessive contact pressure can lead to damage to the working part.
- Excessive heat generation! Excessive contact pressure and/or worn or blunt instruments can cause hazardous heat generation which increases the risk of thermal necrosis (damage to tissue). Excessive contact pressure therefore has to be avoided. Damaged sonic tips (changed shape, corroded surfaces, void spots or defect thread etc.) have to be replaced. The sonic tips are medical devices and as such subject to the applicable legal provisions of the relevant country. In compliance with such regulations, the sonic tips may only be used by properly qualified operators within the scope of the indicated application and according to the valid industrial Health & Safety regulations and observing the valid rules for accident prevention. The provisions further stipulate that all equipment used must be free of defects and suitable for the intended application. Furthermore, the patient and third persons must be protected from hazards. Possible contamination caused by the product must be ruled out.
Safety control
Fracture may occur due to the tips being subjected to continuous or sudden stress (for example if the tips are dropped on the floor or through a mechanical change of the original shape). Check the tips prior to each use to make sure they are safe to use. To do this, apply slight pressure with your thumb and index finger (fig. 1). For additional safety, load the tips mechanically with approx. 10N (1kg), without function (fig. 2). Always wear protective gloves when checking, inserting or removing the sonic tips to avoid risks of injury and infection. We recommend that the tips be exchanged regularly.

mounting wrench to ensure proper reprocessing of both components.
3. Check the patency of the cooling passage by guiding the nozzle cleaner (accessory of the sonic hand piece) through the passage. Replace instrument if the passage is obstructed.

Mechanical reprocessing (validated):
1. Rinse sonic tip/cooling adapter thoroughly under running water to prevent any residues of the detergent/disinfectant from getting into the washer/disinfectant.
2. Unscrew the nozzle on the injector bar of the washer/disinfectant or pull it out. Then screw or push the external thread of the rinse adapter SF1978 in the silicon holders on the injector bar. Now screw the sonic tip in the internal thread at the top of the rinse adapter. When using the rinse adapter SF1979, proceed as follows: Unscrew the nozzle on the injector bar of the washer/disinfectant. Then screw the large external thread of the rinse adapter SF1977 in the silicon holders on the injector bar. Now screw the cooling adapter in the external thread at the top of the rinse adapter. Observe the directions in the instructions for use provided with the rinse adapter.
3. Put chemical detergent into the washer/disinfectant, following the indications on the label and the instructions of the manufacturer of the washer/disinfectant.
4. Start the Vario TD programme including thermal disinfection. Thermal disinfection takes place allowing for the Aa value and observing national provisions (prEN/ISO 15883).
5. On completion of the cycle remove sonic tip/cooling adapter from the washer/disinfectant and dry (preferably with compressed air as recommended by the Commission for Hospital Hygiene and Infectious Disease Prevention of the Robert Koch Institute). When drying the inside of the instruments, please make sure to hold the compressed air gun close to the cooling passage, in order to guarantee that sufficient air flows through the sonic tip.
6. Visual examination to ensure that the instrument is clean and undamaged. If after mechanical reprocessing there are still visible residues of contamination, repeat the cleaning and disinfecting process until no visible contamination is left.

4. Appropriate use
The tips are supplied non-sterile and have to be prepared prior to first use. The sonic tips for oral surgery made by Komet and developed by Dr. Ivo Agabiti (www.sonosurgery.it) can be used:
- in the Komet sonic hand piece SF1LM/SF1LS
- in the SONICflex™ hand piece made by the co. KaVo Series 2003N/L/X/LX
The above listed sonic hand pieces are provided with a MultiFlex™ coupling.
The tip changer with torque SF1975 is used for inserting and removing the sonic tips. Screw the chosen sonic tip into the hand piece with the mentioned tip changer, turning the tool in clockwise direction. To remove the tip, place tip changer on the tip and turn anticlockwise.
Note: Please observe the directions in the instructions supplied with your sonic hand piece.

5. Recommendations for use
The sonic tips may only be used with the authorised power systems and permitted power settings. The power settings 2 (gentle) and 3 (standard) may be used. This applies to both the Komet sonic hand piece SF1LM and the KaVo SONICflex™ 2003. The power setting is regulated by means of a regulating ring. Please observe the recommendations in the instructions for use provided with the sonic hand pieces. The sonic tips are at their most effective when used with low contact pressure. Excessive contact pressure may lead to undesirable effects (see section 3 "warnings"). Given that these sonic tips may fracture without warning even when used according to instruction and their designated use, always plan two sonic tips for the operation, i.e. always have a reserve at hand so that the planned operation can still be finished even if one tip has broken.
Provision of cooling
The surgical sonic tips are cooled with at least 50 ml/min of an externally fed, sterile cooling liquid. Make sure to turn off the water supply of the dental treatment unit. It is absolutely essential to thoroughly and continuously irrigate the site at all times during the operation. This can be ensured by means of the cooling adapter SF1979. One end of the cooling adapter is screwed into the sonic hand piece and the sonic tip is screwed into the other end of the cooling adapter. The flexible tube supplying the cooling agent is placed onto the small tube of the cooling adapter. The sonic tip can be removed from the cooling adapter by means of the mounting wrench 566. Controlled supply of the cooling agent by means of a suitable surgical power system is preferable, for example with the KaVo INTRASurg 1000 for oscillating work. There are various compression systems available on the market that ensure pressure supported transport of the cooling agent:
- OMNIA system: manual supply of a sterile cooling agent by means of a compression pouch.
Article no. H2200 Manual cooling system with compression pouch (Omnia no. 30E0008.00), article no. H2204 Manual cooling system with flux control, to be attached to the hand piece (Omnia no. 32F0157). For further information, visit www.omniasrl.com.
- Aqua-Press sterile irrigation system, made by the company Ronvig Dental A/S, Gl. Vejlevej 59, DK-8721 Dagaard, tel. 0045 70 23 34 11, fax: 0045 76 74 07 98, www.ronvig.com

Manual reprocessing (alternative):
1. Place sonic tip/cooling adapter onto a 10 ml syringe and rinse with detergent/disinfectant. In case of visible contamination emerging from the perforations, rinse again with detergent/disinfectant.
2. Remove stubborn contamination with a nylon brush under running water, turning the sonic tip constantly.
3. Rinse sonic tip/cooling adapter thoroughly with running water.
4. Visual examination to ensure that the instrument is clean. If there are still visible residues of contamination, repeat the cleaning process until no visible contamination is left.
5. Place sonic tip/cooling adapter in a suitable container (e.g. Komet 9952) into the ultra sonic device or instrument bath filled with detergent/disinfectant.
6. During chemical disinfection in the ultrasonic device or instrument bath, observe the instructions of the manufacturer regarding concentration and immersion time. Be sure to observe the full correct immersion time which does not start until the last tip has been immersed. Attention: do not exceed 45°C (risk of protein coagulation)!
7. On completion of the immersion time, rinse sonic tip/cooling adapter thoroughly with suitable water (preferably with demineralised water to avoid residues of lime). Rinse the perforation of the sonic tip with at least 10 ml of demineralised water by means of a cannula, to ensure that no residues of the detergent are left in the cooling passage.
8. Dry sonic tip (preferably with compressed air as recommended by the Commission for Hospital Hygiene and Infectious Disease Prevention of the Robert Koch Institute). When drying the inside of the instruments, please make sure to hold the compressed air gun close to the cooling passage, in order to guarantee that sufficient air flows through the sonic tip.
9. Visual examination to ensure that the instrument is clean and undamaged. If there are still visible residues of contamination, repeat the cleaning and chemical disinfecting process until no visible contamination is left.
According to the recommendations of the Commission for Hospital Hygiene and Infectious Disease Prevention of the Robert Koch Institute, reprocessing should preferably be carried out mechanically. Information on the reprocessing of instruments in compliance with the standard DIN EN ISO 17664, the present instructions for use, as amended, and more detailed information can be downloaded from our website www.kometdental.de or requested directly from the manufacturer Gebr. Brasseler.

6. Maintenance/care
Place sonic tips in a cleaning/disinfection tank (Fräsator) filled with a suitable detergent/disinfectant (e.g. Komet DC1) immediately after use in the mouth to prevent drying of residues on the instruments (protein fixation). Hold sonic tips at an angle when placing them into the liquid to ensure that all hollow spaces are reached. It is recommended to reprocess the instruments within one hour of use at the very latest. The instruments should be in the cleaning/disinfection tank (Fräsator) when transported to the site where the reprocessing is to take place.
7. Reprocessing
The tips are supplied non-sterile and have to be prepared prior to first use.
Manual pre-cleaning (prior to mechanical and manual reprocessing):
1. Remove sonic tip and cooling adapter from cleaning/disinfection tank (Fräsator) immediately before reprocessing.
2. Separate the sonic tip from the cooling adapter before the start of the cleaning process with the tip changer or

5. Recommendations for use
The sonic tips may only be used with the authorised power systems and permitted power settings. The power settings 2 (gentle) and 3 (standard) may be used. This applies to both the Komet sonic hand piece SF1LM and the KaVo SONICflex™ 2003. The power setting is regulated by means of a regulating ring. Please observe the recommendations in the instructions for use provided with the sonic hand pieces. The sonic tips are at their most effective when used with low contact pressure. Excessive contact pressure may lead to undesirable effects (see section 3 "warnings"). Given that these sonic tips may fracture without warning even when used according to instruction and their designated use, always plan two sonic tips for the operation, i.e. always have a reserve at hand so that the planned operation can still be finished even if one tip has broken.
Provision of cooling
The surgical sonic tips are cooled with at least 50 ml/min of an externally fed, sterile cooling liquid. Make sure to turn off the water supply of the dental treatment unit. It is absolutely essential to thoroughly and continuously irrigate the site at all times during the operation. This can be ensured by means of the cooling adapter SF1979. One end of the cooling adapter is screwed into the sonic hand piece and the sonic tip is screwed into the other end of the cooling adapter. The flexible tube supplying the cooling agent is placed onto the small tube of the cooling adapter. The sonic tip can be removed from the cooling adapter by means of the mounting wrench 566. Controlled supply of the cooling agent by means of a suitable surgical power system is preferable, for example with the KaVo INTRASurg 1000 for oscillating work. There are various compression systems available on the market that ensure pressure supported transport of the cooling agent:
- OMNIA system: manual supply of a sterile cooling agent by means of a compression pouch.
Article no. H2200 Manual cooling system with compression pouch (Omnia no. 30E0008.00), article no. H2204 Manual cooling system with flux control, to be attached to the hand piece (Omnia no. 32F0157). For further information, visit www.omniasrl.com.
- Aqua-Press sterile irrigation system, made by the company Ronvig Dental A/S, Gl. Vejlevej 59, DK-8721 Dagaard, tel. 0045 70 23 34 11, fax: 0045 76 74 07 98, www.ronvig.com

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