	<p>Komet Dental Gebr. Brasseler GmbH & Co. KG Trophagener Weg 25 32657 Lemgo · Germany</p> <p>Telefon+49 (0) 526 1 701-0 Telefax+49 (0) 526 1 701-289 info@kometdental.de www.kometdental.de</p>
	

Schallspitzen SFS109/SFS109F

© 10/2013, Gebr. Brasseler GmbH & Co. KG
310454 - Rev 10/2013

DE EN FR ES IT

rostfreier Stahl - stainless steel - acier inoxydable - acero inoxidable - acero inossidabile


Abb./Fig. 1

Abb./Fig. 2

DE

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EN

Surgical sonic tips for the creation of a bone window as part of an open sinus lift, developed by Dr. Ivo Agabiti

1. Indication

The sonic tips SFS109 and SFS109F are designed for preparing a bone window as part of a surgical method for carrying out an open sinus lift. They mainly serve to definitely detach the bone plate, to smooth the bone margins, to prepare the bone plate or to complement the work of the rotary instruments.

The correct use of the sonic tips reduces the risk of an inadvertent perforation of the mucosa in the maxillary sinus, which would cause severe problems to the surgeon and the patient. The sonic tips may only be used by dentists and medical personnel that are sufficiently trained in the relevant fields of dental surgery.

2. Contra-indications

All general contraindications for dental surgery must be taken into account. All commonly known contraindications in oral maxillo facial surgery (e.g. implantation, jaw surgery) also have to be considered.

3. Warnings

The useful life of the sonic tips is limited. They therefore have to be used with the utmost care.

Attention:

- **Risk of fracture!** Spontaneous and premature fracture must be expected. The sonic tips may fracture unexpectedly and without warning. In order to reduce the risk of fracture, make sure to stick to the correct power setting and observe all recommendations for use. Never allow the instruments to jam. Do not use instruments for leverage. Avoid excessive contact pressure. Use sonic tips with the utmost care as even minor mechanical damage can provoke the fracture of the instruments.

- **Risk of ingestion!** Fracture of the sonic tip leads to the risk of the tip being accidentally swallowed. Take all appropriate steps to prevent this.

- **Premature wear!** Work on hard dental substance can cause wear of the diamond coating. Likewise, excessive contact pressure can lead to damage to the working part of the diamond coating. Inspect instruments after every use for signs of wear or damage. Exchange if necessary.

- **Excessive heat generation!** Excessive contact pressure and/or the use of worn/blunt instruments can cause hazardous heat generation which increases the risk of thermal necrosis. Excessive contact pressure therefore has to be avoided. Damaged sonic tips (changed shape, corroded surfaces, void spots or defect thread etc.) have to be replaced. The sonic tips are medical devices and as such subject to the applicable legal provisions of the relevant country. In compliance with such regulations, the sonic tips may only be used by properly qualified operators within the scope of the indicated application and according to the valid industrial safety regulations and observing the valid rules for accident prevention. The provisions further stipulate that all equipment used must be free of defects and suitable for the intended application. Furthermore, the patient and third persons must be protected from hazards. Possible contamination caused by the product must be ruled out.

Safety control:

Fracture may occur due to the tips being subjected to continuous or sudden stress (for example if the tips are dropped on the floor or through a mechanical change of the original shape). Check the tips prior to each use to make sure they are safe to use. To do this, apply slight pressure with your thumb and index finger. For additional safety, load the tips mechanically with approx. 10 N (1kg), without function. Always wear protective gloves when checking, inserting or removing the sonic tips to avoid risks of injury and infection. We recommend that the tips be exchanged regularly.

4. Appropriate use

The tips are supplied non-sterile and have to be prepared prior to first use. The Komet® sonic tips SFS109 and SFS109F can be used in the following power systems:

- in the Komet sonic hand piece SF11M

- in the SONICflex® hand piece made by the co. KaVo Series 2003N/L/X/LX

The above listed sonic hand pieces are provided with a MULTIFlex™ coupling. The tip changer with torque SF1975 is used for inserting and removing the sonic tips or cooling adapter (fig. 1). Screw the chosen sonic tip into the hand piece with the mentioned tip changer, turning the tool in a clockwise direction. To remove the tip, place tip changer on the tip and turn anticlockwise. Note: Please observe the directions in the instructions supplied with your sonic hand piece.

5. Recommendations for use

The sonic tips may only be used with the authorised power systems and permitted power settings. The power setting 1-3 may be used. This applies to both the Komet sonic hand piece SF11M and the KaVo SONICflex® 2003N/L/X/LX. In both power systems, the power setting is regulated by means of a regulating ring. Please observe the recommendations in the instructions for use provided with the sonic hand pieces. The sonic tips are at their most effective when used with low contact pressure. Excessive contact pressure may lead to undesirable effects (see section 3 "warnings"). Given that these sonic tips may fracture without warning even when used according to instruction and their designated use, always plan two sonic tips for the operation, i.e. always have a reserve at hand so that the planned operation can still be finished if one tip has broken.

Provision of cooling

The surgical sonic tips are cooled with at least 50 ml/min of an externally fed, sterile cooling liquid. Make sure to turn off the water supply of the dental treatment unit. It is absolutely essential to thoroughly and continuously irrigate the site at all times during the operation. This can be ensured by means of the cooling adapter SF1979. One end of the cooling adapter is screwed into the sonic hand piece and the sonic tip is screwed into the other end of the cooling adapter. The flexible tube supplying the cooling agent is placed onto the small tube of the cooling adapter. The sonic tip can be removed from the cooling adapter by means of the mounting wrench 566. Controlled supply of the cooling agent by means of a suitable surgical power system is preferable, for example with the KaVo INTRASurg 1000 for oscillating work. There are various compression systems available on the market that ensure pressure supported transport of the cooling agent:

- OMNIA system: manual supply of a sterile cooling agent by means of a compression pouch.

Article no. H2200 Manual cooling system with compression pouch (Omnia no. 30E0008.00), article no. H2204 Manual cooling system with flux control, to be attached to the hand piece (Omnia no. 32F0157). For further information, visit www.omniasrl.com.

- Aqua-Press sterile irrigation system, made by the company Ronvig Dental A/S, Gl. Vejlevej 59, DK-8721 Daugaard, tel. 0045 70 23 34 11, fax: 0045 76 74 07 98, www.ronvig.com.

6. Maintenance/care

Place sonic tips and cooling adapter in a cleaning/disinfection tank (Fräsator) filled with a suitable detergent/disinfectant (e.g. Komet DC1™) immediately after use in the mouth to prevent drying of residues on the instruments (protein fixation). Hold sonic tips and cooling adapter at an angle when placing them into the liquid to ensure that all hollow spaces are reached. It is recommended to reprocess the instruments within one hour of use at the very latest. The instruments should be in the cleaning/disinfection tank (Fräsator) when transported to the site where the reprocessing is to take place.

7. Reprocessing

The tips are supplied non-sterile and have to be prepared prior to first use.

Manual pre-cleaning (prior to mechanical and manual reprocessing):

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