Komet Dental
 Gebr. Brasseler GmbH & Co. KG
 Trophagener Weg 25
 32657 Lemgo - Germany
 Telefon+49 (0)5261 701- 0
 Telefax+49 (0)5261 701- 289
 info@brasseler.de
 www.brasseler.de

AlphaKite
 © Gebr. Brasseler GmbH & Co. KG
3 10467 - Rev 2018-12

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AlphaKite System for the mechanical preparation of root canals

1. Indication

The instruments are intended for the mechanical removal of diseased tissue and root filling materials within the course of a root canal
 treatment.

2. Contraindication

The contra-indications are (e. g.):

- pronounced abrupt curves
- severely calcified root canals
- preparation of milk teeth

3. Proper use

The general recommendations for a mechanical root canal preparation with NiTi files apply.
 - Creation of a straight canal access
 - Adherence to the Crown-Down technique or Double-Flare technique
 - Gentle work in intermittent mode, without applying pressure
 - Only leave instrument in the canal for a short period of time (5-8 seconds at most)
 - To be used in lubricated canal only
 - The instrument must be in rotary mode when inserted into the canal
 - Lubricant used must contain EDTA
 - Always rinse canal thoroughly (NaOCl, citric acid, chlorhexidine)
 - Check patency by inserting a hand instrument, for example a K-file of size 010 or 015
 - Monitor and record carefully how many times each instrument has been used (5-8 canals at most, when used in extremely curved canals,
 discard after one use, if necessary)
 - Use in a torque-limited endodontic crown-angle, e.g. EndoPilot mobil (Komet), W&H ENTRAN or an endodontic motor, e.g. E-Drive (Komet),
 KaVo ENDOadvance or Sironiti/Sironiti Air+
 - The use of a suitable torque-limited power system guarantees sufficient safety from torsional fracture. Do not exceed the file's load limit
 or its maximum fracture moment. Check instruments visually prior to use. Discard damaged or bent instruments. Prior to first use on a
 patient, the user should acquaint himself with the system by practicing on extracted teeth.
 - A straight access cavity of sufficient size should be created for a root canal preparation. The canal orifice is prepared up to the curve with
 the access reamer AK10 (L15 or L19),035 at a torque of 2,0 Nm and an optimum speed of \varnothing_{max} 500 rpm (max. \varnothing_{max} 800 rpm). This
 instrument can be identified by a green spiral. Never insert this file as far as the curve of the canal (risk of fracture or incorrect preparation).
 The preparation can for example be carried out according to the Crown-Down technique. Subsequent reaming is carried out with a file with
 a taper of 08 and/ or 06, followed by a file with a taper of 04. The preparation in apical direction is carried out with an 02 tapered file. The
 apical area is then enlarged by one size with a file with a taper of 02. The apex can be further enlarged with an 04 file if this is required by
 the filling technique used.
 - Instrument sequences
 When used according the crown-down technique:
 Narrow canals, yellow AK files 020: AK10L19.035 green, AK08L25.020 yellow, AK06L25.020 yellow, AK04L25.020 yellow, AK02L25.020
 yellow, AK02L25.025 red and optionally once again AK04L25.020 yellow to enlarge to taper 04.

Average canals, red AK files 025: AK10L19.035 green, AK08-06-04-02L25.025 red, AK02L25.030 blue, optionally: AK04L25.025 red.
 Wide canals, blue AK files 030: AK10L19.035 green, AK08-06-04-02L25.030 blue, AK02L25.035 green, optionally: AK04L25.030 blue
 The files can be used for other techniques as well, for example the Double-Flare technique. A decisive factor for the success of any
 preparation is the constant probing with a manual file of size 010 or 015, thus not only creating a smooth path for the rotating instruments
 but also preventing possible blockages during preparation. The canal should also be rinsed with an adequate quantity of rinsing solution
 after each change of file.

4. Reprocessing

Remove silicone stopper from the file prior to reprocessing. Disinfect the AlphaKite files with cleaning and disinfecting agents for rotary
 instruments (e.g. Komet DC1). For recommendations for use (immersion time, concentration, suitability) of the cleaning and disinfecting
 agents please refer to the instructions provided by the manufacturer of the agents used. The instruments can be reprocessed in the thermal
 disinfecter, provided that the agent used is suitable for rotary instruments (mildly alkaline detergent). Cleaning and disinfecting agents
 must be rinsed off thoroughly with water and the instruments must be dried carefully (preferably with compressed air according to the
 recommendations of the Commission for Hospital Hygiene and Infectious Disease Prevention of the Robert Koch Institute). Do not store
 instruments for a longer period in wet or humid condition. Take great care to remove all traces of NaOCl and limit contact time to a few
 minutes only. The drills must not come in contact with each other during the cleaning in the ultrasonic bath. Carry out a visual control of the
 instruments. Damaged or blunt instruments must be rejected and their use discontinued.

Sterilization:

The sterilization of the instruments is carried out in the autoclave at 134°C. Observe the instructions provided by the manufacturer. After
 sterilization, check instruments for surface defects (corrosion). Do not reuse corroded instruments. The operator of medical products is
 responsible for ensuring that proper treatment is carried out by qualified personnel, using suitable materials and suited equipment, as
 recommended by the Commission for Hospital Hygiene and Infectious Disease Prevention of the Robert Koch Institute. Information on
 proper reprocessing of instruments according to DIN EN ISO 17664, the present instructions for use, as amended, and more detailed
 information can be downloaded from our web site www.kometdental.de or requested from the manufacturer Gebr. Brasseler.

5. Storage

Store packed instruments away from UV radiation and heat in a clean and dry environment. Do not store in the same room as solvents or
 chemicals.

6. Disposal

The instruments have to be disposed of in tight, non-breakable and non-pierceable containers (protection from contamination).

7. Safety and liability

See also "Instructions for use and safety recommendations" in the current Komet Dental catalogue. The user is responsible for checking
 the products prior to use as to whether they are suited for the intended purpose. The user is responsible for the appropriate use of the
 instruments. In case of contributory negligence by the user, Gebr. Brasseler partially or totally declines liability for all resulting damages,
 particularly if these are due to nonobservance of our recommendations for use or warnings as well as inadvertent misuse by the user.

8. Procedure in case of serious adverse events

In case of serious adverse events, please proceed in line with the legal regulations in your country.