	Komet Dental
	Gebr. Brasseler GmbH & Co. KG
	Trophagener Weg 25
	32657 Lemgo - Germany
	Telefon+49 (0)5261 701-0 Telefax+49 (0)5261 701-289 info@kometdental.de www.kometdental.de

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Sonic tips for retrograde endodontic treatments

1. Indication
 These sonic tips are designed for use during apicoectomies, when cysts or apical periodontitis can no longer be treated with traditional treatment methods. In the front an lateral region, sonic tips for retrograde endodontic treatments permit the retrograde axial preparation of the apical part of the root canal while preserving a maximum of bone substance. The tips are used in three chronological steps, each instrument stage covering a special indication.
 SF56 and SF57 are used for localizing the canal entrance. Thanks to their fine tips, they can create an isthmus, if required. SF16 and SF17 are used during the actual preparation of the root canal. Due to their specially adapted curvatures extending in different directions, they can be used in every quadrant.
 Thanks to their T-shaped heads, the SF20 and SF21 can create retentive undercuts by performing circular movements in order to prevent loss of the canal filling and to provide the new apex with a bacteria proof seal.
 The special shape of the SF55 tip permits work on roots with apexes strongly curved in an oral direction.

2. Contraindication
 When a patient is considered for treatment, all general contraindications for dental surgery must be taken into account. All commonly known contraindications in oral maxillo facial surgery (e.g. implantation, jaw surgery) also have to be considered.
3. Warnings
 The useful life of the sonic tips is limited. They therefore have to be used with the utmost care.
 Attention:
 - Risk of fracture! Spontaneous and premature fracture must be expected. The sonic tips may fracture unexpectedly and without warning. In order to reduce the risk of fracture, make sure to stick to the correct power setting (see instructions for use) and observe all recommendations for use. Never allow the instruments to jam. Do not use instruments for leverage. Avoid excessive contact pressure. Use sonic tips with the utmost care as even minor mechanical damage can provoke the fracture of the instruments.
 - Risk of ingestion! Fracture of the sonic tip leads to the risk of the tip being accidentally swallowed. Take all appropriate steps to prevent this (for example by using a rubber dam).
 - Premature wear! Work on hard dental substance can cause premature blunting of the blades and/or wear of the diamond coating. Likewise, excessive contact pressure can lead to damage to the working part, the blades or the diamond coating. Inspect instruments after every use for signs of wear or damage. Exchange if necessary.
 - Excessive heat generation! Excessive contact pressure and/or worn or blunt instruments can cause hazardous heat generation which increases the risk of thermal necrosis (damage to tissue). Excessive contact pressure therefore has to be avoided. Damaged sonic tips (changed shape, corroded surfaces, void spots or defect thread etc.) have to be replaced.
 Incorrect use:
 - In scalers made by ca. W&H (e. g. series Proxeo® ZA-55/L/LM/M/LS and ST ZE-55RM/BC, series Synea® ZA-55/L/LM/M or Alegre® ST ZE-55RM/BC)
 - In the SIRONIA L provided by co. Dentsply Sirona

The tip changer for SF1975 is used for inserting and removing the sonic tips. Screw the chosen sonic tip into the hand piece with the mentioned tip changer, turning the tool in clockwise direction. To remove the tip, place tip changer on the tip and turn anticlockwise. Note: Please observe the directions in the instructions supplied with your sonic hand piece.
5. Proper use
 Prior to the use of the sonic tips, the following surgical procedure is recommended:
 - Siting of the mucosa according to the habitual surgical technique.
 - Creation of a mucoperiosteal flap to expose the apical operative site.
 - Work on the bone in the region around the apex. It is recommended to prepare a bone window with a repositionable lid (approx. 6 mm lateral length), for example with the Komet sonic tip SFS100/101/102.
 - After that, the root apex is removed in a right angle to the tooth axis, for example with the tungsten carbide bone cutter Komet H254E.
 - To enable the retrograde preparation, remove any diseased tissue by means of a sharp spoon curette and stop the resulting bleeding. This is followed by the actual preparation of the root canal according to the steps described under paragraph 1 "indication". Start up sonic tip prior to application to the tooth, providing plenty of spray cooling. For the treatment, the power setting 1 is authorised for the SF1LM and SF1LS and the approved power systems made by W&H and Dentsply Sirona, and power setting 2 for the SONICflex™ 2003 (this corresponds to 100% capacity of the SONICflex™ 2000).
 The use of power setting 2 is not permitted for the SF1LM, SF1LS as well as the approved power systems made by W&H and Dentsply Sirona, and power level 3 is prohibited with any of these power systems when working with sonic tips (see warnings).
 Apply low contact pressure only, excessive contact pressure reduces the vibration of the sonic tip and consequently its effectiveness. During the treatment, the removed material can be evacuated by the cooling liquid and by performing gentle up-and-down movements with the instrument (Attention: this does not apply to the tips SF20 and SF21).
 After that, the root canal is filled with a material suitable for retrograde fillings.
 Finally, the wound is closed by pressing the mucoperiosteal flap down for 1 minute. Larger defects might have to be filled with bone material.
 Cooling:
 The surgical sonic tips for retrograde endodontic treatments are cooled with at least 50 ml/min of an externally fed, medical cooling liquid. Make sure to turn off the water supply of the dental treatment unit. It is absolutely essential to thoroughly and continuously irrigate the site at all times during the operation. This can be ensured by means of the cooling adapter SF1979. One end of the cooling adapter is screwed into the sonic hand piece and the sonic tip is screwed into the other end of the cooling adapter. The flexible tube supplying the cooling agent is placed onto the small tube of the cooling adapter (fig. 3). The sonic tip can be removed from the cooling adapter by means of the mounting wrench 566. Controlled supply of the cooling agent by means of a suitable surgical power system is preferable, for example with the KaVo INTRASung 1000 for oscillating work.
 There are various manual systems available on the market that ensure the external supply of the cooling agent, e.g.:
 System Aqua-Press sterile irrigation system made by the company Rovig Dental A/S, Gl. Vejlevej 59, DK-8721 Dagaard. Tel.: 0045-70 23 34 11, Fax: 0045-76 74 07 98; www.rovig.com.
 The company Implants Deutschland sells a range of suitable manual cooling systems.
 Germany: (0180) 5 90 75 90, free fax: (00800) 05 05 51 23
 Austria: (0049 180) 5 90 75 90, free fax: (00800) 05 05 51 23
 To order online, contact www.implants.eu. Choice of articles:
 H2200 Manual cooling system with compression pouch (Omnia-Nr. 30E0008.00)
 H2204 Manual cooling system with flux control (Omnia-Nr. 32FO157)

6. Maintenance
 Place sonic tip and/or cooling adapter in a cleaning/disinfection tank (Fräsator) filled with a suitable detergent/disinfectant (e.g. Komet DC1) immediately after use in the mouth to prevent drying of residues on the instruments (protein fixation). Hold sonic tip and/or cooling adapter at an angle when placing them into the liquid to ensure that all hollow spaces are reached. It is recommended to reprocess the instruments within one hour of use at the very latest. The instruments should be in the cleaning/disinfection tank (Fräsator) when transported to the site where the reprocessing is to take place.

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7. Reprocessing

Sonic tips supplied non-sterile have to be prepared prior to first use.
 Manual pre-cleaning (prior to manual and mechanical reprocessing):
 1. Remove sonic tip and cooling adapter from cleaning/disinfection tank (Fräsator) immediately before reprocessing.
 2. Separate the sonic tip from the cooling adapter before the start of the cleaning process with the tip changer to ensure proper reprocessing of both components.
 3. Remove stubborn contamination with a nylon brush under running water, turning the sonic tip and the rinse adapter constantly.
 4. Check the patency of the cooling passage by guiding the nozzle cleaner (accessory of the sonic hand piece) through the passage. Replace instrument if the passage is obstructed.
 Mechanical reprocessing (validated):
 1. Rinse sonic tip/cooling adapter thoroughly under running water to prevent any residues of the detergent/disinfectant from getting into the washer/disinfecter.
 2. Unscrew the nozzle on the injector bar of the washer/disinfecter. Then screw the external thread of the rinse adapter SF1978/L in the silicon holders on the injector bar. Now screw the sonic tip in the internal thread at the top of the rinse adapter (Fig. 4 depicts the rinse adapter SF1978/L). Observe the directions in the instructions for use provided with the rinse adapter SF1978/L. When using the rinse adapter SF1979, proceed as follows: Unscrew the nozzle on the injector bar of the washer/disinfecter. Then screw the large external thread of the rinse adapter SF1977 in the silicon holders on the injector bar. Now screw the cooling adapter in the small external thread at the top of the rinse adapter. Observe the directions in the instructions for use provided with the rinse adapter.
 3. Put chemical detergent into the washer/disinfecter, following the indications on the label and the instructions of the manufacturer of the washer/disinfecter.
 4. Start the Vario TD programme including thermal disinfection. Thermal disinfection takes place allowing for the Av value and observing national provisions (pREN/ISO 15883).
 5. On completion of the cycle remove sonic tip/cooling adapter from the washer/disinfecter and dry (preferably with compressed air according to the recommendations of the Commission for Hospital Hygiene and Infectious Disease Prevention of the Robert Koch Institute). When drying the inside of the instruments, please make sure to hold the compressed air gun close to the cooling passage, in order to guarantee that sufficient air flows through the sonic tip or rinse adapter.
 6. Visual examination to ensure that the instrument is clean and undamaged. If after mechanical reprocessing there are still visible residues of contamination, repeat the cleaning and disinfecting process until no visible contamination is left.
 Manual reprocessing (alternative):
 1. Place sonic tip/cooling adapter onto a 10 ml syringe and rinse with detergent/disinfectant. In case of visible contamination emerging from the perforations, rinse again with detergent/disinfectant. Remove stubborn contamination with a nylon brush under running water, turning the sonic tip/cooling adapter constantly.
 2. Remove stubborn contamination with a nylon brush under running water, turning the sonic tip/cooling adapter constantly.
 3. Rinse sonic tip/cooling adapter thoroughly with running water.
 4. Visual examination to ensure that the instrument is clean. If there are still visible residues of contamination, repeat the cleaning process until no visible contamination is left.
 5. Place sonic tip/cooling adapter in a suitable container (e.g. Komet 9952) into the ultra sonic device or instrument bath filled with detergent/disinfectant.
 6. During chemical disinfection in the ultrasonic device or instrument bath, observe the instructions of the manufacturer regarding concentration and immersion time. Be sure to observe the full correct immersion time which does not start until the last tip has been immersed. Attention: do not exceed 45°C (risk of protein coagulation!)
 7. On completion of the immersion time, rinse sonic tip/cooling adapter thoroughly with suitable water (preferably with demineralised water to avoid residues of lime). Rinse the perforation of the sonic tip with at least 10 ml of demineralised water by means of a cannula, to ensure that no residues of the detergent are left in the cooling passage.
 8. Dry sonic tip (preferably with compressed air according to the recommendations of the Commission for Hospital Hygiene and Infectious Disease Prevention of the Robert Koch Institute). When drying the inside of the instruments, please make sure to hold the compressed air gun close to the cooling passage, in order to guarantee that sufficient air flows through the sonic tip or rinse adapter.
 9. Visual examination to ensure that the instrument is clean and undamaged. If there are still visible residues of contamination, repeat the cleaning and chemical disinfecting process until no visible contamination is left.
 According to the recommendations of the Commission for Hospital Hygiene and Infectious Disease Prevention of the Robert Koch Institute, the further reprocessing should be carried out mechanically. Information on the reprocessing of instruments in compliance with the standard DIN EN ISO 17664, the present instructions for use and more detailed information, as amended, can be downloaded from our website www.kometdental.de or requested directly from the manufacturer Gebr. Brasseler.
 Sterilisation in the autoclave:
 Make sure that the packaging is suitable for sonic tips classed as risk group critical B, their accessories and the chosen method of sterilisation.
 Single pack: The packaging must be large enough to ensure that there is no pressure on the seal. In the sterilisation container: Place sonic tips or rinse adapter into a suitable sterilisation container (e.g. 9952). For sterilisation, the container has to be sealed in suitable packaging. For this sterilisation, the following instructions have to be observed:
 Steam sterilisation using a vacuum process at 134°C in a device that complies with the provisions of DIN EN 13060; with validated processes:
 -fractionated pre-vacuum (type B)
 -sterilisation temperature: 134°C
 -hold time: at least 5 minutes (full cycle)
 -drying time: at least 10 minutes

In order to prevent staining and corrosion, the steam must be free of particles. The recommended limits for particle contents in feed water and condensed steam are defined by standard DIN EN 13060. Make sure not to exceed the maximum capacity of the sterilizer when sterilizing several instruments. Follow the instructions of the device manufacturer.
 Universally valid notes:
 Observe the legal provisions regarding the reprocessing of medical products valid in your country (e.g. www.rki.de). The manufacturer confirms that the above detailed reprocessing methods are suitable for preparing the above named instrument group to enable their reuse. The person carrying out the reprocessing is responsible for ensuring that the applied method is carried out with appropriate equipment, materials and trained personnel at the reprocessing site and that it actually achieves the desired result. To guarantee this, routine controls of the validated mechanical preparation methods are necessary. Any deviation from the above detailed process must be carefully checked by the operator to ensure effectiveness and to avoid possible adverse consequences.
 Accessories:
 Komet tip changer SF1975
 Komet cooling adapter SF1979 + mounting wrench 566 (in the set 4602)
 For preparation in the Miele washer/disinfecter: Komet rinse adapter SF1978 for Komet sonic tips and Komet rinsing adapter SF1977 for Komet cooling adapter SF1979

8. Storage
 Store packed instruments away from UV radiation and heat in a clean and dry environment. Do not store in the same room as solvents or chemicals.

9. Disposal
 The instruments have to be disposed of in tight, non-breakable and non-pierceable containers (protection from contamination).

10. Safety and liability
 See also "Instructions for use and safety recommendations" in the current Komet Dental catalogue. The user is responsible for checking the products prior to use as to whether they are suited for the intended purpose. The user is responsible for the appropriate use of the instruments. In case of contributory negligence by the user,

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© Gebr. Brasseler GmbH & Co. KG