**Ceradrill K210**  
 © Gebr. Brasseler GmbH & Co. KG  
 310569 - Rev 2016-07


Implantologie, Implantology, Implantologie, Implantología, Implantologia

**STERILE/O** Sterilisation durch Ethylenoxid (Sterile Instrumente sind auf dem Verpackungsetikett entsprechend gekennzeichnet), Sterilisation with ethylene oxide (Sterile instruments are marked as such on the label of the packaging), Stérilisation avec oxyde d'éthylène (Les instruments stériles sont marqués comme tels sur l'étiquette de l'emballage), Esterilización con óxido de etileno (Los instrumentos estériles son marcados como tales sobre la etiqueta del embalaje), Sterilizzazione con ossido di etilene (Gli strumenti sterili sono marcati come tali sull'etichetta dell'imballaggio)

Verwendbar bis (nur zutreffend bei sterilen Instrumenten), Use by (only applies to sterile instruments), A utiliser avant le (ne s'applique qu'aux instruments stériles), Fecha de caducidad (se aplica solamente a instrumentos estériles), Data di scadenza (riguarda soltanto strumenti sterili)

Bei beschädigter Verpackung nicht verwenden (nur zutreffend bei sterilen Instrumenten), Do not use in case of damaged packaging (only applies to sterile instruments), Ne pas utiliser si l'emballage est endommagé (ne s'applique qu'aux instruments stériles), No utilizar si el embalaje está dañado (se aplica solamente a instrumentos estériles), Non utilizzare se l'imballo è danneggiato (riguarda soltanto strumenti sterili)

○<sub>min</sub> 800 - 1.000 min<sup>-1</sup>/rpm ○<sub>max</sub> 6.000 min<sup>-1</sup>/rpm

**DE**  
 Ceradrill - Implantologievorböhrer aus Keramik

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**EN**  
 Ceradrill - Ceramic pilot drill for use in implantology

- System overview**  
 See table 1.
- Indication**  
 - For the preparation of the implant site for the most commonly used implants. If a template is used make sure that the instrument does not jam in order to prevent possible breakage of the instrument.  
 - Determination of position, direction and depth of the implant.
- Clinical sequence**  
 - Select the required size of the pilot drill and determine the drilling depth.  
 - Determine the size and the exact position of the implant with the help of radiographs (bitewing radiograph and orthopantomogramme).  
 - In order to establish the exact position and the depth of the borehole, a computer tomographic examination in addition to the bitewing radiograph and the orthopantomogramme is recommended.  
 - Examine the adjacent structures  
 - In order to prevent any danger to adjacent structures, carefully check the area surrounding the operation site.  
 - Administer a local anaesthesia to the site where the instrument is to operate.  
 - Cut the mucosa and periost to the bone and expose the operation site.  
 - Mark the exact position with the help of a round bur, for example K 160A.205.023.  
 - Apply the pilot bur at the chosen position.  
 - Drilling should be carried out in intermittent mode using constant exterior cooling with sterile physiological saline solution to prevent thermal damage to the bone tissue and to ensure proper chip removal. Make sure that the drill neither jams nor gets stuck during use (increased risk of fracture). Prepare to the desired depth at a speed of 800 - 1.000 rpm (green contra-angle, reduction at least 10:1), exerting low pressure. The indicated speed should be observed in order to avoid fracture of the instrument. Do not exceed a torque of 30 Ncm.
- Further treatment**  
 - The subsequent instrumentation is used to enlarge the implant site to the defined implant diameter.
- Recommendations for maintenance**  
 The instruments are supplied non-sterile and have to be prepared prior to first use.  
 - Do not allow residues from the operation (blood, secretions, tissue) to dry.  
 - Immerse the drills in disinfecting solution immediately after the operation.  
 - Disinfect the drills with disinfecting and cleaning agents for rotary instruments (for example Komet DC1).

For recommendations for use (immersion time, concentration, suitability) of the disinfecting and cleaning agents please refer to the manufacturers' instructions. Only use brushes with metal-free bristles for pre-cleaning, as these do not leave discolorations caused by abrasion on the instruments. The instruments can be reprocessed in the thermo disinfectant, provided that the agent used is suitable for rotary instruments (mildly alkaline detergent). Disinfecting and cleaning agents must be rinsed off thoroughly with water and the instruments must be dried carefully (preferably with compressed air according to the recommendations of the Commission for Hospital Hygiene and Infectious Disease Prevention of the Robert Koch Institute). Do not store instruments in wet or damp condition for a longer period of time. The drills should not be allowed to come in contact with each other during the cleaning in the ultrasonic bath. Examine the cleaned instruments visually. Damaged or blunt instruments must be rejected and their use discontinued. The same applies to instruments that were subjected to sudden stress (for example by getting dropped off the table or similar). Do not continue using these instruments, as there is an increased risk of breakage! Sterilization is carried out in the autoclave (134°C) following the usual methods. Make sure to observe the instructions provided by the manufacturer of the devices. After sterilization, the instruments have to be examined for surface damages (cracks). Discontinue the use of damaged instruments (increased risk of breakage). The operator of medical products is responsible for seeing that proper treatment is carried out by qualified personnel, using the appropriate materials and suited equipment, as recommended by the Commission for Hospital Hygiene and Infectious Disease Prevention of the Robert Koch Institute. Work instructions with regard to proper reprocessing of instruments according to DIN EN ISO 17664 can be downloaded from our web site [www.kometdental.de](http://www.kometdental.de) or requested from the manufacturer Gebr. Brasseler.

**6. Storage**  
 Store sterile packed instruments in their original packaging, away from light and heat in a clean and dry environment. Do not store in the same room as solvents or chemicals.

**7. Safety and liability**  
 Due to the cutting of hard bone material and tooth substance premature blunting of the blades may occur. Therefore, the drills have to be checked for blunt blades or other damages after each use. Exchange if necessary. Be sure to observe the indicated speed in order to avoid instrument breakage. Attention: Risk of injury caused by the sharp blades of the drill! Risk of injury due to jamming and slipping of the drill! Prior to application, the user undertakes to check whether the product is suited for the intended use. The user is responsible for the application of the drill. In case of contributory negligence by the user, Gebr. Brasseler partially or totally declines liability for all resulting damages, particularly if these are due to non-observance of our recommendations for use or warnings as well as inadvertent misuse by the user.

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