Kavitätenpräparation, Cavity preparation, Préparation cavitaire, Preparación de cavidades, Preparazione di cavità

Verpackungseinheit, Packing unit, Unité d'emballage, Unidad de embalaje, Unità d'imballo

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
**Instructions for use**  
Please read carefully and retain for future reference.

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**Sonic tips Micro/Bevel**

- 1. Indication**  
Micro/Bevel sonic tips are used for defect orientated preparations of micro lesions without damaging the adjacent teeth. The sonic tips come in special shapes that can be used for beveling tiny lesions.
- 2. Contraindication**  
There are no known contra-indications related to the instruments.
- 3. Warnings**  
The useful life of the sonic tips is limited. They therefore have to be used with the utmost care  
Attention:  
- Risk of fracture! The sonic tips may fracture unexpectedly and without warning. In order to reduce the risk of fracture, make sure to stick to the correct power setting (see instructions for use) and observe all recommendations for use. Never allow the instruments to jam. Do not use instruments for leverage. Avoid excessive contact pressure. Use sonic tips with the utmost care as even minor mechanical damage can provoke the fracture of the instruments.  
- Risk of ingestion! Fracture of the sonic tip leads to the risk of the tip being accidentally swallowed. Take all appropriate steps to prevent this.  
- Premature wear! Work on hard dental substance can cause premature wear of the diamond coating. Likewise, excessive contact pressure can lead to damage to the working part. Therefore, inspect instruments after every use for signs of wear or damage. Exchange if necessary.  
- Excessive heat generation! Excessive contact pressure and/or worn or blunt instruments can cause hazardous heat generation which increases the risk of thermal necrosis (damage to tissue). Excessive contact pressure therefore has to be avoided. Damaged sonic tips (changed shape, corroded surfaces, void spots or defect thread etc.) have to be replaced. Incorrect use can lead to the above detailed undesirable effects as well as to injury to the user, the patients or third persons.  
Increased oscillation:  
Please note that, depending on the chosen power level, the oscillation of sonic tips intended for this indication might exceed 200 µm. The increased oscillations do not affect the safety of the tips, provided that the power settings indicated in paragraph 5 (Recommendations for use) are observed and the threshold value regarding the air supply of the authorized power system used are observed.  
Safety control:  
Fracture may occur due to the tips being subjected to continuous or sudden stress (for example if the tips are dropped on the floor or through a mechanical change of the original shape). Check the tips prior to each use to make sure they are safe to use. To do this, apply slight pressure with your thumb and index finger (fig. 1). For additional safety, load the tips mechanically with approx. 10 N (1 kg), without function (fig. 2). Always wear protective gloves when checking, inserting or removing the sonic tips to avoid risks of injury and infection. We recommend that the tips be exchanged regularly.
- 4. Start-up**  
These sonic tips that are supplied non-sterile have to be prepared prior to first use.  
The sonic tips can be used in the following power systems:  
- in the Komet sonic hand pieces SF1LM and SF1LS  
- in the SONICflex™ Komet hand piece of the Fa. KaVo (Series 2000N/L/LX or series 2003N/L/X/LX)  
- in scalers made by ca. W&H (Series Synea® ZA-55/L/LM/M or Alegria® ST ZE-55RM/BC) (Series Proxo® ZA-55/L/LM/M or Proxo® ST ZE-55RM/BC)  
- in the SIROAIR I provided by co. Dentsply Sirona  
The required tip is screwed into the hand-piece with the tip changer SF1975 which is supplied with the hand piece or alternatively with the SF4887 in a clockwise direction until a clicking sound can be heard. When using the tip changer, make sure that the tip faces inwards, i.e. towards the semi-radius, during insertion. To remove the tip, place the tip changer onto the tip and remove turning in an anti-clockwise direction.  
Note: Please observe the directions in the instructions supplied with your sonic hand piece.
- 5. Proper use**  
All Micro/Bevel sonic tips are coated on one side only to prevent damage to the adjacent teeth. The sonic tips SF30M (mesial) and SF30D (distal) in the sizes 016 and 024 are intended for the preparation of tiny lesions. The tips are suitable for both occlusal and lateral access to the cavity. The further treatment of the prepared tiny lesions is carried out with the instruments SF58M and SF58D in order to bevel the margins with precision to an angle of 45°. Always work with low contact pressure, excessive contact pressure reduces the vibration of the sonic tip and consequently its effectiveness.  
Start sonic tip with plenty of water spray before applying the tip to the tooth. These sonic tips may be used with the above mentioned sonic hand-pieces with the power level 1. For safety reasons, do not use power levels 2 and 3 when working with these tips. This does not apply to the SONICflex™ power systems listed under point 4, where the power settings 1 and 2 are permitted (i.e. the SONICflex™ 2000 works to 75% or 100% capacity). Non-compliance with these instructions can lead to undesirable effects as mentioned in the chapter Warnings.
- 6. Maintenance**  
Place sonic tips in a cleaning/disinfection tank filled with a suitable detergent/disinfectant (e.g. Komet DC1) immediately after use in the mouth to prevent drying of residues on the instruments (protein fixation). Hold sonic tips at an angle when placing them into the liquid to ensure that all hollow spaces are reached. It is recommended to reprocess the instruments within one hour of use at the very latest. The instruments should be in the cleaning/disinfection tank when transported to the site where the reprocessing is to take place.
- 7. Reprocessing**  
The sonic tips are supplied non-sterile and have to be processed prior to first use. Manual pre-cleaning (prior to mechanical and manual reprocessing):  
1. Remove sonic tip from the cleaning /disinfection tank immediately before reprocessing.  
2. Remove surface contamination thoroughly under running water. Completely remove stubborn contamination with a nylon brush below water level, turning the tip constantly. Rinse again under running water.  
3. Check patency of the cooling passage by guiding the nozzle cleaner (accessory of the sonic and piece) through the passage. Replace instrument if the passage is obstructed.  
Mechanical reprocessing (validated):  
1. Rinse sonic tip thoroughly under running water to prevent any residues of the detergent/disinfectant from getting into the washer/disinfectant.  
2. Unscrew the nozzle on the injector bar of the washer/disinfectant or pull it out. Then screw or push the external thread of the rinse adapter SF1978/L in the silicon holders on the injector bar. Now screw the sonic tip in the internal thread at the top of the rinse adapter

Fig. 3 shows the rinse adapter SF1978/L. Observe the directions in the instructions for use provided with the SF1978/L.  
3. Put chemical detergent into the washer/disinfectant, following the indications on the label and the instructions of the manufacturer of the washer/disinfectant.

4. Start the Vario TD programme including thermal disinfection. Thermal disinfection takes place allowing for the A<sub>0</sub> value and observing national provisions (prEN/ISO 15883).
5. On completion of the cycle remove sonic tip from the washer/disinfectant and dry (preferably with compressed air as recommended by the Commission for Hospital Hygiene and Infectious Disease Prevention of the Robert Koch Institute). When drying the inside of the instruments, please make sure to hold the compressed air gun close to the cooling passage, in order to guarantee that sufficient air flows through the sonic tip.
6. Visual examination to ensure that the instrument is clean and undamaged. If after mechanical reprocessing there are still visible residues of contamination, repeat the cleaning and disinfecting process until no visible contamination is left.  
Manual reprocessing (alternative):  
1. Place sonic tip onto a 10 ml syringe and rinse with detergent/disinfectant. In case of visible contamination emerging from the perforation, rinse again with detergent/disinfectant.  
2. Remove stubborn contamination with a nylon brush under running water, turning the sonic tip constantly.  
3. Rinse sonic tip thoroughly with running water.  
4. Visual examination to ensure that the instrument is clean. If there are still visible residues of contamination, repeat the cleaning process until no visible contamination is left.
5. Place sonic tip in a suitable container (e.g. Komet 9952) into the ultrasonic device or instrument bath filled with detergent disinfectant.
6. During chemical disinfection in the ultrasonic device or instrument bath, observe the instructions of the manufacturer regarding concentration and immersion time. Be sure to observe the full correct immersion time which does not start until the last tip has been immersed. Attention: do not exceed 45°C (risk of protein coagulation)!
7. On completion of the immersion time, rinse sonic tip thoroughly with suitable water (preferably with demineralised water to avoid residues of lime). Rinse the perforation of the sonic tip with at least 10 ml of demineralised water by means of a cannula, to ensure that no residues of the detergent are left in the cooling passage.
8. Dry sonic tip (preferably with compressed air as recommended by the Commission for Hospital Hygiene and Infectious Disease Prevention of the Robert Koch Institute). When drying the inside of the instruments, please make sure to hold the compressed air gun close to the cooling passage, in order to guarantee that sufficient air flows through the sonic tip.
9. Visual examination to ensure that the instrument is clean and undamaged. If there are still visible residues of contamination, repeat the cleaning and chemical disinfecting process until no visible contamination is left. According to the recommendations of the Commission for Hospital Hygiene and Infectious Disease Prevention of the Robert Koch Institute, the further reprocessing should be carried out mechanically. Work instructions with regard to proper reprocessing of instruments according to DIN EN ISO 17664, the present instructions for use and more detailed information, as amended, can be downloaded from our web site www.kometdental.de or requested from the manufacturer Gebr. Brasseler.  
Thermal disinfection in the autoclave:  
Classed as risk group semi-critical B, these sonic tips are suitable for thermal disinfection in a steam sterilizer, in an unwrapped condition inside a suitable container (e.g. 9952).  
For this procedure, the following instructions have to be observed:  
Steam sterilisation using a vacuum process at 134°C in a device that complies with the provisions of DIN EN 13060; with validated processes.  
- fractionated pre-vacuum (type B)  
- sterilisation temperature: 134°C  
- hold time: at least 5 minutes (full cycle)  
- drying time: at least 10 minutes  
In order to prevent staining and corrosion, the steam must be free of particles. The recommended limits for particle contents in feed water and condensed steam are defined by standard DIN EN 13060. Make sure not to exceed the maximum capacity of the sterilizer when sterilizing several instruments. Follow the instructions of the device manufacturer.  
Universally valid notes:  
Observe the legal provisions regarding the reprocessing of medical products valid in your country (e.g. www.rki.de). The manufacturer confirms that the above detailed reprocessing methods are suitable for preparing the above named instrument group to enable their reuse. The person carrying out the reprocessing is responsible for ensuring that the applied method is carried out with appropriate equipment, materials and trained personnel at the reprocessing site and that it actually achieves the desired result. To guarantee this, routine controls of the validated mechanical and/or manual preparation methods are necessary. Any deviation from the above detailed process must be carefully checked by the operator to ensure effectiveness and to avoid possible adverse consequences.

- Accessories:  
Komet tip changer SF4887 (supplied with the hand piece)  
Komet tip changer SF1975  
For reprocessing in the Miele washer/disinfectant: Komet rinse adapter SF1978/L for Komet sonic tips  
Komet sterilization container 9952
- 8. Storage**  
Store packed instruments away from UV radiation and heat in a clean and dry environment. Do not store in the same room as solvents or chemicals.
- 9. Disposal**  
The instruments have to be disposed of in tight, non-breakable and non-pierceable containers (protection from contamination).
- 10. Safety and liability**  
Only use the instruments according to their intended purpose. The non-observance of these safety instructions can lead to injuries or damage to the power system. It is the responsibility of the user to check the products prior to use to ensure that they are suited for the intended purpose. In case of contributory negligence by the user, Gebr. Brasseler partially or totally declines liability for all resulting damages, particularly due to non-observance of our recommendations for use or warnings as well as inadvertent misuse by the user.
- 11. Procedure in case of serious adverse events**  
In case of serious adverse events, please proceed in line with the legal regulations in your country.

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