


Komet Dental  
 Gebr. Brasseler GmbH & Co. KG  
 Trophagener Weg 25  
 32657 Lemgo · Germany  
 Telefon+49 (0)5261 701- 0  
 Telefax+49 (0)5261 701- 289  
 info@kometdental.de  
 www.kometdental.de

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**Ultrasonic tips for prophylaxis and periodontal treatments**

**1. Indication**

- Ultrasonic scaler tips are used in prophylactic treatments. They are intended for the supra and subgingival scaling of calculus, removal of soft and hard concretions.
- Ultrasonic periodontal tips are used for the removal of soft and hard dental plaque from deeper periodontal pockets.

**2. Contraindication**

- Cardiac pacemakers and defibrillators: Patients wearing pacemakers or defibrillators are excluded from treatments with ultrasonic tips because the ultrasonic vibrations of the instruments and/or hand pieces may impair the proper function of these devices.

**3. Warnings**

- The useful life of ultrasonic tips is limited. They therefore have to be used with the utmost care.
- **Attention:**
- **Risk of fracture!** Spontaneous and premature fracture must be expected. The ultrasonic tips may fracture unexpectedly and without warning. In order to reduce the risk of fracture, make absolutely sure to stick to the correct power setting (see proper use) and observe all recommendations for use. Never allow the instruments to jam. Do not use instruments for leverage. Avoid excessive contact pressure. Use ultrasonic tips with the utmost care as even minor mechanical damage can provoke the fracture of the instruments.
- **Risk of ingestion!** Fracture of the ultrasonic tip leads to the risk of the tip being accidentally swallowed. Take all appropriate steps to prevent this.
- **Premature wear** Due to the nature of their range of application, the functional parts of these ultrasonic tips will always be subject to moderate wear, even when the tips are applied correctly and according to their intended use. A test card (412830) can be used to determine the extent of wear of the tips (3L, 4L, PE1, 4PS, PE9,SI1).
- **Excessive heat generation!** Excessive contact pressure and/or worn or blunt instruments can cause hazardous heat generation which increases the risk of thermal necrosis (damage to tissue). Excessive contact pressure therefore has to be avoided. Damaged ultrasonic tips (changed shape, corroded surfaces or defect thread etc.) have to be replaced.
- **Excessive heat generation!** The use of a high volume evacuator might lead to an excessive evacuation of the cooling liquid, in which case the dental surface or the instrument will not be cooled sufficiently. The use of a high volume evacuator should therefore be avoided.
- **Excessive heat generation!** The tip, the treated tooth or the drive can heat up very strongly and quickly in case of an insufficient flow rate of the cooling liquid through the ultrasonic drive. To prevent this, always ensure a flow rate of at least 30 ml/min. of water through the drive. Never use the tips without irrigation.
- **Non-sterile instruments:** Instruments that have not been sterilised before the first use or not correctly reprocessed after each subsequent use can cause infection. The instruments therefore have to be cleaned and sterilised prior to each use in compliance with the reprocessing directions indicated in the present instructions for use.
- **Damage to hard tissues:** Excessive contact pressure might damage the hard dental substance and/or root surfaces. Excessive contact pressure therefore has to be avoided. This applies in particular to the periodontal tips which may only be used with low contact pressure.
- **Damage to hard tissues:** Ultrasonic tips perform longitudinal oscillations relative to the longitudinal axis of the instrument. In order to prevent damage to the hard dental substance, avoid frontal positioning of the instrument on the dental surface.
- **Ultrasonic tips may not be used on metal or ceramic restorations or dentures as these might get damaged by the use these instruments.**
- **Do not use the hand piece in areas that are exposed to explosion hazards, especially in the vicinity of flammable anaesthetics or gases.**
- **The ultrasonic tips may only be used by properly trained and qualified personnel.**
- **Always wear protective gloves when checking, inserting or removing the ultrasonic tips to avoid risks of injury and infection.** Exchange tips regularly.
- **Risk of injury:** When not in use, place tip changer onto the screwed-in ultrasonic tip, this will serve as a tip protector.
- **Risk of injury:** Protect soft tissue from contact with the active tip to prevent burns.
- **Risk of fracture:** Subject the tip to axial loads prior to each use.
- **Incorrect use** can lead to the above detailed undesirable effects as well as to injury to the user, the patient or third persons.
- **To ensure risk-free use, always follow the instructions and safety and maintenance recommendations provided by the manufacturer of the device.**

**4. Start-up**

- The ultrasonic tips delivered in a non-sterile condition have to be cleaned, disinfected and sterilised prior to first use.
- The ultrasonic tips can be used in the following hand-pieces:
  - PerioSonic for PerioScan.
  - SIROSONIC/L and SIROSONIC TL, used in the treatment units TENEO and SINIUS.
  - SIROSON S and SIROSON C8, used in the treatment unit M1.
- To insert or remove the ultrasonic tips, use the enclosed torque wrench only. Screw the chosen tip into the hand piece, turning the torque wrench in a clockwise direction. To fasten the tip correctly, turn the torque wrench another quarter-turn after the tip has been completely screwed in. To remove the tip, place the torque wrench onto the tip and unscrew in an anticlockwise direction.

**5. Proper use**

- **Start-up ultrasonic tip** prior to application to the tooth, supplying an adequate amount of spray cooling. All ultrasonic tips described in the present instructions require plenty of cooling liquid, at least 30 ml/min. Recommended use of each individual ultrasonic tip:
  - Group of tips used for scaling:
    - 1L.SI.1.: Coarse supragingival scaling in all quadrants. The tip can be used from lingual or buccal direction. Not suitable for use in interproximal areas. Can be used at a power setting of up to 100%.
    - 3L.SI.1.: Removal of supragingival and subgingival concretions. Versatile bilateral use thanks to fine instrument tip. Suitable for interdental use. Can be used at a power setting of up to 80%.
    - 4L.SI.1.: Removal of supragingival or subgingival concretions in all quadrants. Particularly suitable for tooth necks and interproximal regions. Easy access to rear molars thanks to its curved shape of the tip. Can be used at a power setting of up to 70%.
  - Group of tips used for periodontal treatments:
    - PE1.SI.1.: Removal of subgingival deposits mainly in the anterior region, suitable for cleaning deep periodontal pockets and dental roots. Can be used at a power setting of up to 70%.
    - PE2.SI.1. + PE3.SI.1.: Curved versions of the PE1, bent to the left or right. Suitable for cleaning roots and rinsing deep periodontal pockets. Mainly used in lateral teeth. Can be used at a power setting of up to 100%.
    - 4PS.SI.1.: For cleaning roots and rinsing deep periodontal pockets. Can be used at a power setting of up to 80%.
    - PE5.SI.1. and PE6.SI.1.: Ultrasonic tips with curve to the left or right for cleaning roots and rinsing deep periodontal pockets in the lateral region. Can be used at a power setting of up to 100%.
    - PE9.SI.1.: The scope of application of this tip is situated between scaling and periodontology, which means that it is suitable for removing concretions in the supragingival and subgingival region of all quadrants. Can be used at a power setting of up to 80%.
- **Maintenance**

Place ultrasonic tips in a cleaning/disinfection tank (Fräsator) filled with a suitable detergent/disinfectant (e.g. Komet DC1) immediately after use in the mouth to prevent drying of residues on the instruments (protein fixation). Hold ultrasonic tips at an angle when placing them into the liquid to ensure that all hollow spaces are reached. It is recommended to reprocess the instruments within one hour of use at the very latest. The instruments should be in the cleaning/disinfection tank (Fräsator) when transported to the site where the reprocessing is to take place.

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**Instructions for use**

Please read carefully and retain for future reference.

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**7. Reprocessing**

- Manual pre-cleaning (prior to mechanical and manual reprocessing):
  - Remove ultrasonic tip from the cleaning/disinfection tank (Fräsator) immediately before reprocessing.
  - Check patency of the cooling passage by guiding a fine needle or a cleaning wire (e.g. Komet 97509) through the passage. Replace instrument if the passage is obstructed.
  - Thoroughly rinse under running water to remove surface contamination. Remove all stubborn contamination under water level with a nylon brush, constantly turning the ultrasonic tip.
- Mechanical reprocessing (validated):
  - Rinse ultrasonic tip thoroughly under running water to prevent any residues of the detergent/disinfectant from getting into the washer/disinfector.
  - Unscrew the nozzle on the injector bar of the washer/disinfector. Then screw the external thread of the rinse adapter 1978.SI1 or 1978L.SI1 in the silicon holders on the injector bar. Now screw the ultrasonic tip in the internal thread at the top of the rinse adapter. Observe the directions in the instructions for use provided with the 1978.SI1 or 1978L.SI1.
  - Put chemical detergent into the washer/disinfector, following the indications on the label and the instructions of the manufacturer of the washer/disinfector.
  - Start the Vario TD programme including thermal disinfection. Thermal disinfection takes place allowing for the As value and observing national provisions (pREN/ISO 15883).
  - On completion of the cycle remove bar block and inserts from the washer/disinfector and dry (preferably with compressed air www.rki.de). Afterwards the Commission for Hospital Hygiene and Infectious Disease Prevention of the Robert Koch Institute). When drying the inside of the instruments, please make sure to hold the compressed air gun close to the cooling passage, in order to guarantee that sufficient air flows through the ultrasonic tip.
  - Visual examination to ensure that the instrument is clean and undamaged. If after mechanical reprocessing there are still visible residues of contamination, repeat the cleaning and disinfecting process until no visible contamination is left.
- Manual reprocessing (alternative):
  - Place ultrasonic tip onto a 10 ml syringe and rinse with detergent/disinfectant. In case of visible contamination emerging from the perforations, rinse again with detergent/disinfectant.
  - Remove stubborn contamination with a nylon brush under running water, turning the ultrasonic tip constantly.
  - Rinse ultrasonic tip thoroughly with running water.
  - Visual examination to ensure that the instrument is clean. If there are still visible residues of contamination, repeat the cleaning process until no visible contamination is left.
  - Place ultrasonic tip in a suitable container into the ultrasonic device or instrument bath filled with detergent/disinfectant.
  - During chemical disinfection in the ultrasonic device or instrument bath, observe the instructions of the manufacturer regarding concentration and immersion time. Be sure to observe the full correct immersion time which does not start until the last tip has been immersed. Attention: do not exceed 45°C (risk of protein coagulation)!
  - On completion of the immersion time, rinse ultrasonic tip thoroughly with suitable water (preferably with demineralised water to avoid residues of lime). Rinse the perforation of the ultrasonic tip with at least 10 ml of demineralised water by means of a cannula, to ensure that no residues of the detergent are left in the cooling passage.
  - Dry ultrasonic tip (preferably with compressed air as recommended by the Commission for Hospital Hygiene and Infectious Disease Prevention of the Robert Koch Institute). When drying the inside of the instruments, please make sure to hold the compressed air gun close to the cooling passage, in order to guarantee that sufficient air flows through the ultrasonic tip.
  - Visual examination to ensure that the instrument is clean and undamaged. If there are still visible residues of contamination, repeat the cleaning and chemical disinfecting process until no visible contamination is left.
  - According to the recommendations of the Commission for Hospital Hygiene and Infectious Disease Prevention of the Robert Koch Institute, reprocessing should preferably be carried out mechanically. Information on proper reprocessing of instruments according to DIN EN ISO 17664, the present instructions for use, as amended, and more detailed information can be downloaded from our web site www.kometdental.de or requested from the manufacturer Gebr. Brasseler.
- Sterilisation in the autoclave:
  - Make sure that the packaging is suitable for ultrasonic tips classed as risk group critical B and the chosen method of sterilisation.
  - Single pack: The packaging must be large enough to ensure that there is no pressure on the seal. Multipack: Place ultrasonic tips inside the torque wrench into a suitable sterilisation container (e.g. 97607). For sterilisation, the container has to be sealed in suitable packaging.
  - For this procedure, the following instructions have to be observed:
    - Steam sterilisation using a vacuum process at 134°C in a device that complies with the provisions of DIN EN 13060; with validated processes.
    - fractionated pre-vacuum (type B)
    - sterilisation temperature: 134°C
    - hold time: at least 5 minutes (full cycle)
    - drying time: at least 10 minutes
  - In order to prevent staining and corrosion, the steam must be free of particles. The recommended limits for particle contents in filter water and condensed steam are defined by standard DIN EN 13060. Make sure not to exceed the maximum capacity of the sterilizer when sterilizing several instruments. Follow the instructions of the device manufacturer.
  - Universally valid notes:
    - Observe the legal provisions regarding the reprocessing of medical products valid in your country (e.g. www.rki.de). The manufacturer confirms that the above detailed reprocessing methods are suitable for preparing the above named instrument group to enable their reuse. The person carrying out the reprocessing is responsible for ensuring that the applied method is carried out with appropriate equipment, materials and trained personnel at the reprocessing site and that it actually achieves the desired result. To guarantee this, routine controls of the validated mechanical preparation methods are necessary. Any deviation from the above detailed process must be carefully checked by the operator to ensure effectiveness and to avoid possible adverse consequences.

**8. Storage**

Store packed instruments away from UV radiation and heat in a clean and dry environment. Do not store in the same room as solvents or chemicals.

**9. Disposal**

The instruments have to be disposed of in tight, non-breakable and non-pierceable containers (protection from contamination).

**10. Safety and liability**

Only use the instruments according to their intended purpose. The non-observance of these safety instructions can lead to injuries or damage to the power system. It is the responsibility of the user to check the products prior to use to ensure that they are suited for the intended purpose. In case of contributory negligence by the user, Gebr. Brasseler partially or totally declines liability for all resulting damages, particularly due to non-observance of our recommendations for use or warnings as well as inadvertent misuse by the user.

**11. Procedure in case of serious adverse events**

In case of serious adverse events, please proceed in line with the legal regulations in your country.

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© Gebr. Brasseler GmbH & Co. KG