Komet Dental
 Gebr. Brasseler GmbH & Co. KG
 Trophagener Weg 25
 32657 Lemgo - Germany
 Telefon+49 (0)5261 701- 0
 Telefax+49 (0)5261 701- 289
 info@brasseler.de
 www.brasseler.de

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F6 SkyTaper System of NITI files with a taper of .06 for the mechanical preparation of root canals

- Indication**
 The instruments are intended for the mechanical removal of diseased tissue and root filling materials within the course of a root canal treatment.
- Contraindication**
 The contra-indications are (e. g.):
 - pronounced abrupt curves
 - severely calcified root canals
 - preparation of milk teeth
- Proper use**

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- Probe/glide path ISO 025 : Use F6 SkyTaper in size 035
- Probe/glide path ISO 030 : Use F6 SkyTaper in size 040
- Insert the files into the canal in rotary mode. The instruments are used in "pecking motion", i.e. the files are moved up and down the canal wall during a period of max. 5 seconds, according to the single-length principle.
- Clean the instruments periodically during the treatment, for example in the interim storage block 595. Rinse the root canal thoroughly in order to flush out all debris.

- Further treatment**
 After the administration of a medicament, if necessary, the prepared root canal can be filled according to the usual methods, for example with condensed gutta-percha such as Komet F6 SkyTaper gutta-percha points GPF06 a carrier based filling system, e.g. Komet F360 Fill or a thermo plastic procedure.

- Reprocessing**
 The present instruments are supplied in sterile packaging. They can therefore be used without previous processing, provided that the packaging is intact. Do not use instruments inside damaged packaging because they are no longer sterile. The sterile packed F6 SkyTaper files are not designed for reprocessing, i.e. cleaning/disinfection and sterilisation to enable their reuse. The instruments can only be used once (one patient only). Used instruments have to be discarded because a risk-free reuse of the products cannot be guaranteed. Instruments in damaged packaging are no longer sterile and should not be used. Information on proper reprocessing of instruments according to DIN EN ISO 17664, the present instructions for use, as amended, and more detailed information can be downloaded from our web site www.kometdental.de or requested from the manufacturer Gebr. Brasseler.

- Storage**
 Store packed instruments away from UV radiation and heat in a clean and dry environment. Do not store in the same room as solvents or chemicals.

- Disposal**
 The instruments have to be disposed of in tight, non-breakable and non-pierceable containers (protection from contamination).

- Safety and liability**
 See also "instructions for use and safety recommendations" in the current Komet Dental catalogue. The user is responsible for checking the products prior to use as to whether they are suited for the intended purpose. The user is responsible for the appropriate use of the instruments. In case of contributory negligence by the user, Gebr. Brasseler partially or totally declines liability for all resulting damages, particularly if these are due to nonobservance of our recommendations for use or warnings as well as inadvertent misuse by the user.

- Procedure in case of serious adverse events**
 In case of serious adverse events, please proceed in line with the legal regulations in your country.

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© Gebr. Brasseler GmbH & Co. KG