 Komet Dental
 Gebr. Brasseler GmbH & Co. KG
 Trophagener Weg 25
 32657 Lemgo - Germany
 Telefon+49 (0)5261 701-0
 Telefax+49 (0)5261 701-289
 info@kometdental.de
 www.kometdental.de

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Instructions for use
 Please read carefully and retain for future reference.

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Ultrasonic tips for prophylaxis and periodontal treatments

1. Indication
 - Ultrasonic scaler tips are used in prophylactic treatments. They are intended for the sub and supragingival scaling of calculus, removal of soft and hard concretions.
 - Ultrasonic periodontal tips are used for the removal of soft and hard dental plaque from deeper periodontal pockets. They are also suited for the finishing of root surfaces after cleaning as part of periodontal flap surgery.

2. Contraindication
 - Cardiac pacemakers and defibrillators: Patients wearing pacemakers or defibrillators are excluded from treatments with ultrasonic tips because the ultrasonic vibrations of the instruments and/or hand pieces may impair the proper function of these devices.

3. Warnings
 The useful life of ultrasonic tips is limited. They therefore have to be used with the utmost care. Attention:
 - Fracture hazard! Spontaneous and premature fracture must be expected. The ultrasonic tips may fracture unexpectedly and without warning. In order to reduce the risk of fracture, make sure to stick to the correct power setting (see instructions for use) and observe all recommendations for use. Never allow the instruments to jam. Do not use instruments for leverage. Avoid excessive contact pressure. Use ultrasonic tips with the utmost care as even minor mechanical damage can provoke the fracture of the instruments.
 - Risk of ingestion! Fracture of the ultrasonic tip leads to the risk of the tip being accidentally swallowed. Take all appropriate steps to prevent this.
 - Premature wear! Due to the nature of their range of application, the functional parts of these sonic tips will always be subject to moderate wear, even when the tips are applied correctly and according to their intended use. A test card (413586) can be used to determine the extent of wear of tips (LU8 - LU11, LU63). Excessive contact pressure might cause damage to the working part. Therefore, the instruments have to be checked after each use. Damaged or worn instruments have to be replaced.
 - Excessive heat generation! Excessive contact pressure and/or worn or blunt instruments can cause excessive heat generation which increases the risk of thermal necrosis (damage to tissue). Excessive contact pressure therefore has to be avoided. Damaged ultrasonic tips (changed shape, corroded surfaces or defect thread etc.) have to be replaced.
 - Excessive heat generation! The use of a high volume evacuator might lead to an excessive evacuation of the cooling liquid, in which case the dental surface or the instrument will not be cooled sufficiently. The use of a high volume evacuator should therefore be avoided.
 - Excessive heat generation! The tip, the treated area on the tooth or the hand piece can heat up very strongly and quickly in case of an insufficient flow rate of the cooling liquid through the ultrasonic drive. The quantity of water supplied can be adjusted by means of a regulating ring on the ultrasonic hand piece. The quantity of water is correctly set when the oscillating ultrasonic tip provides an even jet of water spray. All instruments have to be used with a high flow rate. Do not use any circumstances work without cooling.
 - Non-sterile instruments: Instruments that have not been sterilised before the first use or not correctly reprocessed after each subsequent use can cause infection. The instruments therefore have to be cleaned and sterilised prior to each use in compliance with the reprocessing directions indicated in the present instructions for use.
 - Damage to hard tissues: Excessive contact pressure might damage the hard dental substance and/or root surfaces. Excessive contact pressure therefore has to be avoided. This applies in particular to the periodontal tips which may only be used with low contact pressure.
 - Damage to hard tissues: Ultrasonic tips perform longitudinal oscillations relative to the longitudinal axis of the instrument. In order to prevent damage to the hard dental substance, avoid frontal positioning of the instrument on the dental surface.
 - Ultrasonic tips may not be used on metal or ceramic restorations or dentures as these might get damaged by the use these instruments.
 - Do not use the hand piece in the vicinity of flammable anaesthetics or gases.
 - The ultrasonic tips may only be used by properly trained and qualified personnel.
 - Always wear protective gloves when checking, inserting or removing the ultrasonic tips to avoid risks of injury and infection. Exchange tips regularly.
 - Risk of injury: When not in use, place tip changer onto the screwed-in ultrasonic tip, this will serve as a tip protector.
 - Risk of injury: Protect soft tissue from contact with the active tip to prevent burns.
 - Risk of fracture: Subject the tip to axial loads prior to each use.
 - Incorrect use can lead to the above detailed undesirable effects as well as to injury to the user, the patient or third persons.
 - To ensure risk-free use, always follow the instructions and safety and maintenance recommendations provided by the manufacturer of the device.

4. Start-up
 Ultrasonic tips delivered in a non-sterile condition have to be cleaned, disinfected and sterilised prior to first use. The ultrasonic tips can be used in the following hand-pieces:
 - In the hand-piece PIEZOlux of the company KaVo.
 - In the hand-piece SONOsoft and SONOsoft LUX of the company KaVo.
 In treatment units set to level "E", the oscillation intensity of the hand piece can only be regulated by actuating the foot pedal. Having selected memory locations 1 - 3, the oscillating intensity can be set in steps of 0.25, from minimum 1 to maximum 10. Scalers and periodontal tips can be used up to power level 10, with the exception of the LU8 which can only be used up to power level 6.
 To insert or remove the ultrasonic tips, use the torque wrench SF4887 or the corresponding torque wrench of the system used. Screw the chosen tip into the hand piece, turning the torque wrench in a clockwise direction. Correct fastening of the ultrasonic tip can be recognised by the typical audible and perceivable locking of the ratchet. Further turning of the torque ratchet does not increase the retention force of the ultrasonic tip in the handpiece. To remove the tip, place the torque wrench onto the tip and apply pressure in an anticlockwise direction.

5. Proper use
 Start up ultrasonic tip prior to application to the tooth, supplying an adequate amount of spray cooling. All ultrasonic tips described in the present instructions have to be used with a high flow rate (see item 3). Recommended use of each individual ultrasonic tip:
 Group of tips used for scaling:
 LU9.KA1: Removal of large amounts of supragingival calculus in all quadrants.
 LU10.KA1: Removal of supragingival calculus in all quadrants.
 LU11.KA1: Removal of large amounts of calculus in all quadrants.
 LU8.KA1: This ultrasonic tip can be used both for scaling and periodontal treatment, in particular removal of subgingival and supragingival concretions in all quadrants.
 Group of tips used for periodontal treatment:
 LU63.KA1: The following tips are suitable both for initial and follow-up treatments.
 LU63.KA1: Subgingival removal of plaque and deposits on root surfaces in all quadrants, in particular for the anterior region and smooth surfaces. This tip is also suitable for polishing root surfaces after cleaning.
 LU64.KA1: Subgingival removal of plaque and deposits on root surfaces, in particular in the left mandibular vestibular, right mandibular lingual, right maxillary vestibular and left maxillary palatal regions.
 LU65.KA1: Subgingival removal of plaque and deposits on root surfaces, in particular in the left mandibular lingual, right mandibular vestibular, right maxillary palatal and left maxillary vestibular regions.

6. Maintenance
 Place ultrasonic tips in a cleaning/disinfection tank (Fräsator) filled with a suitable detergent/disinfectant (e.g. Komet DC1) immediately after use in the mouth to prevent drying of residues on the instruments (protein

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© Gebr. Brasseler GmbH & Co. KG