 Komet Dental
Gebr. Brasseler GmbH & Co. KG
Trophagener Weg 25
32657 Lemgo · Germany
Telefon +49 (0) 5261 701-0
Telefax +49 (0) 5261 701-289
info@kometdental.de
www.kometdental.de

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Sonic tips for gentle surgery, developed by Dr. Agabiti

1. Indications

The Komet sonic tips with Quick connection SFSQ 100, 101 and 102 are designed for oscillating work on bones in addition to rotary work with hand pieces and contra angles. These sonic tips are particularly suited for the fragmentation of bones, alveolotomy, osteotomy – in particular for exposure of impacted teeth -, orthodontic corticotomy and for bone splitting. The saw blade of an SFSQ100/101/102 tip is provided with sharp blades. These sonic tips may only be used by dentists or medical personnel who are sufficiently qualified in dental implantology and the relevant fields of surgery.

2. Contraindications

All general contraindications for dental surgery must be taken into account. All commonly known contraindications in oral maxillo facial surgery (e.g. implantation, jaw surgery) also have to be considered.

3. Warnings

The useful life of the sonic tips is limited. They therefore have to be used with the utmost care. Attention:
- Risk of fracture! The sonic tips may fracture unexpectedly and without warning. In order to reduce the risk of fracture, make sure to stick to the correct power setting (see instructions for use) and observe all recommendations for use. Never allow the instruments to jam. Do not use instruments for leverage. Avoid excessive contact pressure. Use sonic tips with the utmost care to even minor mechanical damage can provoke the fracture of the instruments.
- Risk of ingestion! Fracture of the sonic tip leads to the risk of the tip being accidentally swallowed. Take all appropriate steps to prevent this.

- Premature wear! Inspect instruments after every use for signs of wear or damage. Exchange if necessary. Work on hard dental substance can cause moderate wear. Likewise, excessive contact pressure can lead to damage to the working part.

- Excessive heat generation! Excessive contact pressure and/or worn or blunt instruments can cause hazardous heat generation which increases the risk of thermal necrosis (damage to tissue). Excessive contact pressure therefore has to be avoided. Damaged sonic tips (changed shape, corroded surfaces, void spots or defect thread etc.) have to be replaced. The sonic tips are medical devices and as such subject to the applicable legal provisions of the relevant country. In compliance with such regulations, the sonic tips may only be used by properly qualified operators within the scope of the indicated application and according to the valid industrial Health& Safety regulations and observing the valid rules for accident prevention. The provisions further stipulate that all equipment used must be free of defects and suitable for the intended application. Furthermore, the patient and third persons must be protected from hazards. Possible contamination caused by the product must be ruled out. Safety control Fracture may occur due to the tips being subjected to continuous or sudden stress (for example if the tips are dropped on the floor or through a mechanical change of the original shape). Check the tips prior to each use to make sure they are safe to use. To do this, apply slight pressure with your thumb and index finger (fig. 1). For additional safety, load the tips mechanically with approx. 10N (1kg), without function (fig. 2). Always wear protective gloves when checking, inserting or removing the sonic tips to avoid risks of injury and infection. We recommend that the tips be exchanged regularly.

4. Appropriate use

Sonic tips supplied in non-sterile condition have to be prepared prior to first use. The sonic tips for oral surgery made by Komet and developed by Dr. Ivo Agabiti (www.sonosurgery.it) including the letter "Q" in the first part of the figure number, can be used:

- In the Komet sonic handpiece Quick SFQ2008L and

- SONICflex™ Quick Series 2008/S/L/LS made by the co. KaVo

The above listed sonic hand pieces are provided with a MULTIflex® quick coupling or with a Sirona® coupling (2008S/LS). To insert the sonic tip, screw the cooling adapter SFQ1979 manually on the chisel tip. Together, the products are screwed in the handpiece using the tip changer SF4887 (enclosed with the handpieces), or alternatively with the SF1975, turning the tool in clockwise direction. Correct fastening of the tip can be recognised by the typical audible locking of the tip changer. When using the tip changer, make sure that the tip faces towards the semi-radius during insertion (fig. 3). To remove the tip, place tip changer on the tip and turn anticlockwise. The cooling adapter is separated from the tip with the mounting wrench 566 (fig. 4).

Note: Please observe the directions in the instructions supplied with your sonic hand piece.

5. Recommendations for use

The sonic tips may only be used with the authorised power systems and permitted power settings. The power settings 2 (gentle) and 3 (standard) may be used. This applies to both the Komet sonic hand piece SFQ2008L and the KaVo SONICflex™ 2008. The power setting is regulated by means of a regulating ring. Please observe the recommendations in the instructions for use provided with the sonic hand pieces. The sonic tips are at their most effective when used with low contact pressure. Excessive contact pressure may lead to undesirable effects (see section 3 "warnings"). Given that these sonic tips may fracture without warning even when used according to instruction and their designated use, always plan two sonic tips for the operation, i.e. always have a reserve at hand so that the planned operation can still be finished even if one tip has broken.

Provision of cooling

The surgical sonic tips are cooled with at least 50 ml/min of an externally fed, sterile cooling liquid. Make sure to turn off the water supply of the dental treatment unit. It is absolutely essential to thoroughly and continuously irrigate the site at all times during the operation. The sterile cooling agent can be supplied manually, for example externally by means of a large syringe with a blunt needle held by the dental chair assistant. The sterile cooling can also be provided by the cooling adapter SFQ1979. The flexible tube supplying the cooling agent is placed onto the small tube of the cooling adapter. Controlled supply of the cooling agent by means of a suitable surgical power system is preferable, for example with the KaVo INTRASurg 1000 for oscillating work. There are various compression systems available on the market that ensure pressure supported transport of the cooling agent, for example:

- OMNIA system: manual supply of a sterile cooling agent by means of a compression pouch, available in Germany from the company Implantis Disposable Medical Devices (www.implantis.eu). They can be contacted as follows: From Germany: (0180) 5 90 75 90, free fax (00800) 05 05 51 23

From Austria: (0049 180) 5 90 75 90, free fax (00800) 05 05 51 23

Article no. H2204 Manual cooling system with compression pouch (Omnia no. 30E0008.00), article no. H2204 Manual cooling system with flux control, to be attached to the hand piece (Omnia no. 32F0157). For further information, visit www.omniasrl.com.

6. Maintenance/care

Place sonic tips and cooling adapter in a cleaning/disinfection tank (Fräsator) filled with a suitable detergent/disinfectant (e.g. Komet® DC17) immediately after use in the mouth to prevent drying of residues on the instruments (protein fixation). Hold sonic tips at an angle when placing them into the liquid to ensure that all hollow spaces are reached. It is recommended to reprocess the instruments within one hour of use at the very latest. The instruments should be in the cleaning/disinfection tank (Fräsator) when transported to the site where the reprocessing is to take place.

7. Reprocessing

Sonic tips delivered in non-sterile condition have to be prepared prior to first use.

Manual pre-cleaning (prior to mechanical and manual reprocessing):

1. Remove sonic tip and cooling adapter from cleaning/disinfection tank (Fräsator) immediately before reprocessing.

2. Separate the sonic tip from the cooling adapter before the start of the cleaning process.

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© 12/2015, Gebr. Brasseler GmbH & Co. KG