Komet Dental  
Gebr. Brasseler GmbH & Co. KG  
Trophagener Weg 25  
32657 Lemgo · Germany  
Telefon +49 (0) 5261 701-0  
Telefax +49 (0) 5261 701-289  
info@kometdental.de  
www.kometdental.de




ER TitanPost X Coated 6 mm  
© Gebr. Brasseler GmbH & Co. KG  
311797 - Rev 2016-07




 Stiftsysteme, post systems, systèmes de tenons, sistemas de pernos, sistemi di perni

Größe Size Taille Tamaño Misura	Farbe: orange Colour: orange Couleur: orange Color: naranja Colore: arancione	Größe Size Taille Tamaño Misura	Farbe: rot Colour: red Couleur: rouge Color: rojo Colore: rosso	Werkstoff Material Matériau Material Materiale
	<b>070</b>		<b>090</b>	

 ② TPXCL6.000	 ② TPXCL6.000	Reintitan Pure titanium Titane pur Titanio puro Titanio puro
--	---	--

 183LB.204	 183LB.204	Rostfreier Stahl Stainless Steel Acier inoxydable Acero inoxidable Acciaio inossidabile
 1.000 - 5.000 min <sup>-1</sup>  20.000 min <sup>-1</sup>	 1.000 - 5.000 min <sup>-1</sup>  20.000 min <sup>-1</sup>	

 196S.204	 196S.204	Rostfreier Stahl Stainless Steel Acier inoxydable Acero inoxidable Acciaio inossidabile
 1.000 - 2.000 min <sup>-1</sup>  20.000 min <sup>-1</sup>	 1.000 - 2.000 min <sup>-1</sup>  20.000 min <sup>-1</sup>	

 196DS.644	 196DS.644	Rostfreier Stahl · diamantiert Stainless Steel · diamond coated Acier inoxydable · diamanté Acero inoxidable · diamantado Acciaio inossidabile · diamantato
---	--	---

Tab. 1

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ER TitanPost X Coated 6 mm for build-ups with plastic materials, made of pure titanium, pre-fabricated with SiO<sub>2</sub> and SiH<sub>4</sub>, and provided with an adhesion enhancing layer of polymer for an adhesive bond.

1. System overview  
See table 1

2. Indication  
Prosthetic restoration of teeth with partially or deeply damaged crown using the adhesive technique.

Basic requirements:

- Remaining substance of 2,0 mm in height around the crown (ferrule design) or surgical crown extension.
- Successfully concluded endodontic treatment with remaining apical root filling of 4,0 mm.
- Wall adaptation of the post along its entire length.
- The lower retention element of the retention head (up to the first retention groove) can be countersunk in the root cross-section.

3. Contraindications

- Funnel-shaped cervical canal lumen, because the post is not adapted to the canal wall along its entire length.
- Subgingival caries.
- Circular isogingival destruction.
- Restorations where the adhesive technique cannot be used for fixing the post / building up the core.
- Traditional cementing of the coated posts might not result in the required adhesion which might lead to a loose fit of the post.

4. General preparatory procedure

- The instruments which are supplied non-sterile have to be prepared prior to first use.
- 4.1 The tooth has been prepared according to the indication and received endodontic treatment. In case of teeth with deeply damaged crown, reduce the stump to the sound tooth substance and prepare the surface to a concave shape by using a round diamond instrument, for example 801 (fig. 1).
- 4.2 Choose suitable post
- Select post on the basis of the radiograph according to the anatomic condition of the tooth to be restored. For stability reasons, the post has to fit closely to the canal walls and be anchored in the root canal as deeply as possible.
- 4.3 Removal of the root filling to the depth required with the pilot bur 183LB, according to the chosen post size.
- Mark the required preparation length on the pilot drill with a stopper.
- 4.4 Enlarge the root canal according to the chosen post.
- With the corresponding reamer 196S which is provided with a step to simultaneously prepare the central inlay cavity.
- Prepare inlay cavity to a slightly oval shape by means of a cylindrical diamond grinding instrument, e.g. 835KR, to protect the post/stump build-up against rotation.
- 4.5 Check proper fit of the prepared post size by inserting the selected post. Eliminate any imprecisions with the reamer 196S.
- 4.6 If necessary, take a radiograph for control.

5. Insertion of post and reconstruction

- 5.1 Mechanical conditioning of the canal wall with roughening instrument 196DS.
- Rotate the instrument manually 2 – 3 times in the root canal without applying pressure.
- 5.2 Rinse root canal and inlay cavity for example with 1% NaOCl, followed by distilled water. Dry with paper points.
- 5.3 Clean surface of the post: degrease with medical alcohol and dry with airstream.
- Note: Do not immerse pin in alcohol, just disinfect by wiping. The post must not be sterilized. Do not touch the coated part of the post (other than the handling element) with your fingers after disinfection and before inserting the post into the root canal.
- 5.4 The posts are designed for adhesive fixation.
- Chemical conditioning of enamel and dentin. Follow the manufacturer's instructions for the Composite System.
- Cover post with fixing composite.
- Follow the manufacturer's instructions for the Composite system.

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**Instructions for use**

Please read carefully and retain for future reference.

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- Introduce post slowly into the root canal, gently rotating the post during insertion. Fixate by applying pressure with your thumb.
- Allow the composite to polymerize and remove excess material, if necessary.
- Important: For adhesive fixation of the post do not use purely light-curing composite or a purely light-curing bonding agent.
- 5.5 After insertion of the root post, carefully snap off the handling element.
- 5.6 Conclusion of the definite preparation.
- Adaptation of the plastic build-up material. Follow the manufacturer's instructions for the Composite system. After this, the coronal section of the root canal should be completely covered with composite.

6. Further treatment

For the definite preparation, please observe the following rules:

- The crown stump has to be surrounded in a circumferential manner.
- For all-ceramic crowns we recommend shoulder preparations.

The crown has to surround the post reinforced root in an at least 2 mm wide parallel zone so that forces exerted are distributed over the entire tooth.

- Impression and temporary restoration in the usual manner.
- Manufacture and insert crown.
- Due to the biological, mechanical and static aspects of the definitive prosthesis, the long-term success of the restoration can be improved by replacing the temporary appliance with the definitive prosthesis as soon as possible.

7. Recommendations for maintenance

The instruments which are supplied non-sterile have to be prepared prior to first use. Root posts should be cleaned and disinfected with medical alcohol. Do not use sterile packed instruments inside damaged packaging. In case of sterile packed instruments, observe the date on the "use by" symbol. Root posts are single-use products. The reuse of these products poses a risk of infection and/or the safety of the products can no longer be guaranteed. The instruments are to be disinfected with anticorrosive disinfecting and cleaning agents for rotary instruments (e.g. Komet DC1). Regarding recommendations for use (immersion time, concentration, suitability) of the disinfecting or cleaning agents see manufacturers' instructions. The instruments can be reprocessed in the thermo disinfecter, provided that the agent used is suitable for rotary instruments (mildly alkaline detergent). Disinfecting and cleaning agents must be rinsed off thoroughly with water and the instruments must be dried carefully (preferably with compressed air according to the recommendations of the Commission for Hospital Hygiene and Infectious Disease Prevention of the Robert Koch Institute). Do not store instruments in wet or humid condition for a prolonged period of time. Make sure that instruments do not come in contact with each other during ultrasonic cleaning. Inspect cleaned instruments visually. Do not reuse damaged or blunt instruments. These have to be discarded.

Sterilisation:

Sterilisation of the instruments is carried out in the autoclave at 134°C. Observe instructions provided by the manufacturer. After sterilization, check instruments for defects on the surface (corrosion). Corroded instruments must not be used any longer. The operator of medical products is responsible for seeing that proper treatment is carried out by qualified personnel, using the appropriate materials and suited equipment, as recommended by the Commission for Hospital Hygiene and Infectious Disease Prevention of the Robert Koch Institute. Work instructions with regard to proper reprocessing of instruments according to DIN EN ISO 17664, the present instructions for use and more detailed information, as amended, can be downloaded from our web site [www.kometdental.de](http://www.kometdental.de) or requested from the manufacturer Gebr. Brasseler.

8. Storage

Store sterile packed instruments in their original packaging, away from light and heat in a clean and dry environment. Do not store in the same room as solvents or chemicals.

9. Safety and liability

See also "Instructions for use and safety recommendations" in the current Komet catalogue. The user is responsible for checking the products prior to use as to whether they are suited for the intended purpose. In case of contributory negligence by the user, Gebr. Brasseler partially or totally declines liability for all resulting damages, particularly if these are due to non-observance of our recommendations for use or warnings as well as inadvertent misuse by the user.