

 KFO, Orthodontics, Traitement orthodontique, Gnato - ortopedia, Ortodonzia

 Kavitätenpräparation, Cavity preparation, Préparation cavitaire, Preparación de cavidades, Preparazione di cavità

 Füllungsbearbeitung, Working on fillings, Préparation des obturations, Trabajo en obturaciones, Lavorazione di otturazioni

 Kronenpräparation, Crown preparation, Préparation coronaire, Preparación de coronas, Preparazione di corone

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EN
Sonic tips with Quick connection for stripping and shaping

1. Indication
Orthodontics:
- Stripping of front and lateral teeth during interproximal enamel reduction (IER). Use the sonic tips for stripping that are coated on the flat side.
- Final shaping of the front and lateral teeth with the sonic tips for shaping that are coated on the convex side.
Interproximal preparation of cavities:
- Smoothing of transitions at the interproximal preparation limit
Work on composite fillings:
- Anatomical shaping of the interproximal surfaces with the sonic tips for shaping that are coated on the convex side.

Crown core preparation:
- Elimination of the interproximal contact (separation) prior to the preparation with the tips for stripping that are coated on the flat side, in order to prevent damage to the adjacent teeth.

2. Contra-indications
There are no known contra-indications.

3. Warnings
The useful life of the sonic tips is limited. They therefore have to be used with the utmost care.
Attention:
- Risk of fracture! Spontaneous and premature fracture must be expected. The sonic tips may fracture unexpectedly and without warning. In order to reduce the risk of fracture, make sure to stick to the correct power setting (see instructions for use) and observe all recommendations for use. Never allow the instruments to jam. Do not use instruments for leverage. Avoid excessive contact pressure. Use sonic tips with the utmost care as even minor mechanical damage can provoke the fracture of the instruments.
- Risk of ingestion! Fracture of the sonic tip leads to the risk of the tip being accidentally swallowed. Take all appropriate steps to prevent this.
- Premature wear! Work on hard dental substance can cause premature wear of the diamond coating. Likewise, excessive contact pressure can lead to damage to the working part or the diamond coating. Inspect instruments after every use for signs of wear or damage. Exchange if necessary.
- Excessive heat generation! Excessive contact pressure and/or worn or blunt instruments can cause hazardous heat generation which increases the risk of thermal necrosis (damage to tissue). Excessive contact pressure therefore has to be avoided. Damaged sonic tips (changed shape, corroded surfaces, void spots or defect thread etc.) have to be replaced. Incorrect use can lead to the above detailed undesirable effects as well as to injury to the user, the patient or third persons.

Safety control:
Fracture may occur due to the tips being subjected to continuous or sudden stress (for example if the tips are dropped on the floor or through a mechanical change of the original shape). Check the tips prior to each use to make sure they are safe to use. To do this, apply slight pressure with your thumb and index finger (fig. 1). For additional safety, load the tips mechanically with approx. 5N (0,5 kg), without function (fig. 2). Always wear protective gloves when checking, inserting or removing the sonic tips to avoid risks of injury and infection. We recommend that the tips be exchanged regularly.

4. Appropriate use
The sonic tips are supplied in a non-sterile condition and have to be prepared prior to first use.
Sonic tips with a "Q" in the first part of their reference number can be used in the following power systems:
- In the Komet sonic hand piece Quick SFQ2008L
- In the SONICflex™ Quick sonic hand piece made by the co. KaVo Series 2008/S/L/LS
The above listed sonic hand pieces are provided with a quick MULFflex® or Sirona® coupling (2008S/LS). The required tip is screwed into the hand-piece with the tip changer SF4887 which is supplied with the hand-piece or alternatively with the SF1975 in a clockwise direction until a clicking sound can be heard. When using the tip changer, make sure that the tip faces inwards, i.e. towards the semi-radius, during insertion. To remove the tip, place the tip changer onto the tip and remove turning in an anti-clockwise direction.
Note: The above observations the directions in the instructions supplied with your sonic hand piece.

5. Recommendations for use
- Start up sonic tip prior to application to the tooth, providing plenty of spray cooling.
- During stripping and shaping, the above listed sonic hand pieces have to be used at power setting 1
- The use of power settings 2 and 3 not permitted when working with the above described tips (see "warnings").
The below described sonic tips are designed for use in the front and lateral region:
The tips SFQD1F and SFQM1F - marked "strip" - are intended for use in orthodontics and crown core preparation.
The tips SFQD2F and SFQM2F - marked "shape" - are intended for use in orthodontics and work on composite fillings.
Work on composite fillings:
- Apply low contact pressure only. Excessive contact pressure reduces the vibration of the sonic tip and consequently its effectiveness.
- After conclusion of the treatment with the sonic tips, continue the therapy in the usual manner.

6. Maintenance/care
Place sonic tips in a cleaning/disinfection tank (Fräsator) filled with a suitable detergent/disinfectant (e.g. Komet® DC1“) immediately after use in the mouth to prevent drying of residues on the instruments (protein fixation). Hold sonic tips at an angle when placing them into the liquid to ensure that all hollow spaces are reached. It is recommended to reprocess the instruments within one hour of use at the very latest. The instruments should be in the cleaning/disinfection tank when transported to the site where the reprocessing is to take place.

7. Reprocessing
The sonic tips are supplied non-sterile and have to be processed prior to first use.
Manual pre-cleaning (prior to mechanical and manual reprocessing):
1. Remove sonic tip from the cleaning /disinfection tank immediately before reprocessing.
2. Remove stubborn contamination with a nylon brush under running water, turning the sonic tip constantly.
3. Check patency of the cooling passage by guiding the nozzle cleaner (accessory of the sonic and piece) through the passage. Replace instrument if the passage is obstructed.

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Instructions for use
Please read carefully and retain for future reference.

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Mechanical reprocessing (validated):

- Rinse sonic tip thoroughly under running water to prevent any residues of the detergent /disinfectant from getting into the washer/disinfectant.
- Unscrew the nozzle on the injector bar of the washer/disinfectant. Then screw the external thread of the rinse adapter SFQ1978 in the silicon holders on the injector bar. Now screw the sonic tip in the internal thread at the top of the rinse adapter. Observe the directions in the instructions for use provided with the rinse adapter SFQ1978 (fig. 3).
- Put chemical detergent into the washer/disinfectant, following the indications on the label and the instructions of the manufacturer of the washer/disinfectant.
- Start the Vario TD programme including thermal disinfection. Thermal disinfection takes place allowing for the A₀ value and observing national provisions (prEN/ISO 15883).
- On completion of the cycle remove sonic tip from the washer/disinfectant and dry (preferably with compressed air as recommended by the Commission for Hospital Hygiene and Infectious Disease Prevention of the Robert Koch Institute). When drying the inside of the instruments, please make sure to hold the compressed air gun close to the cooling passage, in order to guarantee that sufficient air flows through the sonic tip.
- Visual examination to ensure that the instrument is clean and undamaged. If after mechanical reprocessing there are still visible residues of contamination, repeat the cleaning and disinfecting process until no visible contamination is left.

Manual reprocessing (alternative):

- Place sonic tip onto a 10ml syringe and rinse with detergent/disinfectant. In case of visible contamination emerging from the perforation, rinse again with detergent/disinfectant.
- Remove stubborn contamination with a nylon brush under running water, turning the sonic tip constantly.
- Rinse sonic tip thoroughly with running water.
- Visual examination to ensure that the instrument is clean. If there are still visible residues of contamination, repeat the cleaning process until no visible contamination is left.
- Place sonic tip in a suitable container (e.g. Komet 9952) into the ultrasonic device or instrument bath filled with detergent disinfectant.
- During chemical disinfection in the ultrasonic device or instrument bath, observe the instructions of the manufacturer regarding concentration and immersion time. Be sure to observe the full correct immersion time which does not start until the last tip has been immersed. Attention: do not exceed 45°C (risk of protein coagulation)!
- On completion of the immersion time, rinse sonic tip thoroughly with suitable water (preferably with demineralised water to avoid residues of lime). Rinse the perforation of the sonic tip with at least 10 ml of demineralised water by means of a cannula, to ensure that no residues of the detergent are left in the cooling passage.
- Dry sonic tip (preferably with compressed air as recommended by the Commission for Hospital Hygiene and Infectious Disease Prevention of the Robert Koch Institute). When drying the inside of the instruments, please make sure to hold the compressed air gun close to the cooling passage, in order to guarantee that sufficient air flows through the sonic tip.
- Visual examination to ensure that the instrument is clean and undamaged. If there are still visible residues of contamination, repeat the cleaning and chemical disinfecting process until no visible contamination is left. According to the recommendations of the Commission for Hospital Hygiene and Infectious Disease Prevention of the Robert Koch Institute, the further reprocessing should be carried out mechanically. Information on the reprocessing of instruments in compliance with the standard EN ISO 17664 can be downloaded from our website www.kometdental.de or requested directly from the manufacturer Gebr. Brasseler.

Thermal disinfection in the autoclave:

Classed as risk group semi-critical B, these sonic tips are suitable for thermal disinfection in a steam sterilizer, in an unwrapped condition inside a suitable container (e.g. Komet 9952).
For this procedure, the following instructions have to be observed:
Steam sterilisation using a vacuum process at 134°C in a device that complies with the provisions of DIN EN 13060; with validated processes.
- Fractionated pre-vacuum (type B)
- Sterilisation temperature: 134°C
- Hold time: at least 5 minutes (full cycle)
- Drying time: at least 10 minutes

In order to prevent staining and corrosion, the steam must be free of particles. The recommended limits for particle contents in feed water and condensed steam are defined by standard DIN EN 13060. Make sure not to exceed the maximum capacity of the sterilizer when sterilizing several instruments. Follow the instructions of the device manufacturer.

Universally valid notes:

Observe the legal provisions regarding the reprocessing of medical products valid in your country (e.g. www.rki.de). The manufacturer confirms that the above detailed reprocessing methods are suitable for preparing the above named instrument group to enable their reuse. The person carrying out the reprocessing is responsible for ensuring that the applied method is carried out with appropriate equipment, materials and trained personnel at the reprocessing site and that it actually achieves the desired result. To guarantee this, routine controls of the validated mechanical and/or manual preparation methods are necessary. Any deviation from the above detailed process must be carefully checked by the operator to ensure effectiveness and to avoid possible adverse consequences.

Accessories:

Komet tip changer SF4887 (supplied with the hand-piece)
Komet tip changer SF1975

For reprocessing in the Miele washer/disinfectant: Komet rinse adapter SFQ1978 for Komet sonic tips Komet sterilization container 9952

8. Safety and liability:

Only use the instruments according to their intended purpose. The non-observance of these safety instructions can lead to injuries or damage to the power system. It is the responsibility of the user to check the products prior to use to ensure that they are suited for the intended purpose. In case of contributory negligence by the user, Gebr. Brasseler partially or totally declines liability for all resulting damages, particularly due to non-observance of our recommendations for use or warnings as well as inadvertent misuse by the user.

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