Komet Dental
Gebr. Brasseler GmbH & Co. KG
Trophagener Weg 25
32657 Lemgo · Germany
Telefon +49 (0) 5261 701-0
Telefax +49 (0) 5261 701-289
info@kometdental.de
www.kometdental.de

EndoExplorer
© Gebr. Brasseler GmbH & Co. KG
312021 - Rev 2020-06

Wurzelkanalaufbereitung, Root canal preparation,
Traitement endodontique, Preparación del canal
radicular, Preparazione canalare

STERILE strahlensterilisiert (Sterile Instrumente sind auf dem
Verpackungsetikett entsprechend gekennzeichnet),
Sterilized using irradiation (Sterile instruments are
marked as such on the label of the packaging),
Stérilisation par irradiation (Les instruments
stériles sont marqués comme tels sur l'étiquette
de l'emballage), Esterilización por radiación (Los
instrumentos estériles son marcados como tales
sobre la etiqueta del embalaje), Sterilizzazione per
mezzo di radiazione (Gli strumenti sterili sono
marcati come tali sull'etichetta dell'imballaggio)

Verpackungseinheit, Packing unit, Unité d'emballage,
Unidad de embalaje, Unità d'imballo

verwendbar bis (nur zutreffend bei sterilen
Instrumenten), Use-by date (only applies to sterile
instruments), A utiliser avant le (ne s'applique qu'aux
instruments stériles), Utilizar antes de (se aplica
solamente a instrumentos estériles), Data di scadenza
(riguarda soltanto strumenti sterili)

bei beschädigter Verpackung nicht verwenden
(nur zutreffend bei sterilen Instrumenten),
Do not use in case of damaged packaging
(only applies to sterile instruments),
Ne pas utiliser si l'emballage est endommagé
(ne s'applique qu'aux instruments stériles),
No utilizar si el embalaje está dañado
(se aplica solamente a instrumentos estériles),
Non utilizzare se l'imballo è danneggiato
(riguarda soltanto strumenti sterili)

EndoExplorer				
				Hartmetall Tungsten carbide Carbure de tungstène Carburo de tungsteno Carburo di tungsteno
EX1S.204.007 ○ _{opt.} 2.000 min ⁻¹ ○ _{max.} 40.000 min ⁻¹	EX1.204.007 ○ _{opt.} 2.000 min ⁻¹ ○ _{max.} 40.000 min ⁻¹	EX1L.204.007 ○ _{opt.} 2.000 min ⁻¹ ○ _{max.} 40.000 min ⁻¹	EX1.310.007 ○ _{opt.} 10.000 min ⁻¹ ○ _{max.} 160.000 min ⁻¹	
EX2S.204.011 ○ _{opt.} 2.000 min ⁻¹ ○ _{max.} 40.000 min ⁻¹	EX2.204.011 ○ _{opt.} 2.000 min ⁻¹ ○ _{max.} 40.000 min ⁻¹	EX2L.204.011 ○ _{opt.} 2.000 min ⁻¹ ○ _{max.} 40.000 min ⁻¹	EX2.310.011 ○ _{opt.} 10.000 min ⁻¹ ○ _{max.} 160.000 min ⁻¹	

Tab. 1

EndoExplorer

1. Indikation

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EndoExplorer

1. Indication

Removal of dentin during endodontic access preparation.

Indications + : Preferably - : Not suitable	EX1	EX2
Straightening of the walls of access cavities		+
Exposure of the anatomy of the pulp chamber	+	
Conservative opening of root canal entrances	+	
Removal of dentin overhangs in the region of the canal entrance		+
Exposure of obliterated root canals	+	
Exposure of posts and instrument fragments	+	
Removal of gutta-percha and sealer residues from the walls of the cavity after completed root filling		+

2. Contra indications

Isthmuses situated too deeply in case of narrow coronal canals, isthmuses or obliterations in curved canal sections.

3. Appropriate use

Instruments that are supplied non-sterile have to be prepared prior to use. Capable of cutting in both axial and radial direction, the high-efficiency EndoExplorers EX1 and EX2 were developed for the exposure of all anatomical structures at the bottom of the pulp chamber that are relevant for the endodontic treatment. They are used practically without pressure and without water spray at an optimum speed of 2,000 rpm (contra-angle) or 10,000 rpm (FG). When removing dentin, make sure to proceed step by step. Prepare the canal access to a tapered shape to avoid possible fracture caused by the instrument jamming in the canal. The long, slim instrument neck guarantees almost unimpaired vision on the operation site. We recommend working with a surgical microscope in order to ensure optimum orientation and precise substance removal and to minimize the risk of an iatrogenic perforation.

When using the EndoExplorer, bear in mind/make sure that
- the instrument neck or the working part do not jam in the canal
- the instrument is used practically without contact pressure
- the instrument is not used for leverage
- the instrument does not jam
- the rotary instrument is only used for eliminating tissue during the removal of fragments
- unobstructed view of the tip of the working part is guaranteed at all times during substance removal
- the instrument is used with visual aids (ideally with a surgical microscope)
Non-observance of these recommendations can lead to excessive wear of the instrument (increased risk of fracture).

4. Further treatment

After initial opening of the root canal with the EndoExplorer, reaming with the Opener OP10L19.204.030 and following the creation of a glide path, for example mechanically with the PathGlider PG03L25.204.015, the root canal can be prepared with the NiTi file systems F360, F6 SkyTaper or R6 ReziFlow.

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Instructions for use

Please read carefully and retain for
future reference.

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5. Preparation

The instruments come in sterile packaging and do not require processing prior to first use, provided that the sterile packaging is intact. Komet Dental regrets being unable to guarantee sterility once the instruments have undergone a reprocessing cycle or in case of damaged packaging.

Manual pre-cleaning of the instruments

After use on gutta-percha:

Immerse the instrument in "Guttasolv" and leave to react for 2 minutes. Clean instrument with a cloth. Re-immers the instrument in "Guttasolv" and allow to react for 1 further minute. Clean thoroughly with a nylon brush 9873 and rinse under running water. Control the instrument visually to make sure that it is clean and undamaged. In case of residual contamination, repeat the above described cleaning steps until all visual contamination has been removed.

The instruments have to be treated with anti-corrosive cleaning and disinfection agents suited for rotary instruments (e. g. Komet DC1). Regarding recommendation for use (immersion time, concentration, suitability) of the cleaning and disinfection agents see manufacturers' instructions. cleaning and disinfection agents must be rinsed off thoroughly with water and the instruments must be dried carefully (preferably with compressed air according to the recommendations of the Commission for Hospital Hygiene and Infectious Disease Prevention of the Robert Koch Institute). Do not store instruments in a wet or humid condition for a prolonged period of time. Make sure that instruments do not come in contact with each other during ultrasonic cleaning. The instruments can be reprocessed in the thermos disinfectant, provided that a mildly alkaline cleaning and disinfection agent is used that is suitable for rotary instruments. Inspect cleaned instruments visually. Do not reuse damaged or blunt instruments. These have to be discarded.

Sterilization:

Sterilization of the instruments is carried out in the autoclave at 134°C. Observe instructions provided by the manufacturer. After sterilization, check instruments for defects on the surface (corrosion). Do not reuse corroded instruments. The operator of medical products is responsible for seeing that proper treatment is carried out by qualified personnel, using the appropriate materials and suited equipment, as recommended by the Commission for Hospital Hygiene and Infectious Disease Prevention of the Robert Koch Institute. Information with regard to proper reprocessing of instruments according to DIN EN ISO 17664, the present instructions for use, as amended, and more detailed information can be downloaded from our web site www.kometdental.de or requested from the manufacturer Gebr. Brasseler.

6. Storage

Store sterile packed instruments in their original packaging, away from UV radiation and heat in a clean and dry environment. Do not store in the same room as solvents or chemicals.

7. Disposal

The instruments have to be disposed of in tight, non-breakable and non-pierceable containers (protection from contamination).

8. Safety and liability

See also "Instructions for use and safety recommendations" in the current Komet dental catalogue. It is the responsibility of the user to check the products prior to use as to whether they are suited for the intended purpose. The user is responsible for the correct use of the instruments. In case of contributory negligence by the user, Gebr. Brasseler GmbH & Co. KG partly or totally declines liability for all resulting damages, particularly due to non-observance of our instructions for use as well as misuse by the user.

9. Procedure in case of serious adverse events

In case of serious adverse events, please proceed in line with the legal regulations in your country.

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