Komet Dental
 Gebr. Brasseler GmbH & Co. KG
 Trophagener Weg 25
 32657 Lemgo - Germany
 Telefon+49 (0)5261 701- 0
 Telefax+49 (0)5261 701- 289
 info@brasseler.de
 www.brasseler.de

Endo ReStart
 © Gebr. Brasseler GmbH & Co. KG
3 12700 - Rev 2019-02

Wurzelkanalaufbereitung, Root canal preparation,
 Traitement endodontique, Preparación de canal
 radicular, Strumentazione canalare

nur zum Einmalgebrauch, For single use only, Ne pas
 réutiliser (usage unique), No reutilizar (un solo uso),
 Esclusivamente monouso

Verpackungseinheit, Packing unit, Unité d'emballage,
 Unidad de embalaje, Unità d'imballo

strahlensterilisiert (Sterile Instrumente sind auf dem Verpackungsetikett entsprechend gekennzeichnet), Sterilized using irradiation (Sterile instruments are marked as such on the label of the packaging), Stérilisation par irradiation (Les instruments stériles sont marqués comme tels sur l'étiquette de l'emballage), Esterilización por radiación (Los instrumentos estériles son marcados como tales sobre la etiqueta del embalaje), Sterilizzazione per mezzo di radiazione (Gli strumenti sterili sono marcati come tali sull'etichetta dell'imballaggio)

verwendbar bis (nur zutreffend bei sterilen Instrumenten), Use-by date (only applies to sterile instruments), A utiliser avant le (ne s'applique qu'aux instruments stériles), Utilizar antes de (se aplica solamente a instrumentos estériles), Data di scadenza (riguarda soltanto strumenti sterili)

bei beschädigter Verpackung nicht verwenden (nur zutreffend bei sterilen Instrumenten), Do not use in case of damaged packaging (only applies to sterile instruments), Ne pas utiliser si l'emballage est endommagé (ne s'applique qu'aux instruments stériles), No utilizar si el embalaje está dañado (se aplica solamente a instrumentos estériles), Non utilizzare se l'imballo è danneggiato (riguarda soltanto strumenti sterili)

Endo ReStart NITI-Revisionsfeilensystem

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medicament, if required, and filled according to the established methods, for example condensed gutta-percha such as F360 (GPF04), F6 SkyTaper (GPF06) or R6 ReziFlow (GPR06), gutta-percha points, a carrier based filling system – e.g. F360 Fill – or a thermoplastic method.

Endo ReStart System of NITI files for retreatment

- Indication**
The instruments are intended for the mechanical removal of diseased tissue and root filling materials within the course of a root canal treatment.
- Contraindication**
The contraindications are (e.g.):
 - pronounced abrupt curves
 - severely calcified root canals
 - preparation of milk teeth
 Endo ReStart Opener:
 - Do not use in the apical and middle third of the root canal
- Proper use**
General information:
 - The Endo ReStart Opener is intended for repeated use. Control and document how many times the instrument has been used. Do not use for preparing more than 8 canals.
 - The Endo Restart files are intended for single use (one patient). Up to 5 canals (of the same patient) can be treated with one file.
 - Observe the recommendations regarding the use of mechanical NITI instruments and the cleaning and disinfection of the root canal.
 - Keep an eye out for plastic deformation of the instruments during treatment and intermediate cleaning, for example untwisted windings in the region of the working part. Discard affected instruments.
 Prior to use of the files:
 - It is recommended to determine the canal morphology beforehand with the help of radiographs.
 - Insert a solvent, if necessary, for example eucalyptus oil to soften the gutta-percha.
 Correct use of the files:
 - Use the instruments RE10L15 (Endo ReStart Opener), RE05L21 & RE05L25 (Endo ReStart) in a torque limited contra-angle at a speed of 300 rpm (max. speed 500 rpm) and with a torque of 1.8 Ncm.
 - Removal of the root filling in the coronal third of the canal with the Endo ReStart Opener:
 - The Opener is provided with a cutting tip to enable it to penetrate even hard root fillings
 - The toothly only stretches over a short part of the file because its use is limited to the coronal part of the root canal only.
 - The remaining root filling is removed with the Endo ReStart files following the single-length method.
 - The file is provided with a scraping rather than a cutting tip, to enable it to further penetrate hard root fillings.
 - The toothling of the file is provided with pronounced winding to ensure that the file can work its way into the root filling (screw-in).
 - Start up file and insert into the canal. The instruments are used in "picking motion", i.e. in short, constant up-and-down movements. To ensure a long durability of the Endo ReStart files, please make sure that the gutta-percha root filling is removed step by step, applying the file with a minimum of axial force.
 - Clean the instruments mechanically from time to time, for example in the interim bur block 595. Rinse the root canal thoroughly to completely flush out all debris.
- Further treatment**
After the successful removal of the root filling, the root canal can be enlarged manually or mechanically, for example with the NITI file systems F360, F6 SkyTaper or R6 ReziFlow for the mechanical preparation of root canals. After that, the root canal can be provided with a

- Reprocessing**
Endo ReStart files:
The files are sterile packed and do not need to be prepared, i.e. cleaned/disinfected and sterilized, prior to use. The instruments were developed for single use (on just one patient) and have to be discarded after use. The reuse of these products poses a risk of infection and/or the safety of the products can no longer be guaranteed. Do not use instruments inside damaged packaging because they are no longer sterile.
Endo ReStart Opener:
The present instruments are supplied in sterile packaging. They can therefore be used without previous processing, provided that the packaging is intact. The instruments are to be treated with anti-corrosive disinfecting and cleaning agents suited for rotary instruments (e.g. Komet DC1). Regarding recommendation for use (immersion time, concentration, suitability) of the cleaning or disinfection agents see manufacturers' instructions. Cleaning and disinfection agents must be rinsed off thoroughly with water and the instruments must be dried carefully (preferably with compressed air according to the recommendations of the Commission for Hospital Hygiene and Infectious Disease Prevention of the Robert Koch Institute). Do not store instruments in wet or humid condition for a prolonged period of time. Make sure to remove all traces of NaOCl and limit the contact period to no more than a few minutes. Make sure that instruments do not come in contact with each other during ultrasonic cleaning. The instruments can be reprocessed in the thermo disinfecter, provided that the agent used is suitable for rotary instruments (mildly alkaline detergent). Inspect cleaned instruments visually. Do not reuse damaged or blunt instruments. These have to be discarded.
Sterilization:
The sterilization of the instruments is carried out in the autoclave at 134°C. Observe the instructions provided by the manufacturer. After sterilization, check instruments for surface defects (corrosion). Do not reuse corroded instruments. The operator of medical products is responsible for ensuring that proper treatment is carried out by qualified personnel, using suitable materials and suited equipment, as recommended by the Commission for Hospital Hygiene and Infectious Disease Prevention of the Robert Koch Institute. Work instructions with regard to the correct reprocessing of instruments according to DIN EN ISO 17664 can be downloaded from our web site www.kometdental.de or requested from the manufacturer Gebr. Brasseler.
- Storage**
Store packed instruments away from UV radiation and heat in a clean and dry environment. Do not store in the same room as solvents or chemicals.
- Disposal**
The instruments have to be disposed of in tight, non-breakable and non-pierceable containers (protection from contamination).
- Safety and liability**
See also "Instructions for use and safety recommendations" in the current Komet Dental catalogue. The user is responsible for checking the products prior to use as to whether they are suited for the intended purpose. The user is responsible for the appropriate use of the instruments. In case of contributory negligence by the user, Gebr. Brasseler partially or totally declines liability for all resulting damages, particularly if these are due to nonobservance of our recommendations for use or warnings as well as inadvertent misuse by the user.
- Procedure in case of serious adverse events**
In case of serious adverse events, please proceed in line with the legal regulations in your country.

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