Komet Dental  
Gebr. Brasseler GmbH & Co. KG  
Trophagener Weg 25  
32657 Lemgo · Germany  
Telefon +49 (0) 5261 701-0  
Telefax +49 (0) 5261 701-289  
info@kometdental.de  
www.kometdental.de

**Implantologievorböhrer 210IK16/19**  
© Gebr. Brasseler GmbH & Co. KG  
**313417 - Rev 2017-10**

Implantologie, Implantology, Implantologie,  
Implantología, Implantologia

Strahlensterilisiert (Sterile Instrumente sind auf dem Verpackungsetikett entsprechend gekennzeichnet), Sterilized using irradiation (Sterile instruments are marked as such on the label of the packaging), Stérilisation par irradiation (Les instruments stériles sont marqués comme tels sur l'étiquette de l'emballage), Esterilización por radiación (Los instrumentos estériles son marcados como tales sobre la etiqueta del embalaje), Sterilizzazione per mezzo di radiazione (Gli strumenti sterili sono marcati come tali sull'etichetta dell'imballaggio)


Verwendbar bis (nur zutreffend bei sterilen Instrumenten), Use by (only applies to sterile instruments), A utiliser avant le (ne s'applique qu'aux instruments stériles), Fecha de caducidad (se aplica solamente a instrumentos estériles), Data di scadenza (riguarda soltanto strumenti sterili)

Bei beschädigter Verpackung nicht verwenden (nur zutreffend bei sterilen Instrumenten), Do not use in case of damaged packaging (only applies to sterile instruments), Ne pas utiliser si l'emballage est endommagé (ne s'applique qu'aux instruments stériles), No utilizar si el embalaje está dañado (se aplica solamente a instrumentos estériles), Non utilizzare se l'imballo è danneggiato (riguarda soltanto strumenti sterili)

**Pilotbohrer mit Innenkühlung · Pilot drill with internal cooling · Foret pilote avec refroidissement interne · Fresa piloto con refrigeración interna · Fresa pilota con raffreddamento interno**

⌀<sub>int.</sub> 800 - 1.000 min<sup>-1</sup>/rpm ⌀<sub>max.</sub> 6.000 min<sup>-1</sup>/rpm

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Stainless steel pilot drill for use in implantology (with internal cooling)

**1. System overview**

See table 1  
The working parts of instruments with internal cooling are provided with a bore hole with orifice for feeding the cooling liquid. A Teflon seal prevents the reflux of the liquid.

**2. Indication**

- Creation of a pilot hole during the preparation of the implant site. If a template is used make sure that the instrument does not jam in order to prevent possible fracture of the instrument.  
- Initial determination of the position, direction and depth of the implant.

**3. Contraindications**

- Poor bone quality

**4. Clinical sequence**

- Select the required size of the pilot drill and determine the drilling depth.  
- Determine the size and the exact position of the implant with the help of radiographs (bitewing radiograph and orthopantomogramme).  
- In order to establish the exact position and the depth of the borehole, a computer tomographic examination in addition to the bitewing radiograph and the orthopantomogramme is recommended.  
- Examine the adjacent structures.  
- In order to prevent any danger to adjacent structures, carefully check the area surrounding the operation site.  
- Anaesthesia the site where the instrument is to operate.  
- Cut the mucosa and periost to the bone and expose the operative site.  
- Optional: Mark the exact position with the help of a round bur, for example H141A.205.023.  
- Apply the pilot drill at the chosen position.  
- Drilling should be carried out in intermittent mode, providing constant exterior cooling with sterile physiological saline solution to prevent thermal damage to the bone tissue and to ensure proper chip removal. Make sure that the drill neither jams nor gets stuck during use (increased risk of fracture). Prepare to the desired depth at a speed of 800-1.000 rpm (green contra-angle, reduction at least 10:1), exerting low pressure. The indicated speed should be observed in order to avoid fracture of the instrument. Do not exceed a torque of 30 Nm.

**5. Further treatment**

The instrumentation of the implant system used is subsequently used to enlarge the implant site to the defined implant diameter.

**6. Recommendations for maintenance**

Instruments supplied in a non-sterile condition have to be disinfected, cleaned and sterilized prior to first use.  
- Do not allow residues from the operation (blood, secretions, tissue) to dry.  
- Immerse the drills in disinfecting solution immediately after the operation.

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**Instructions for use**

Please read carefully and retain for future reference.

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- Disinfect the drills with cleaning and disinfecting agents for rotary instruments (e.g. Komet DC1).

For recommendations for use (immersion time, concentration, suitability) of the cleaning and disinfecting agents please refer to the instructions provided by the manufacturer of the agents used. The instruments can be reprocessed in the thermodesinfector, provided that the agent used is suitable for rotary instruments (mildly alkaline detergent). Prior to mechanical preparation, remove instrument from cleaning/disinfection tank and check the patency of the cooling orifices by guiding the cleaning wire through the passages. Place instrument on a 10 ml cannula and rinse with demineralised water. Cleaning and disinfecting agents must be rinsed off thoroughly with water and the instruments must be dried carefully (preferably with compressed air according to the recommendations of the Commission for Hospital Hygiene and Infectious Disease Prevention of the Robert Koch Institute). Do not store instruments in wet or humid condition for a prolonged period of time. Make sure that instruments do not come in contact with each other during ultrasonic cleaning. Inspect cleaned instruments visually. Do not reuse damaged or blunt instruments. These have to be discarded due to an increased risk of fracture

The sterilization of the instruments is carried out in the autoclave at 134°C according to the known methods. Observe the instructions provided by the manufacturer. The operator of medical products is responsible for ensuring that proper treatment is carried out by qualified personnel, using suitable materials and suited equipment, as recommended by the Commission for Hospital Hygiene and Infectious Disease Prevention of the Robert Koch Institute. Work instructions with regard to proper reprocessing of instruments according to DIN EN ISO 17664, the present instructions for use and more detailed information, as amended, can be downloaded from our web site [www.kometdental.de](http://www.kometdental.de) or requested from the manufacturer Gebr. Brasseler.

**7. Storage**

Store sterile packed instruments in their original packaging, away from light and heat in a clean and dry environment. Do not store in the same room as solvents or chemicals.

**8. Safety and liability**

The cutting of hard bone material may lead to premature blunting of the blades. The drills therefore have to be checked for blunt blades or other damage after each use. Exchange if necessary. Be sure to observe the indicated speed in order to avoid instrument fracture. Attention: Risk of injury due to the sharp blades of the bur! Risk of injury caused by jamming and slipping of the bur! The user is responsible for checking the products prior to use as to whether they are suited for the intended purpose. In case of contributory negligence by the user, Gebr. Brasseler partially or totally declines liability for all resulting damages, particularly if these are due to non-observance of our recommendations for use or warnings as well as inadvertent misuse by the user.