	<p>Komet Dental Gebr. Brasseler GmbH & Co. KG Trophagener Weg 25 32657 Lemgo · Germany Telefon +49 (0) 5261 701-0 Telefax +49 (0) 5261 701-289 info@kometdental.de www.kometdental.de</p>
---	---

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EN

Sonic tips with Quick connection for the treatment of furcations

1. Indication

These special sonic tips are used for the gentle cleaning and smoothing of roots and furcations as part of treatments both with and without the creation of a flap. Depending on the severity of the deposits on the furcation in the upper or lower jaw, these sonic tips can be used on their own or in combination with other instruments for work on root surfaces during various treatment options or for removing dental plaque (biofilms) as part of conservative treatments.

2. Contra-indication

These sonic tips for the treatment of furcations are not suitable for use on dental crowns (enamel surfaces), dental restorations and implant surfaces.

3. Warnings

The useful life of the sonic tips is limited. They therefore have to be used with utmost care.

Attention:

- Risk of fracture! Spontaneous and premature fracture must be expected. The sonic tips may fracture unexpectedly and without warning. In order to reduce the risk of fracture, make sure to stick to the correct power setting (see instructions for use) and observe all recommendations for use. Never allow the instruments to jam. Do not use instruments for leverage. Avoid excessive contact pressure. Use sonic tips with the utmost care as even minor mechanical damage can provoke the fracture of the instruments.

- Risk of ingestion! Fracture of the sonic tip leads to the risk of the tip being accidentally swallowed. Take all appropriate steps to prevent this.

- Premature wear! Due to the nature of their range of application, the functional parts of these sonic tips will always be subject to moderate wear, even when the tips are applied correctly and according to their intended use. Excessively worn tips have to be exchanged.

- Excessive heat generation! Excessive contact pressure and/or worn or blunt instruments can cause hazardous heat generation which increases the risk of thermal necrosis (damage to tissue). Excessive contact pressure therefore has to be avoided. Damaged sonic tips (changed shape, corroded surfaces or defect thread etc.) have to be replaced.

- Damage to hard tissues! Excessive contact pressure might damage the root surfaces. Excessive contact pressure therefore has to be avoided. This applies in particular to the periodontal tip which may only be used with low contact pressure (0.5 N).

- Incorrect use can lead to the above detailed undesirable effects as well as to injury to the user, the patient or third persons.

Safety control:

Fracture may occur due to the tips being subjected to continuous or sudden stress (for example if the tips are dropped on the floor or through a mechanical change of the original shape). Check the tips prior to each use to make sure they are safe to use. To do this, apply slight pressure with your thumb and index finger (fig. 1). For additional safety, load the tips mechanically with approx. 10 N (1 kg), without function (fig. 2). Always wear protective gloves when checking, inserting or removing the sonic tips to avoid risks of injury and infection. We recommend that the tips be exchanged regularly.

Increased oscillation: Please note that the oscillations of these sonic tips can exceed 200µm when used at the authorized power level. The increased oscillations do not constitute a safety relevant aspect, as long as the power settings specified in paragraph 5 (recommendations for use) are observed and the limit value concerning the air supply provided by the authorized power system used is not exceeded.

4. Appropriate use

Sonic tips delivered non-sterile have to be prepared prior to first use. For treatments with and without flap, the sonic tips can be used in the following power systems:

- In the Komet sonic hand piece Quick SFQ2008L (MULTIFlex™ coupling) and SFQ2008/LS (Sirona® coupling)

- In the SONICflex™ Quick sonic hand piece made by the Co. KaVo Series 2008/S/L/LS

The above listed sonic hand pieces are provided with a quick MULTIFlex™ or Sirona® coupling (2008S/LS).

The required tip is screwed into the hand-piece with the tip changer SF4887 which is supplied with the hand-piece or alternatively with the SF1975 in a clockwise direction until a clicking sound can be heard. When using the tip changer, make sure that the tip faces inwards, i.e. towards the semi-radius, during insertion. To remove the tip, place the tip changer onto the tip and remove turning in the anti-clockwise direction. Please observe the directions in the instructions supplied with your sonic hand piece.

5. Recommendations for use

- The accessibility of the furcation should be checked according to the severity of the case. To this end, the entrance of the furcation should be slightly enlarged and finely smoothed. The instrument has to be guided across the root surface without pressure (choose the lowest power level). This applies equally to therapies with and without the creation of a flap.

- If an open treatment after the creation of a flap is intended, it is absolutely essential to supply sterile cooling agent. This is facilitated by the use of the cooling adapter SF1979 (see chapter "Cooling"). If you are planning a closed treatment, i.e. without flap, the cooling adapter can be used to prevent secondary infections. In this case, the use is optional.

- Start the sonic tip prior to application to the tooth, providing plenty of spray cooling.

- Power level 1 is the standard choice for all of the above listed power systems. The sonic tip is not authorized for use with power levels 2 and 3 (see warnings).

- Use with low contact pressure (0.5 N) and proceed parallel to the root surface. Excessive contact pressure reduces the vibration of the sonic tip and consequently its effectiveness. Non-compliance can lead to negative effects as described in the warnings.

Cooling:

The sonic tips are cooled with at least 50 ml/min of an externally fed, medically cooling liquid. Make sure to turn off the water supply of the dental treatment unit. It is absolutely essential to thoroughly and continuously irrigate the site at all times during the operation. The sterile cooling liquid can be fed manually, for example by the assistant at the treatment chair by means of a larger piston syringe with a blunt needle. We recommend feeding the sterile cooling liquid through the cooling adapter SFQ1979 (fig. 3). One end of the cooling adapter is screwed into the sonic hand piece and the sonic tip is screwed into the other end of the cooling adapter. The flexible tube supplying the cooling agent is placed onto the small tube of the cooling adapter (fig. 3). The sonic tip can be removed from the cooling adapter by means of the mounting wrench 566.

Controlled supply of the cooling agent by means of a suitable surgical power system is preferable, for example with the KaVo INTRASurg 1000. There are a number of compression systems for pressure supported transport of the cooling agent available on the market, for example:

System OMNIA, manual supply of sterile cooling agents by means of a compression pump, available in Germany from the company Implantis Medical Devices. They can be contacted as follows:
From Germany: (0180) 5 90 75 90, free fax: (00800) 05 05 51 23.

From Austria: (0049 180) 5 90 75 90, free fax: (00800) 05 05 51 23.

H2200 Manual cooling system with compression pump (Omnia no. 30E0008.00); H2204 Manual cooling system with flux control, to be attached to the hand piece (Omnia no. 32F0157), www.omniasr.com

System Aqua-Press sterile irrigation system made by the company Ronvig Dental A/S, Gl. Vejlevej 59, DK-8721 Dagaard. Tel.: 0045-70 23 34 11, Fax: 0045-76 74 07 98, www.ronvig.com.

6. Maintenance/care

Place unscrewed sonic tips and cooling adapter in a cleaning/disinfection tank (Fräsator) filled with a suitable detergent/disinfectant (e.g. Komat DC1) immediately after use in the mouth to prevent drying of residues on the instruments (protein fixation). Hold sonic tips at an angle when placing them into the liquid to ensure that all hollow spaces are reached. It is recommended to reprocess the instruments within one hour of use at the very latest. The instruments should be in the cleaning/disinfection tank (Fräsator) when transported to the site where the reprocessing is to take place.

7. Reprocessing

Sonic tips delivered non-sterile have to be prepared prior to first use.

Manual pre-cleaning (prior to mechanical and manual reprocessing):

1. Remove sonic tip with or without cooling adapter from cleaning/disinfection tank (Fräsator) immediately before reprocessing.

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© Gebr. Brasseler GmbH & Co. KG